

Product datasheet for **SC315810**

FBXO46 (NM_001080469) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FBXO46 (NM_001080469) Human Untagged Clone
Tag:	Tag Free
Symbol:	FBXO46
Synonyms:	20D7-FC4; Fbx46; FBXO34L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC315810 representing NM_001080469.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGACCGTGGGAGCCTCTGCCCTCCAGCTATGGTGCCCCGGCCCTTTGGCACCTACTCACAGAAC
CAGCCACGCCCGCCTTCTGCGGCCCTCAAGCCATCAGCCTGCCCTGAGCCTGGTGGCGGGCCGAGCCA
GACCATGGGCCTGCCACTCAGAGAACACACCACCCGCCTTGGCCACTGAGGTCCCTGCCTCCAGCCG
GCTCCGCTCCTCTCAGCAGCAGCTGCTGGTGATGAGGGTCGAGTCTCCTGGACACGTGGTATGTATC
AAGCCCGGAATACAAAGGAGAAGGTGGCCTTCTTTGTGGCCACCAGTGTGGTGGGGCAGCCGGGCC
AGCTCCATGAAGGTCAAGGGGCACTGGGGCAGCGATAGCTCCAAGGCCAAGCGCGGAGGCGCTGTCTT
GACCCACCAAGGCTCCTCCTGACCCAGGGGGCCGGGAGGGCCCCCTGCTGCTGAGGAGGGCCCCGCC
TCAGCCGGTGAGGACGTGGACCTGCTCTGTGGCCGAGATGGTGGCCCTGGTGGAAACAGCGGGCAGCC
CTGGCCCTGCAGAGCTACCCACGACCGACCACCCAGCGCCTGTAGTCTTTGTGTCGGCCGAGCAAGGT
GGACCGCCAAGGGGGTGGGATCTGAACGGCGGTCCGGTGGTGGGACTGCAGCCGTGTAGCCGAGGCC
GTGGCCCACTTTGAAGCGCAGAGGGACAGCCCTCCACCAAGGGCCTCCGCAAGGAAGAGCGGCCGGG
CCAGGCCCTGGGGAGGTGCGCATCGCCTTCCGCATCTCCAACGGCCGGGAGCCCCGTGCACCAGACAGC
GGCCTGCCAGTGGGGCGGGGGCAGGCCCTGGTTGTGCCTACCCTGGCAGCCAGGTCTGGGGCCCGA
GCCAAGGACAAGATCACATGTGACTTATACCAGCTCATCAGCCCTCGCGGGATGCCCTCCCCAGCAAC
GTGGAGTTCCTGCTGGCCAGGGCGGATGAGGCCAGTGAGGGTGACAGCCCGGCACCCGCCAGGCCTGAG
GACACTCCCCGGCGCCCCCTCCGCCCCCTGCCGGGACTGCGGAGCGTCAGGCTTCCACGTGGACGTG
GTGGTGACGGCGTGGTAGATGAGTGCATCTTCTTTGGCAAGGACGGCACCAAGAAGCTGAAGGAGGAG
ACTGTGTGCCTGACGGTCAGCCCGGAGGAACCGCCGCTCCGGGCCAGCTCTTCTTCTCCAGAACCGC
GGCCCGGACGGGCCCGGAGCCACCCCGGCCGACTCCCGGCCACTGCGCCCGGCCAGACGATGCC
GAGGGCACGGCGGACACCTCCCTGTGCCCTTGTACCGGCACGTGTCGCACGACTTCTAGAGATCCGC
TTCAAGATTACGCGGCTGCTGGAGCCGCGACAGTACATGCTGCTGCTGCCGAGCACGTGCTGGTCAAG
ATCTTCAGCTTCCCTGCCACGCGCGCTGGCCGCCCTCAAGTGCACCTGCCACCACTTCAAGGGCATC
ATCGAGGCGTTTGGCGTGGGGCCACAGACTCGCGCTGGAGCCGAGACCCGCTCTACCGGATGATCCG
TGCAAACAGTGCCGAAGAGATACGAGAAGGGCGACGTGTCGCTCTGCCGCTGGCACCCCAAGCCCTAC
CACCATGACCTGCCTTACGGACGTTCTACTGGATGTGCTGCCGTGAGCCGACCCGAGACTCCCGGC
TGCCGCTGGGCCTCCACGATAACAACCTGGGTGCTGCCCTGCAATGGGCCAGGCGGGGGCCGGCCGGC
CGGGAGGAGGGGAGGTGA
ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001080469

Insert Size: 1812 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001080469.1](#)

RefSeq Size: 2990 bp

RefSeq ORF: 1812 bp

Locus ID: 23403

UniProt ID: [Q6PJ61](#)

Cytogenetics: 19q13.32

Protein Families: Druggable Genome

MW: 64.6 kDa

Gene Summary: Members of the F-box protein family, such as FBXO46, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Mar 2008]