

Product datasheet for **SC315820**

ZNF433 (NM_001080411) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF433 (NM_001080411) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF433
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC315820 representing NM_001080411.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGATGTTTCAGGATTCAGTGGCCTTTGAGGATGTGGCTGTGACCTTCACCCAAGAGGAGTGGGCTTTG
CTGGATCCTTCCAGAAAAATCTCTGTAGAGATGTGATGCAAGAAACCTTCAGGAACCTGGCCTCTATA
GGGAAAAAATGGAAACCCAGAACATATATGTAGAGTACGAAAAATCTAAGGAGAAACCTAAGAATTGTG
GGAGAGAGACTCTTTGAAAGTAAAGAAGGTCATCAGCATGGAGAAATTTTGACCCAGGTTCCAGATGAC
ATGCTGAAGAAAAACAACACTGAGTAAATCATGCGAAAGCAGTGTGTATGGAGAAGTAGGCAGTGCT
CATTATCTCTTAATAGGCACATCAGAGATGACACTGGACACAAGGCATATGAGTATCAAGAATATGGA
CAGAAACCATATAAATGTAAATACTGTAAAAACCTTTCAACTGTCTCTCCTCTGTTCCAGACACATGAA
AGGGCTCATAGTGAAGGAACTCTATGTTTGTGAGGAATGCGGAAAAACATTTATTTCCATTCAAAC
CTTCAAAGACACAGGATAATGCACCGTGGAGATGGACCTTATAAGTGAAATTTTGTGGGAAAGCCTTG
ATGTTTCTCAGTTGTATCTTATCCACAAACGAACTCACACTGGAGAGAAACCATATCAATGTAAACAG
TGTGGTAAAGCCTTTAGTCATTCTAGTAGCCTTCGAATACATGAAAGAACTCACACTGGGAGAGAGCCT
TATAAATGTAATGAATGTGGGAAAGCATTCCATAGTTCCACATGCCTTCATGCTCATAAAAGAACTCAC
ACTGGGAGAGCCATATGAATGTAAACAGTGTGGGAAAGCCTTCAGCTCTTCCATTCTTTCAAATA
CATGAAAGAACTCACACGGGGAGAGCCATATGAATGTAAGGAATGTGGAAAAGCATTCAAGTGTCCC
AGTTCTGTTTCGACAGATGAAAGAACCCACTCTAGGAAAAACCCCTATGAATGTAACATTGTGGGAAA
GTATTATCTTATCTTACCAGCTTTCAAACCACTTGGGAATGCACACTGGAGAGATATCTCATAAATGT
AAGATATGTGGGAAAGCCTTTTATTCTCCAGTTCACTTCAAACACATGAAAAAACTCACACTGGAGAG
AAACCTATAAATGCAACCAATGTGGTAAAGCCTTTAATTCTTCCAGTTCTTCCGATATCATGAAAGA
ACTCACACTGGAGAGAAACCTTACGAGTGTAAAGCAATGTGGGAAAGCCTTCAGATCTGCCTCACTCCTT
CAAACACATGGTAGGACTCACACGGGAGAGAAACCCCTATGCATGTAAGGAATGTGGAAAACCATTTAGT
AATTTCTCTTTCTTTCAAATACATGAAAGGATGCACAGAGAAGAGAAGCCGTATGAATGTAAGGGTTAT
GGGAAAACATTAGTTTGGCCAGTTTATTTATAGACATGAAAGGACTCACACTGGAGGAAAAACCTAT
GAATGCAAGCAGTGTGGCAGATCCTTCAACTGTTTCGAGCTCCTTTCGATATCATGGAAGGACTCACACT
GGAGAGAAACCTATGAATGCAAGCAATGTGGAAAAGCCTTCAGATCTGCCTCACAGCTTCAAATTCAT
GGAAGGACTCACACTGGAGAGAAACCTTATGAATGTAAGCAGTGTGGGAAAGCCTTTGGATCTGCCTCA
CACCTTCAAATGCATGGAAGGACTCACACTGGAGAGAAACCTTATGAATGTAAGCAGTGTGGGAAAGTCT
TTTGGATGTGCCTCGCGACTTCAAATGCATGGAAGGACTCACACTGGAGAGAAACCGTATAAATGTAAAG
CAATGTGGGAAAGCTTTTGGATGTCCCTCAAACCTTCGAAGGCATGGAAGGACTCACACTGGAGAGAAA
CCCTATAAATGTAACCAATGTGGTAAAGCTTTAGATGTTCTTCACAACTCAAGTGCATGGAAGGGCT
CACTGCATAGACACCCCAATA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001080411

Insert Size: 2022 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080411.2</u>
RefSeq Size:	2352 bp
RefSeq ORF:	2022 bp
Locus ID:	163059
UniProt ID:	<u>Q8N7K0</u>
Cytogenetics:	19p13.2
Protein Families:	Transcription Factors
MW:	77.2 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest isoform (1).