

Product datasheet for **SC315826**

LRRC4B (NM_001080457) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC4B (NM_001080457) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRC4B
Synonyms:	HSM; LRIG4; NGL-3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC315826 representing NM_001080457.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGC
ATGGCGCGTGCCCGCGGCTCCCGTGCCCCCGCTGCCCGCGGTAGGATGTCCTGGCCCCACGGGGCA
TTGCTCTTCCTCTGGCTCTTCTCCCCACCCCTGGGGCCGGTGGAGGTGGAGTGGCCGTGACGTCTGCC
GCCGGAGGGGGCTCCCCGCCGGCCACCTCTGCCCGTGGCCTGCTCCTGCAGCAACCAGGCCAGCCGG
GTGATCTGCACACGGAGAGACCTGGCCGAGGTCCAGCCAGCATCCCGGTCAACACGCGGTACCTGAAC
CTGCAAGAGAACGGCATCCAGGTGATCCGACGGACACGTTCAAGCACCTGCGGCACCTGGAGATTCTG
CAGCTGAGCAAGAACCTGGTGCGCAAGATCGAGGTGGGCGCCTTCAACGGGCTGCCAGCCTCAACACG
CTGGAGCTTTTGAACCGGCTGACCAGGTGCCACGAGGCCTTCGAGTACCTGTCCAAGCTGCGG
GAGCTCTGGTGCGAACAACCCATCGAGAGCATCCCTCTACGCCTTCAACCGGTGCCCTCGCTG
CGGCGCCTGGACCTGGGCGAGCTCAAGCGGCTGGAATACATCTCGGAGGCGGCCTTCGAGGGGCTGGT
AACCTGCGCTACCTCAACCTGGGCATGTCAACCTCAAGGACATCCCCAACCTGACGGCCCTGGTGCGC
CTGGAGGAGCTGGAGCTGTGGGCAACCGGCTGGACCTGATCCGCCCGGGCTCCTTCCAGGGTCTACC
AGCCTGCGCAAGCTGTGGCTCATGCACGCCAGGTAGCCACCATCGAGCGCAACGCCCTTCGACGACCTC
AAGTCGCTGGAGGAGCTCAACCTGTCCCAACAACCTGATGTCGCTGCCCCACGACCTCTTACGCCC
CTGCACCGCCTCGAGCGGTGCACCTCAACCACAACCCCTGGCATTGCAACTGCGACGTGCTCTGGCTG
AGCTGGTGGCTCAAGGAGACGGTGCCAGCAACACGACGTGCTGCGCCCGCTGTCATGCGCCCGCCGGC
CTCAAGGGGCGCTACATTGGGGAGCTGGACAGTCGCATTTACCTGCTATGCGCCCGTCATCGTGGAG
CCGCCCACGACCTCAACGTACCGAGGGCATGGCTGCCGAGCTCAAATGCCGCACGGGCACCTCCATG
ACCTCCGTCAACTGGCTGACGCCCAACGGCACCTCATGACCCACGGCTCCTACCGGTGCGCATCTCC
GTCTGTCATGACGGCACGCTTAACCTTCAACAACGTACCGTGCAGGACACGGGCCAGTACAGTGCATG
GTGACGAACTCAGCCGGCAACACCACCGCTCGGCCACGCTCAACGTCTCGGCCGTGGACCCCGTGGCG
GCCGGGGGACCGGCAGCGCGGGGGCGGCCCTGGGGGAGTGGTGGTGTGGAGGGGGAGTGGCGGC
TACACCTACTTACCACGGTGACCGTGGAGACCCTGGAGACGACCCGGAGAGGAGGCCCTGCAGCCG
CGGGGGACGGAGAAGGAACCGCCAGGGGCCACGACAGACGGTGTCTGGGGTGGGGCGCGCTGGGGAC
GCGGCCGCGCCTGCCTCGTCTTACCACGGCACCCGCCCCGCGCTCCTCGCGGCCACGGAGAAGGCG
TTCACGGTGCCCATCACGGATGTGACGGAGAAGCCCTCAAGGACCTGGACGACGTATGAAGACCACC
AAAATCATCATCGCTGCTTCTGTGGCCATCACGTTTATGGCCGCGGTGATGCTCGTGGCCTTCTACAAG
CTGCGCAAGCAGCACCAGCTCCACAAGCACCGGCCACGCGCACCGTGGAGATCATCAACGTGGAG
GACGAGCTGCCCGCCGCTCGGCCGTGTCCGTGGCCGCCGCGGCCGCGCTGGCCAGTGGGGTGGTGTG
GGCGGGGACAGCCACCTGGCCCTGCCCGCCCTGGAGCGAGACCACCTCAACCACCACCTACGTGGCT
GCCGCCTTCAAGGCGCACTACAGCAGCAACCCAGCGCGGGGGCTGCGGGGGCAAGGCCCGCCTGGC
CTCAACTCCATCCACGAACCTCTGCTCTTCAAGAGCGGCTCAAGGAGAAGCTGCAAGAGACGCAGATC
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001080457

Insert Size: 2142 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001080457.1
RefSeq Size:	2958 bp
RefSeq ORF:	2142 bp
Locus ID:	94030
UniProt ID:	Q9NT99
Cytogenetics:	19q13.33
Protein Families:	Druggable Genome, Transmembrane
MW:	76.4 kDa
Gene Summary:	Synaptic adhesion protein. Regulates the formation of excitatory synapses. The trans-synaptic adhesion between LRRC4B and PTPRF regulates the formation of excitatory synapses in a bidirectional manner (By similarity).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 both encode the same protein. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.