

Product datasheet for **SC315832**

CACTIN (NM_001080543) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CACTIN (NM_001080543) Human Untagged Clone
Tag:	Tag Free
Symbol:	CACTIN
Synonyms:	C19orf29; fSAPc; NY-REN-24
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001080543, the custom clone sequence may differ by one or more nucleotides

```

ATGGGTCGGGACACACGCTCGCGCTCGCGGTCCGCGGGTCCGCGGGGCCGAAGGCGGCAG
AGTCAGAGCGGGAGCCGAAGTCCGAGCAGGAGCCATGGGCGGGCAAACCGACGGCGCCGG
GAGGACGAGGGACGGCGCAGACGGAGGCGGGAGCCGGGAGCGCAGGTCAAGTTCAGAG
GAAGAGCGGTGGCAGCGCTCAGGGATGCGAAGCCGGAGCCCCCGCGGCCAAGTGGCAC
TCAAGAGATGGGTCTCTCAGTCGGACTCAGGAGAGGAGCAGTACGGGGCCAGTGGGCT
CGCCGGCGACGGCGGCACGCTCGTGGTCTCCTAGCTCCTCAGCATCCAGTCCGCGTCT
CCAGGGCGATCCCAGAGCCCCGGGCGGCGGCTGCCCTGAGCCAGCAGCAGAGCCTG
CAGGAGCGGCTGCGGCTGCGGGAGGAGCGGAAGCAGCAGGAGGAGCTGATGAAGGCCCTC
GAGACGCCGAGGAGAAGCGGCACGGCGGCTGGCCAAGAAGGAGGCCAAGGAGCGCAAG
AAGCGGGAGAAGATGGGCTGGGGTGAAGGAGTACATGGGTACACCAACACCGACAACCC
TTCGAGACAACAACCTGCTGGGCACCTTCATCTGGAATAAGGCCCTGGAGAAGAAGGGG
ATCAGCCACCTGGAGGAGAAGGAGCTGAAGGAGCGGAACAAGAGGATCCAGGAGGACAAC
CGGCTGGAGCTGCAAGAAGGTGAAGCAGCTGCGGCTGGAGCGGGAGCGGGAGAAGGCCATG
CGCGAGCAGGAGCTGGAGATGCTGCAGCGCAGAGAAGGAGGCAGAGCACTTCAAGACATGG
GAGGAGCAGGAGGACAACCTCCACCTCCAGCAGGCCAAGTGCCTTCAAGATCCGCATC
CGGGACGGGCGGGCCAAGCCCATCGACCTGCTGGCCAAGTACATCAGCGCTGAGGATGAC
GATCTGGCCGTGGAGATGCATGAGCCCTACACGTTCTCAACGGCCTCACCGTGGCCGAC
ATGGAGGACCTGCTGGAGGATATCCAGGTCTACATGGAGCTGGAGCAGGGCAAGAAGCC
GACTTCTGGCGGGACATGACCACCATCACCGAGGACGAGATCTCCAAGCTCCGCAAGCTG
GAGGCCCTGGGCAAGGGGCCAGGTGAGCGCCGCGAGGGGGTCAACGCCTCCGTCAAGTCT
GATGTGCAGTCGGTGTTCAGGGGAAGACATAACAACAGCTGCAGTGCATCTTCCAGGGC
ATCGAGGGCAAAATCCGCGCTGGTGGCCCAACCTGGACATGGGCTACTGGGAGAGCCTC
CTGCAGCAGCTTCGTGCCACATGGCACGGGCCCGGCTGCGTGAGCGCCACCAGGACGTG
CTGCGGCAGAAGCTGTACAACTGAAGCAGGAGCAGGGCGTGGAGAGCGAGCCGCTGTTC
CCCATCTCAAGCAGGAGCCCCAGTCCCCAGCCGACGCTGGAGCCTGAGGACGCGGGC
CCCACCCCGCCGGGCCCTCCTCGGAGGGCGGCCCGCGGAGGCCGAGGTGGACGGCGCG
ACCCCGACAGAGGGGACGGCGACGGGGACGGTGAGGGCGAGGGCGAGGGCGAGGCGGTG
CTCATGGAGGAGGACCTGATCCAGCAGAGCCTGGACGACTACGACGCCGGCAGGTACAGC
CCGCGGCTGCTCACGGCGCACGAGCTGCCACTGGACGCGCACGTGCTGGAACCGGATGAG
GACCTGCAGCGCCTGCAGCTCTCGCGCCAGCAGCTCCAGGTACGGGAGACGCCAGCGAG
AGCGCCGAGGACATCTTCTCCGGCGGGCCAAGGAGGGCATGGGCCAGGACGAGGCGCAG
TTCAGCGTGGAGATGCCACTCACCGGAAGGCCTACCTGTGGGCCGACAAGTACCGGCCA
CGCAAGCCGCGCTTCTTCAACCGCGTGCACACGGGCTTCGAGTGGAAACAAGTACAACAG
ACGCACTACGACTTTGACAACCCACCGCCAAGATCGTGCAGGGATACAAGTTCAACATC
TTTACCCCGACCTCATCGACAAGCGCTCCACGCCCGAGTACTTCTGGAGGCTGCGCC
GACAACAAGGATTCGCCATCTGCGCTTCCACGCGGGGCCCTACGAGGACATCGCT
TTCAAGATCGTCAACCGCGAGTGGGAATACTCGCACCGCCACGGCTTCCGCTGCCAGTTT
GCCAACGGCATCTCCAGCTGTGGTTTCACTTCAAGCGCTACCGCTATCGGCGG

```

Restriction Sites: Please inquire

ACCN: NM_001080543

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080543.1</u> , <u>NP_001074012.1</u>
RefSeq Size:	3625 bp
RefSeq ORF:	2277 bp
Locus ID:	58509
UniProt ID:	<u>Q8WUQ7</u>
Cytogenetics:	19p13.3
Gene Summary:	<p>Involved in the regulation of innate immune response. Acts as negative regulator of Toll-like receptor and interferon-regulatory factor (IRF) signaling pathways. Contributes to the regulation of transcriptional activation of NF-kappa-B target genes in response to endogenous proinflammatory stimuli. May play a role during early embryonic development. Probably involved in pre-mRNA splicing.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript. Variants 1 and 2 encode the same protein.</p>