

Product datasheet for **SC315834**

KIAA1751 (CFAP74) (NM_001080484) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIAA1751 (CFAP74) (NM_001080484) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIAA1751
Synonyms:	C1orf222; KIAA1751
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC315834 representing NM_001080484.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAGGATGACGGCAGCCTGCTCCCTGAGGACGAGCTTTGGCCGATGCCCTTCTTTTGAAGATGAA
AGAGATGAATTAGAGGATCCGGAGTTTGACATCAAATGTCTTCTACAGGAAGCTGAGGATGACGTGGAC
CCGGGACACAGCAGCTCAGTGAAAGAACTAGATACTGATGCTGATAAATTGAAGAAGAAAAACAGCTGAG
GACAGAACGCAGGCATTTACCTGCGGCAGAACCTGAGCGCCCTGGATAAGATGCATGAGGAGCAAGAG
CTTTTCACTGAGAAGATGAGAGGGGAGCTGCGCGCTGTGCGCAGAGGCGGGACCTGATCGACAAGCAG
CAGGAGGCTGTGGCAGCCGAGATCGCCACAGAGGAAGAGCGGGCAACATGGCCGCTGTGGGCCGCTC
CAGGCCGTGTCCAGACGCTGTTTGCAGAGCTGGAGAACGAGAGAGACCTGCAGTCGAGGACTGAGGCC
GTGCTGAAGGAGAGCGAGAACCACATGTGCACATCGAGATCCAGGAGGGGCGACTCGAGGCCTCCGG
ACAGCTGACCGTGAGGAGGTGGAGGCCACGGGGCGCGGCTCCAGGTGCGCGCAGCCGAGCAGCTCTGC
AGAGAGCAGGAGGCCCTGGGAAGGTGGAGCGGAACCGACTGCTGAGGATCAGGAAGTCCCTGAACACC
CAGAAGGAGCTCGGGCTCAGGCACCAGAAGCTGCTGGAGGACGCCCGGAAGAACCAAGGTTGCCGTG
CGGTTCTGAAGGCTCCCTGGGAAGAATCCGAGAGCAAGAGAAGAGGAGGATGGAGTGCCACGAG
TACATGAGGCGACGCATGGATGCGGTGGTGGCGCTGAAGGGCAGCATCTCCGCGAACCGGGACACACTG
CGGAAGTTCCAAGCATGGGACCGTGCCAAGGCAGAGCTGGCGGAGCAGAGGTTCCAGGCTGAGAAGAAG
GCGATTCTGGCCAGGGCAGGGATGCATTACGGCACCTTGTCCACCAGCGCCGGCCAGGAGCTGGAG
GCCCAGAAGAGGGCCTTTGAGGAGGAGCAGAAGCTCAGAAAGCAGGAGATCATCAGTCGGATTCTGAAA
GAGGAGGCTGAGGAGGAAAAGAGGAAGAAACAGCATCCCCCACCAGTGCCAGGCACCGGCTGACCTG
AGGGACAAGACCTGGAACCTACATTTCTGACTTTTGCAAGAAGACCACAGTCCCAACCAACAGTACACA
CTGGACTACGAGGCTGCTGCAGGCCCGGCCCTCTCGGTTGCTGGAAGTCGTTTCCAGTGAGCTTATC
CAGGGGGACCCCGGGCCAGCTCAGAGGAGGAAACGTTAGCTGAGCCCGAGATCTCTGGGCTTTGGAAT
GAAGACTACAAGCCATACCAGGTGCCAAGGAGGACGTGGACCGAAAGCCGTGGGCGGGACAAAGATG
GACAAGGACATCCTGGAGCGCACGGTGGAGCGGCTGCGGAGCAGGTTGGTCCACAAGCAGGTGGTGTGG
GGGCGTGAGTTCCAAGGACGCCCTTCAACAGCAAACCGGAGCTCCTCCACTTCCAGGACTTTGATATT
GGCAAAGTGTAAGAAAAAGATCACGTTGGTAAACACCACCTACACGATCAACTACTGCAAGCTGGTG
GGCGTGGAGGAGCACCTCCGGGACTTCATCCACGTTGACTTTGACCCCTGGCCCCCTGTCAGCCGGA
ATGTCCTGTGAAGTGCTGTACCTTCAAGCCGATGATAAACAAGGATCTAGAAGGAAATATCTCATTT
TTGGCTCAGACGGCGAGTTTTAGTTCCACTGAAATGTTCAACAAAGAAATGTTGCTGTCCCTCGAC
AAGGAGCTCATTGACTTCGGCAGCTACGTGGTAGGAGAGACCAGTCTCGGACCATCACGCTGACCAAC
GTTGGGGGCTTGGGCACGACTTTCAAGTTCTGCCAGCTTCAGAGCCCTGTGAGATGGACGACTCCAG
TCTGCCCTGAAATTAAGCAGTCTCCTGACCTACGAAGATAAAGCTTGTATGACAAAGCCGCCACCAGC
TTCTCTGAGCAGCAGCTAGAGGGCAGGAGTCTCCAGGCGGACATGCAGAGCCGAAGGAGCTGGAG
AAGCTGGACAAGGAGCAGGAGGAGGAGCAGCCGAGGTGAGCCAGGCTGCCAGGCCAAGCCAGGTTGG
GCGGGCGGTGAGGGCTCAGGGCAGATTCCAGGGGCTGAGGGGCCAGGGGAGGGCTGGAAGTAGA
GAATTTTCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001080484

Insert Size: 2289 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080484.1</u>
RefSeq Size:	4710 bp
RefSeq ORF:	2289 bp
Locus ID:	85452
Cytogenetics:	1p36.33
MW:	87 kDa
Gene Summary:	As part of the central apparatus of the cilium axoneme may play a role in cilium movement. [UniProtKB/Swiss-Prot Function]