

## Product datasheet for **SC315835**

### SLC45A4 (NM\_001080431) Human Untagged Clone

#### Product data:

|                           |                                             |
|---------------------------|---------------------------------------------|
| Product Type:             | Expression Plasmids                         |
| Product Name:             | SLC45A4 (NM_001080431) Human Untagged Clone |
| Tag:                      | Tag Free                                    |
| Symbol:                   | SLC45A4                                     |
| Mammalian Cell Selection: | Neomycin                                    |
| Vector:                   | pCMV6-Entry (PS100001)                      |
| E. coli Selection:        | Kanamycin (25 ug/mL)                        |



[View online »](#)

**Fully Sequenced ORF:** >SC315835 representing NM\_001080431.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGGAAGCAAGCCCAGCCTCAGGGCTGTCCCGTCCAAAGACCCTTGAGTCCACTGAGGTCAAAT
CAATGGAGTCTTCAGAAAGGCCTTCGGAGCAGTACTACAGCCTCACCTGGTTCCTGAGCCCCATCCTT
GGCCTCATCTTCACACCTCTCATTGGGTCTGCGAGTGACCGGTGCACCCTGAGCTGGGGCCGCGGCGG
CCCTTCATCCTCGCCCTCTGCGTTGGCGTCTCTTTGGCGTTGCACTTTTCTTAACGGCTCTGCCATC
GGTCTGGCCCTCGGCGATGTCCCAACCGGCAGCCCATTTGGCATCGTGCTCACGGTGCTGGGAGTGGTG
GTCCTGGACTTCAGCGCCGATGCCACCGAGGGGCCATCCGTGCCTATCTGCTGGACGTGGTGGACAGC
GAGGAGCAGGACATGGCCCTCAACATCCACGCCTTCTCTGCCGGCCTCGGCGGAGCCATCGGCTACGTG
CTGGGTGGGTGGACTGGACCCAGACCTTCTGGGCAGCTGGTTCGGACCCAGAACCAGGTGCTCTTC
TTCTTTGCCGCCATCATTTACGGTGTCCGTGGCCCTGCACCTGTTGAGCATCGACGAGGAGCAGTAC
AGCCCGCAGCAGGAGCGCAGCGTGAGGAGCCCGCGCCCTGGATGGGGCGAGCCGACGGCGTCCCT
GCCTTCCAGACGAGGTACAGTCGAGACGAGCTGGCCCTGGACTACCCGGACGTGGACATCATGCGC
AGCAAAAGCGACTCGGCATTGCACGTGCCGACACCGCGCTGGACCTGGAGCCCGAGCTGCTGTTCTG
CACGACATCGAGCCCTCCATCTTCCACGACGCCTCTACCCCGCCACCCCGCAGCACCAGCCAGGAG
CTCGCCAAGACCAAGCTGCCCGCCTGGCCACCTTCTCAAGGAAGCCGCAAGGAGGACGAGACCTTG
CTGGATAATCACTTGAATGAAGCTAAAGTCCCAACGGAAGTGGTCCCCCAAAAGACGCCCTCGGC
GGCTACACCAGGTGGACACGAAGCCCTCGGCCACGTGAGCTCCATGCGGGCGGGCGGCACGCGTTC
CGCAGGACGGCCTCCAGACCTTCTCTACTACGGCAAGCTTGGTCCCCTGCTACCGCTACCGGCGC
GCCAACGCCGTGGTCTGATCAAGCCGTGCGCAGCATGAGCGACCTGTACGACATGCAGAAGCGGCAG
CGGCAGCACCCGGCACCGGAACCGAGAGCGGGGCCACCACTCCAGCGGGGACACCGAGAGTGAGGAGGG
GAGGGCGAGACCAGGTGCGCCTGCTGTGGCTCTCCATGCTGAAGATGCCAGGGAGCTGATGCGGCTG
TGCTCTGCCACCTCTCACCTGGTCTCTGTCTCGCCGAGGCCGTGTTCTACACCGACTTCATGGGC
CAGGTCTCTTGAAGGCGACCCCAAGGCCCTCGAACTCGACCGCCTGGCAAGCCTACAACGCCGGG
GTCAAGATGGGTGCTGGGGCTGGTCAATTTATGCCGCCACTGGTGTATTTGTTAGCCCTGTTACAG
AAGTACTTGGACAACACGACCTGAGCGTCAGGGTGATCTACGTGCTGGGACGCTGGGCTTCTGTG
GGCACAGCCGTGATGGCCATGTTTCCCAACGTCTACGTGCGCATGGTACCATCAGCACCATGGGCATC
GTCTCCATGAGCATCTCTACTGCCGTACGCCCTGCTGGGCCAGTACCATGACATCAAGCAGTACATC
CACCACAGCCCCGGAACTCCAAGCGAGGGTTTGGCATAGATTGTGCCATCCTGTCTGCCAAGTGATC
ATCTCGCAGATCCTGGTGGCCTCTGCCCTTGGGGCGTGGTGCAGCGCGTGGGACTGTCCGCGTCATC
CCCATGGTGGCCTCTGTGGGCTCTTTCTGGGCTTCTGACGGCCACATTCTGGTGATCTATCCCAAC
GTGTCAGAGGAGGCCAAGGAGGAGCAGAAAGGCCTGTCTTCCCCGTTGGCCGGCAAGGCAGGGCCGGT
GGGAACAGCGAAAAGCCACCGTGTGAAGTCAACGCGGAAGGAGGGCTGCAGGGACCGGTGGAGACA
GAACGCTGCAGGTATTGACGTCTGTGAGATCTCGTCACATTGGCTGGTGCCGACGTGTTGGAGAGTA
TTTTCTTTATGATTATTTAGAAAAAAATTTCTTTTCCACAATGTGGTCTCTTAGAAGAATGACG
TATCTTCTTTCTCAGCGAGTTGGACACATTGTGCCAGGGCAGCCCTGTCTTGGGCAGCGACCGCA
CACCAAAGCTGGGAGGAGGCTGGTCCGGGGGCTGGGCAGAAGACAGTGA
AGCGGACCGACGCGTACGCGGCCCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

**Restriction Sites:** SgfI-RsrII

**ACCN:** NM\_001080431

**Insert Size:** 2397 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Annotation:</b>        | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.                                                                                                                                                                                                                                                      |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_001080431.2</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>RefSeq Size:</b>           | 6991 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 2397 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 57210                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>UniProt ID:</b>            | <u><a href="#">Q5BKX6</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Cytogenetics:</b>          | 8q24.3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Protein Families:</b>      | Transmembrane                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>MW:</b>                    | 87.6 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |