

Product datasheet for **SC315841**

USP6NL (NM_001080491) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	USP6NL (NM_001080491) Human Untagged Clone
Tag:	Tag Free
Symbol:	USP6NL
Synonyms:	RNTRE; TRE2NL; USP6NL-IT1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC315841 representing NM_001080491.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAGGACTGAAGCTGTTGTTAATTCCTCAAGGCCAGACTGATTATCTCACAAGGATTAGACCAGGAT
GTAGCACTCAAACTTGCCAGGAGCGAGCTGAAATAGTTGCTAAATATGACAGAGGACGAGAAGGTGCA
GAGATTGAACCTTGGAAGATGCTGATTACCTTGTTCACAAAGTCACAGATAGATTGGCTTTTACAT
GAGGAGGAGCTCCAGATCATAATGTGGCTGTGGAACGGCAAAAGCACCTGGAAATTGAAAGAACTACC
AAATGGCTGAAAATGCTGAAAGGATGGGAAAAATACAAGAACTGAAAAGTTTCATAGGCGAATTTAC
AAAGGAATACCACTCCAGCTCAGAGGTGAAGTCTGGGCCCTCTTCTTGAGATCCCTAAATGAAAGAA
GAAACAAGGGACCTGTATAGTAAATTAACACAGAGCACGGGGCTGTTACCTGACATCAGACAAATA
GACCTGGATGTCAACCGCACATTCGGGACCACATTATGTTTAGAGACAGATATGGTGTTAAGCAACAA
TCCTTATTCCATGTGCTTGTCTGCTATTCTATTTATAACACGGAAGTCGGGTATTGTGAGGGGATGAGC
CAGATCACAGCTTACTCCTCATGTATATGAACGAGGAAGATGCCTTCTGGGCCCTGGTCAAACCTCTTC
TCAGGCCCTAAACATGCCATGCATGGCTTTTTTGTCCAAGGTTTTCTAAACTCTTGAGGTTTCAAGAA
CATCATGAAAAATACTGAACAAATTTCTGTCCAAGCTTAAGCAACACTTGGATTCTCAAGAAATCTAC
ACAAGTTTTTACACAATGAAATGGTTTTTTCAGTGTTCCTTGATCGTACTCCCTTACACTAAACCTC
AGAATATGGGATATCTACATCTTTGAAGGAGAAGGAGTCTTACTGCTATGTCTTACACCATCTTAAAA
TTACACAAAAACATCTAATGAAATGTCCATGGAAGAACTTGTAGAATTTTTTCAGGAGACCCTGGCA
AAGGATTTTTTCTTTGAAGATGATTTTGTGATAGAGCAACTTCAGATTTCTATGACAGAACTAAAGCGG
GCAAGTTAGACCTTCCAGAACCTGGTAAAGAGGATGAATATCCAAAGAAGCCCTTGGGGCAGCTTCCA
CCTGAACCTTCAGTCTTGGGGCGTCCATCACTTGAGCAACGGACAGAGGAGCGTGGGCCGGCCGAGCCCG
CTGGCCAGCGGCAGGAGGAGAGCGGGCGCCACAGGAGGCACGAGCACTCCCGCACCCCCAGAGC
AGGACCGGGACGCCGAGAGAGCACAGCCGCCAAGACGGAATCGGTGGAGGAGGAGCAAAAAGCTT
AAAGATGAGGCAGATTTTCAAAGAAAACCTCCATCGGGTCCACAGGACAGTTCAGGCAATATAATCAC
GCAGCTGCCAACCAAAATAGCAACGCCACTTCAAATATCAGGAAGGAGTTTGTGCCCAATGGAATAAA
CCGTGACAGCTCTCAGCTACAGAGAGAACTGCCAAATACACCATGGAAGGCAAGGTGAGCAGCGCAC
CCCGCGCTCGCAGTTACCGTCCCAGGTCTGCCGAGGTGCGGGTGTCAAACGTGCGGCCAAAGATGAAG
GCCCTGGATGCTGAGGACGGGAAGCGGGCTCCACTGCATCGCAGTACGACAACGTGCCAGGCCCGGAG
CTGGACAGCGGCGCTTCCGTGGAGGAGGCGCTGGAAAGGGCTTACTCCCAGAGCCCCGGCATGCCCTT
TACCCTCCTAGCCCGAGAAAGCACGCTGAGCCAAGTTCTAGTCCATCAAAGATATCCAACAAGTTTACT
TTTAAAGTACAGCCTCCAAGTCATGCACGATATCCGTCCCAGCTAGATGGGGAAGCCCGAGGGCTAGCT
CATCCCCCTCTACAGCAATCCCCCGTTTACCACGGAACCTCTCCCAACACTTCCCTACTGCCAAC
AGCAGCTTTGCTTCTCCACAGTTTAGCCCTGGGACTCACTGAATCCTTCCAGGAGACCTCATGGTTCT
ACTCTTCCGTGAGTGCTTCTCCGAGAAATCTTACAGCCGCCAAGCCCTTGTACTGCCGTCTAGT
CGAATAGAAGTCTCCTGTTGACACTGGTGCTGGGGGATTCGGGCAATTCAGGTCACCAAAGAAT
GAAAAATTGATCATTCCACAGTGGATTACTTGCCAGATAACAGAACATGGTCAGAAGTTAGTTATACA
TACAGACCTGAGACGCAGGACAAATCATGGACCGAGATGCTAGCCGTGGCAATTTACAAAATACTCA
GCCTTCAACTCGACCCCTTTCAGGACCATGGCTCCCTGCAGTTTCTGTAGATAGTCCCGTGAATAT
AAAGCTTACCGGCCGAGAAGATGCCAGTCCATCTGGATATCCATATTCAGGGCCCCCGCTCCAGCC
TACCACTACAGGAATCGGGACGGGCTTTCCATCCAAGAGTCAGTGTTGCTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001080491
Insert Size: 2538 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080491.3</u>
RefSeq Size:	4517 bp
RefSeq ORF:	2538 bp
Locus ID:	9712
UniProt ID:	<u>Q92738</u>
Cytogenetics:	10p14
Protein Families:	Druggable Genome
MW:	96 kDa
Gene Summary:	Acts as a GTPase-activating protein for RAB5A and RAB43. Involved in receptor trafficking. In complex with EPS8 inhibits internalization of EGFR. Involved in retrograde transport from the endocytic pathway to the Golgi apparatus. Involved in the transport of Shiga toxin from early and recycling endosomes to the trans-Golgi network. Required for structural integrity of the Golgi complex.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (2) has a longer and distinct N-terminus compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.