

Product datasheet for **SC315842**

KIF18B (NM_001080443) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIF18B (NM_001080443) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIF18B
Synonyms:	kinesin family member 18B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC315842 representing NM_001080443.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCAGTGGAGGACAGCACGCTGCAAGTAGTGGTACGGGTGCGGCCCCACCCTCGGGAGCTGGAC
AGTCAGCGCGCGCCAGTGGTTCAGGTGGTGGACGAGCGGGTGTGGTGTTTAACCTGAGGAGCCGAT
GGAGGGTCCCTGGCCTGAAATGGGGTGGCACCCATGATGGCCCCAAGAAGAAGGGCAAAGACCTGACG
TTTGTCTTTGACCGGGTCTTTGGCGAGGCGGCCACCAACAGGACGTGTTCCAGCACACCCACGCACAGC
GTCCTGGACAGCTTCTCCAGGGCTACAACGCTCAGTGTTCGCTACGGGGCCACCGGGGCTGGGAAG
ACACACACCATGCTGGGAAGGGAGGGGACCCCGGCATCATGTACCTGACCACCGTGGAAGTGTACAGG
CGCCTGGAGGCCCGCCAGCAGGAGAAGCACTTCGAGGTGCTCATCAGCTACCAGGAGGTGTATAATGAA
CAGATCCATGACCTCCTGGAGCCCAAGGGGCCCTTGCCATCCGCGAGGACCCCGACAAGGGGGTGGTG
GTGCAAGGACTTTCTTTCCACCAGCCAGCCTCAGCCGAGCAGCTGCTGGAGATACTGACCAGGGGAAC
CGTAACCGCACGCAGCACCCACTGATGCCAACCGGACTTCTCCGCTCCCATGCCATCTCCAGATC
TTTGTGAAGCAGCAGGACCGGGTCCAGGACTGACCCAGGCTGTCCAGGTGGCCAAGATGAGCCTGATT
GACCTGGCTGGCTCAGAGCGGGCATCCAGCACCCATGCGAAGGGGGAGCGGCTGCGGGAGGGGGCCAAC
ATCAACCGCTCTCTGCTGGCGCTCATCAACGCTCCTCAATGCCCTTGCCGATGCAAAGGGCCGAAGACC
CATGTGCCCTACCGGGACAGCAAAGTACCCGCTGCTCAAAGACTCCCTCGGGGCAACTGCCGCACA
GTGATGATCGCTGCCATCAGCCCTCCAGCCTGACCTACGAGGACACGTACAACACCCCTCAAATATGCC
GACCGGGCAAGGAGATCAGGCTCTCGCTGAAGAGCAATGTGACCAGCCTGGACTGTCACATCAGCCAG
TATGCTACCATCTGCCAACAGCTCCAGGCTGAGGTAGCCGCTCTGAGGAAGAAGTCCAAGTGTATGAG
GGGGGAGGCCAGCCCCACCACAGGACCTCCAGGATCTCCCAAGTCCGGACACCACAGAACCCCTT
CCAGTCCCCCTTGCCACCCACCCCTCCAGCCAGCCCTGCACCCAGAGCTCCCTGCAGGGCCTAGA
GCCCTTCAAGAGGAGAGTCTGGGGATGGAGGCCAGGTGGAGAGGGCCATGGAAGGGAACCTTTCAGAC
CAGGAGCAGTCCCAGAGGATGAGGATGAAGGCCAGCTGAGGAGGTTCCAACCCAGATGCCAGAGCAG
AACCCACACATGCACTGCCAGAGTCCCCTCGCCTGACCTGCAGCCAAAGCCAGTCGTGGGCCACTTC
TCAGCACGGGAAGTGGATGGGGACCGTTCTAAGCAGTTGGCCCTAAAGGTGCTGTGCGTTGCCAGCGG
CAGTACTCCCTGCTCCAAGCAGCAACCTCCTGACGCCGACATGATCACAGAGTTTGGACCCTACAG
CAGCTGGTGAAGAGGAAAAATGAGCCTGGGGCAGAGGCTTGAGGACTTCAGGCCTGGCCAGGGGG
GCACCTCTGGCTCAGGAGCTGTGTTTCAAGTCAAAGCCTCCAGGATACACTGGCCCTGTGACCCGGACT
ATGGCGAGGGCAGTGAAGTGGCCCCCTGCACACCCCTGGGAATCCCGCTGGACCCAACTGCACCCAGCC
CAGGGGTCCCGATGGCCCATGGAGAAGAAGAGGAGGAGACCAAGCGCCTTGAGGGCAGACAGTCCCATG
GCCCAAAGCGGGGACCAAGCGCCAGCGCCAGTCCCTCCTGCCCTGCCTAAGGAGAGGGTCTCTGCCT
GACACCCAACTTCACAGGGGCCAGCACCCCAAGGAGAAAGGGCCTCCTCCCTGCCATTCCTCCCT
CGCGTTTGCCAGCCACAGTCAAAAAGCCGGGTGCCCTGGGCCCTCCGCCATGCAGAACTGCTCC
ACCCCGCTGGCTCGCCACTCGAGACCTCAATGCCACCTTTGATCTCTCTGAGGAGCCTCCCTCAAAG
CCCAGTTTCCATGAATGCATTGGCTGGGACAAAATACCCAGGAGCTGAGCAGGCTGGACCAGCCCTTC
ATCCCCAGGGCACCTGTGCCCTGTTCCACATGAAGGGCCCAAGCCAACATCTTCCCTCCCTGGGACC
TCTGCCTGCAAGAAGAAGCGCGTTGCGAGTTCCTCAGTCTCCATGGCCGAGCCGCATCGCCCGCCTC
CCCAGCAGCACTTTGAAGAGGCCAGCTGGGCCCTTGTACTCCAGAGCTGCCCTTGAGTCCCCTGTGC
CCTAGCAACCGGAGGAATGGAAAGGACCTCATCAGGTTGGGAGAGCACTCTCAGCAGGGAACGGCGTC
ACCAAGGTGTCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001080443

Insert Size: 2568 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001080443.1](#)

RefSeq Size: 2669 bp

RefSeq ORF: 2568 bp

Locus ID: 146909

Cytogenetics: 17q21.31

Protein Families: Druggable Genome

MW: 93.3 kDa

Gene Summary: In complex with KIF2C, constitutes the major microtubule plus-end depolymerizing activity in mitotic cells. Its major role may be to transport KIF2C and/or MAPRE1 along microtubules. [UniProtKB/Swiss-Prot Function]