


```

CAGCTGCCCTGTCCGAAGGGGCAGTTTGGAGAAGGTTGTGCCAGTCGCTGTGACTGTGAC
CACTCTGATGGCTGTGACCCTGTTTCATGGACGCTGTAGTGCCAGGCTGGCTGGATGGGT
GCCCCGCTGCCACCTGTCTGCCCTGAGGGCTTATGGGGAGTCAACTGTAGCAACACCTGC
ACCTGCAAGAATGGGGGCACCTGTCTCCCTGAGAATGGCAACTGCGTGTGTGCACCCGGA
TTCCGGGGCCCCCTCCTGCCAGAGATCCTGTGAGCCTGGCCGCTATGGCAAACGCTGTGTG
CCCTGCAAGTGCCTAACCCTCCTTCTGCCACCCCTCGAACGGGACCTGCTACTGCCTG
GCTGGCTGGACAGGCCCCGACTGCTCCAGCCATGCCCTCCAGGACACTGGGGAGAAAAC
TGTGCCCAGACCTGCCAATGTACCATGGTGGGACCTGCCATCCCCAGGATGGGAGCTGT
ATCTGCCCCCTAGGCTGGACTGGACACCACTGCTTAGAAGGCTGCCCTCTGGGGACATT
GGTGCTAACTGCTCCCAGCCATGCCAGTGTGGTCTGGAGAAAAGTGCCACCCAGAGACT
GGGGCCTGTGTATGTCCCCAGGGCACAGTGGTGACCTTGACAGATTGGAATCCAGGAG
CCCTTTACTGTGATGCCGACCACTCCAGTAGCGTATAACTCGCTGGGTGCAGTGATTGGC
ATTGCAGTGCTGGGGTCCCTTGTGGTAGCCCTGGTGGCACTGTTTCATTGGCTATCGGCAC
TGGCAAAAAGGCAAGGAGCACCACCACCTGGCTGTGGCTTACAGCAGCGGGCGCCTGGAC
GGCTCCGAGTATGTCATGCCAGATGTCCCTCCGAGCTACAGTCACTACTACTCCAACCCC
AGCTACCACACCCGTGTCGAGTGTCCCAAACCCCCACCCCTAACAAAGGTTCCAGGC
CCGCTCTTTGCCAGCCTGCAGAACCCTGAGCGGCCAGGTGGGGCCCAAGGCGATGATAAC
CACACCACCTGCCTGTGACTGGAAGCACCGCCGGGAGCCCCCTCCAGGGCCTCTGGAC
AGGGGGAGCAGCCGCTGGACCGAAGCTACAGCTATAGCTACAGCAATGGCCCAGGCCCA
TTCTACAATAAAGGGCTCATCTCTGAAGAGGAGCTCGGGGCCAGTGTGGCTTCCCTGAGC
AGTGAGAACCCATATGCCACCATCCGGGACCTGCCAGCTTGCCAGGGGGCCCCGGGAG
AGCAGCTACATGGAGATGAAAGGCCCTCCCTCAGGATCTCCCCCAGGCAGCCTCCTCAG
TTCTGGGACAGCCAGAGCGGGCGGCAACCCAGCCACAGAGAGACAGTGGCACCTACGAG
CAGCCCAGCCCCCTGATCCATGACCGAGACTCTGTGGGCTCCCAAGCCCCCTCTGCCTCCG
GGCCTACCCCCCGGCACTATGACTCACCAAGAAGCAGCCACATCCCTGGACATTATGAC
TTGCTCCAGTACGGCATCCCCATCACCTCCACTTCGACGCCAGGACCGTTGA

```

Restriction Sites:

Please inquire

ACCN:

NM_001080471

Insert Size:

3200 bp

OTI Disclaimer:

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080471.1, NP_001073940.1</u>
RefSeq Size:	4866 bp
RefSeq ORF:	3114 bp
Locus ID:	375033
UniProt ID:	<u>Q5VY43</u>
Cytogenetics:	1q23.1
Protein Families:	Transmembrane
Gene Summary:	PEAR1 is a platelet receptor that signals upon the formation of platelet-platelet contacts independent of platelet activation and secondary to platelet aggregation (Nanda et al., 2005 [PubMed 15851471]).[supplied by OMIM, Mar 2008]