

Product datasheet for **SC315878**

KIAA1543 (CAMSAP3) (NM_001080429) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIAA1543 (CAMSAP3) (NM_001080429) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIAA1543
Synonyms:	KIAA1543; NEZHA; PPP1R80
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC315878 representing NM_001080429. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTGGAGGCGGCGCCCCCGGGCCGGGCGCTGCGGAGGACCTTTCTAGTGCCCGAGATCAAGTCG
CTGGACCACTACGATTTCTCGGGGCAAGGCGGCGGCCAGCCTGGCGTGGTGTGCGGGCCGCGTTT
GGGGGCGCAGAGCAGTGCCTCCCGGAGCTGTGGGAGCCCTTCTATACCGACCACTACGCGCAGGAGCAT
GTGAAGCCCCCGGTGACACGGCTGCTGCTCTCAGCCGAGCTCTACTGCAGAGCCTGGCGCCAGGCACTG
CCACAGCTTGAACACCCCCAACCCCTCTGCACTGCTGGCCCTGCTGGCGGGAGGGGCACAGTGCCT
GCTTTGCCGAGCGCCCGGTGCGCGAGGCGCAGCTGAGGCACCAGCCATTCTCATGGGAGCCACCTA
GCTGTCAATTGATGCCCTCATGGCTGCCTTTGCCTTCGAGTGGACAAGACCCTGCCAGGTCCTTGCC
CTGACCAGCTTGAGCACAAGCTGCTTTTCTGGGTGGACACGACCGTCCGGCGGCTGCAGGAGAAGACC
GAGCAGGAAGCGGCCAGCGAGCCTCTCCAGCAGCCCTGCAGACGGGGCGGCCCCGGCGCAGCCCTCG
TGCCCTACGCGCTGGTACTGGAAGCTGGTTCCTCACGCAATTGCCTTCTGTTTGAAGGAGTCGGGGAGC
AAACCCCATGATCCGATACCGCAAGGACCGTGTGGTGGCGGACGCTGCCCTGCTTCCCGACGGTG
ACCAGCTCCAGGACCTGGCCAGTGGGGCCGCGCTGGCCGCCACCATCCACTGCTATTGTCCCAGCTG
CTTCGACTTGAGGAGGTGTGCTTGAAGGACCCATGTCTGTGGCGGACAGCCTGTACAACCTCCAGCTC
GTGCAGGATTTCTGTGCCTCTCGCCTTCCCTGCTGGCTGCCCTGTCCCTTGAGGACTTGTGTACGTC
CCACCGCACTCAAGGTCAACTTGGTGGTGTGCTGGCCGAGTTGTTTGTGTTTGGAGTGTCAAG
CCCGACTTTGTGCAAGTGAAGGACTTGCCCGATGGTACGCTGCCTCCCCCGGGCACTGAGGCCTCC
CCACCTCAGAACAACAGCGGAGTGTCTCCTGTCTTACCTTCCGCCACCCGCTTCTGTATCTGGT
GGCCCCAGTCCCCTCCGAGGATCCACAGGCTCCCTGAAGTCTTCCCGTCCATGTCCCATATGGAG
GCCCTGGCAAGGCTGGAACCGGACGCTCAGCCGTCCCCTCTCCAGGCTGTGCTATTACGACCCCC
TTTGGCTGGACAGCGAGTGGATGTGCTCATGGGAGACCCTGTGCTCCTCCGCTCTGTGAGCTCGGAC
AGCCTGGGCCCCCGCTCCCGCGCCGGCCAGGACCCCAACCAGCCACCCCGGAGCCTGGTACCTG
CCCACCATCGAGGAAGCTCTGCAGATCATCCACAGTGCCGAGCCCGGCTCCTCCAGATGGGGCGGCC
```



[View online >](#)

GACGGCAGCTTCTACCTCCACTCCCCTGAGGGGCCCTCCAAGCCATCCCTGGCCTCCCCCTACCTGCCC
 GAGGGGACCTCCAAACCACTGTCCGACAGGGCCACCAAAGCACCAGTGTACATGCCACACCCGAGACC
 CCCTCGAAACCATCTCCCTGTCTGGTGGGGGAGGCATCGAAACGCCAGCCCATCCGAGGGTCCCCG
 AAGGGGTGGCTTCGTCGCCAGCAGCCACCACTCCGAGGTGAAAATGACCAGCTTTCGAGAACGCAAG
 AAACAGCTGGTGAAGGCAGAGGCTGAGGCCGGAGCGGGTCCCCACGTCCACTCCGGCCCCGCGGAG
 GCCCTGAGCTCGGAGATGAGTGAGCTCAGCGCCGGCTGGAGGAGAAACGCAGAGCCATCGAGGCTCAG
 AAGCGAGGATTGAGCCATATTCGCCAAGCACCGCCAGCGGCTGGGCAAAAGCGCCTTCTGCAGGTG
 CAGCCCGGGAAAGCCTCTGGGGAGGCGGAAGCAGAGGCGGAGGAGCCGATTCCGGTCCAGTCCCTGGT
 GGGGAGCGGCCCGCAGGCGAGGGCCAGGGTGAAGCAACCTCACGGCCAAAGGCAGTGACCTTCTCGCCA
 GACCTGGGCCCGGTGCCACAGAGGGGTGGGGGAATACAATCGAGCGGTGAGCAAGCTGAGTGCCGCC
 TTGAGCTCGTGCAGCGGACATGCAGAGGCTCACGGACCAGCAGCAGCGGCTCTGGCCCCGCGGAG
 GCCCCCGATCCGCCACACCTGCTGCGTGGTTCATCCCTGGCCCCACGACGGGGCCAAAGTGCA
 TCCCCAGCCCCCGCGGAGTCCCGCCACCCGGCGAGCCCTGGGCCGGGCCAGCCAGTCAACC
 CGCAGCCCGAAACACACGCGGCCAGCGGAGCTGCGGCTGGCACCTTGACCAGGGTGCTTACGCCACCC
 CACGACGTAGACAGCCTCCCCACTGCGCAAGTTCTCGCGAGCCAGGTGCCGTGCAGACGCGCTCT
 TCCATCTCTGGCGGAGGAGACGCCCCCGAGGAGCCAGCCGCCGGCCGGCCCTCATCGAGATCCCG
 CTGGGCAGCCTGGCAGATCCCGCCGCGAGGACGAGGGAGACGGGAGCCCGCTGGTGTGAGGATTCC
 TTGAGGAGGAGGCGTCTTCGGAGGGGAGCCCGGTGGGGTGGGGTCTTCTACAAGGATGAAGAC
 AAGCCTGAGGACGAGATGGCCAAAAGCGGGCCAGCCTGCTGGAGCGGCAGCAGCGGCGAGCAGAGGAG
 GCGCGGCGGCGCAAGCAGTGGCAGGAGTGGAGAAGGAACAGCGGAGGGAGGAGCCGCGAGGCTGGCC
 CAAGAGGAGGCCCGGGCCAGCCCCGCTTGTGTCCGAGTCCGATGGCGACTCCAGCCCTGCTGCC
 CGGGCTCCAGCCGAGGAGGAGTGGGCCCGGAAGGGGACTTACGCGGCAGGAGTACGAGCGCCGG
 GCCCAGTGAAGCTGATGGACGACCTCGATAAGGTGCTGCGGCCCGGGCTGCGGGTCCGGGGTCCA
 GGTCCGGGCGGGGAGGAGGCCACCCGGCTCGCTCGGGTTGCTGTGACGACTCAGCCCTGGCACGAAGC
 CCAGCCCGCGGCTGCTGGGCTCTCGGCTGAGCAAAATCTATTCCAGTCCACCCTGTACTGTCCACT
 GTGGCAAACGAGGCCACAATAACCTCGGGTGAAGAGGCCACGTCTCGGGTCCCTCCCCGTAGGT
 CTCATGTCCCAAGCCGCTGCCTGGAAGCCGCAACGGGACTGGGAAAATGGCAGCAATGCCTCTCC
 CCAGCGTCAGTGCCGAGTACACAGTCCACGGCTGTACAAGAACCAGCGCAAGTCCAACAAGTTC
 ATCATCCACAATGCCCTATCACACTGCTGCCTGGCGGCAAGGTGAACGAACCGCAGAAGAATCGCATT
 CTGGAGGAAATTGAGAAAAGCAAGGCCAACCACTTCTGATCCTTTTCGCGACTCGAGCTGCCAGTTC
 CGGGCGCTCTACACGCTGTCGGGGGAGACAGAGGAGCTGTGCGGCTGGCAGGATGGGCCCGGACC
 GTCACGCCGCCATGGTGAAGGCATCTACAAGTACAACCTCGGACCGCAAGCGCTTACCAGATCCCC
 GCCAAGACCATGTCCATGAGCGTCGATGCCTTACCATCCAGGGACACCTCTGGCAGGGCAAGAAACC
 ACCACTCCAAGAAGGGCGGCGGCACCCCCAAATAG
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

- Restriction Sites:** Sgfl-MluI
- ACCN:** NM_001080429
- Insert Size:** 3831 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080429.2</u>
RefSeq Size:	4201 bp
RefSeq ORF:	3831 bp
Locus ID:	57662
UniProt ID:	<u>Q9P1Y5</u>
Cytogenetics:	19p13.2
MW:	137.9 kDa

Gene Summary:

Key microtubule-organizing protein that specifically binds the minus-end of non-centrosomal microtubules and regulates their dynamics and organization (PubMed:19041755, PubMed:23169647). Specifically recognizes growing microtubule minus-ends and autonomously decorates and stabilizes microtubule lattice formed by microtubule minus-end polymerization (PubMed:24486153). Acts on free microtubule minus-ends that are not capped by microtubule-nucleating proteins or other factors and protects microtubule minus-ends from depolymerization (PubMed:24486153). In addition, it also reduces the velocity of microtubule polymerization (PubMed:24486153). Required for the biogenesis and the maintenance of zonula adherens by anchoring the minus-end of microtubules to zonula adherens and by recruiting the kinesin KIFC3 to those junctional sites (PubMed:19041755). Required for orienting the apical-to-basal polarity of microtubules in epithelial cells: acts by tethering non-centrosomal microtubules to the apical cortex, leading to their longitudinal orientation (PubMed:27802168, PubMed:26715742). Plays a key role in early embryos, which lack centrosomes: accumulates at the microtubule bridges that connect pairs of cells and enables the formation of a non-centrosomal microtubule-organizing center that directs intracellular transport in the early embryo (By similarity). Couples non-centrosomal microtubules with actin: interaction with MACF1 at the minus ends of non-centrosomal microtubules, tethers the microtubules to actin filaments, regulating focal adhesion size and cell migration (PubMed:27693509). Plays a key role in the generation of non-centrosomal microtubules by accumulating in the pericentrosomal region and cooperating with KATNA1 to release non-centrosomal microtubules from the centrosome (PubMed:28386021). Through the microtubule cytoskeleton, also regulates the organization of cellular organelles including the Golgi and the early endosomes (PubMed:28089391). Through interaction with AKAP9, involved in translocation of Golgi vesicles in epithelial cells, where microtubules are mainly non-centrosomal (PubMed:28089391).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) encodes the longer isoform (1).