

Product datasheet for **SC315962**

LILRB2 (NM_001080978) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LILRB2 (NM_001080978) Human Untagged Clone
Tag:	Tag Free
Symbol:	LILRB2
Synonyms:	CD85D; ILT-4; ILT4; LIR-2; LIR2; MIR-10; MIR10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_001080978 edited
 ATGACCCCATCGTCACAGTCTGATCTGTCTCGGGCTGAGTCTGGGCCCCAGGACCCAC
 GTGCAGACAGGGACCATCCCAAGCCCACCCTGTGGGCTGAGCCAGACTCTGTGATCACC
 CAGGGGAGTCCCCTCACCTCAGTTGTCAGGGGAGCCTTGAAGCCCAGGAGTACCGTCTA
 TATAGGGAGAAAAATCAGCATCTTGGATTACACGGATACGACCAGAGCTTGTGAAGAAC
 GGCCAGTTCACATCCCATCCATCACCTGGGAACACAGGGCGATATGGCTGTCAAGTAT
 TACAGCCGCGCTCGGTGGTCTGAGCTCAGTGACCCCTGGTGTGATGACAGGAGCC
 TACCCAAAACCCACCCTCTCAGCCAGCCAGCCCTGTGGTGACCTCAGGAGGAAGGGTG
 ACCCTCCAGTGTGAGTACAGGTGGCATTGGCGCTTATTCTGTGTAAGGAAGGAGAA
 GATGAACACCCACAATGCCTGAACTCCCAGCCCCATGCCCGTGGGTCGTCCCGGCCATC
 TTCTCCGTGGGCCCCGTGAGCCGAATCGCAGGTGGTCGCACAGGTGCTATGGTTATGAC
 TTGAACTCTCCATGTGTGGTCTTACCCAGTGTCTCCTGGAGCTCCTGGTCCCAGGT
 GTTTCTAAGAAGCCATCACTCTCAGTGCAGCCGGTCTGTCTGGCCCTGGGAAAGC
 CTGACCCTCCAGTGTCTGTGATGTCGGCTATGACAGATTTGTTCTGTACAAGGAGGGG
 GAACGTGACCTTCGCCAGCTCCCTGGCCGCGAGCCAGGCTGGGCTCTCCAGGCCAAC
 TTCACCCTGGGCCCTGTGAGCCGCTCTACGGGGCCAGTACAGATGCTACGGTGCATAC
 AACCTCTCCTCCGAGTGGTGGGCCCCAGCGACCCCTGGACATCCTGATCACAGGACAG
 ATCCATGGCACACCCCTTCATCTCAGTGCAGCCAGGCCACAGTGGCCTCAGGAGAGAAC
 GTGACCCTGTGTGTCAGTCAATGGCGCAGTTCACACTTTCCTTCTGACCAAGGCGGGA
 GCAGCTGATGCCCCACTCCGTCTAAGATCAATACACGAATATCCTAAGTACCAGGCTGAA
 TTCCCCATGAGTCTGTGACCTCAGCCACGCGGGGACCTACAGGTGCTACGGCTCACTC
 AACTCCGACCCCTACCTGCTGTCTCACCACAGTGGCCCTGGAGCTCGTGGTCTCAGGA
 CCCTCCATGGGTTCCAGCCCCCACCACCGGTCCCATCCACACCTGCAGGCCCTGAG
 GACCAAGCCCTCACCCCACTGGGTCGGATCCCCAAAGTGGTCTGGGAAGGCACCTGGGG
 GTTGTGATCGGCATCTTGGTGGCCGTCGTCTACTGCTCCTCCTCCTCCTCCTCTTTC
 CTCATCTCCGACATCGACGTCAGGGCAAACACTGGACATCAACCCAGAGAAAGGCTGAT
 TTCCAACATCCTGCAGGGGCTGTGGGGCCAGAGCCACAGACAGAGGCCTGCAGTGGAGG
 TCCAGCCAGCTGCCGACGCCAGGAAGAAAACCTCTATGCTGCCGTGAAGGACACACAG
 CCTGAAGATGGGGTGGAGATGGACTCGGGTGTGATCTGAAGCCCCCAGGATGTG
 ACCTACGCCAGCTGCACAGCTTGACCCTCAGACGGAAGGCAACTGAGCCTCCTCCATCC
 CAGGAAGGGGAACCTCCAGCTGAGCCAGCATCTACGCCACCCTGGCCATCCACTAG

Restriction Sites: Please inquire

ACCN: NM_001080978

Insert Size: 1800 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080978.1</u> , <u>NP_001074447.1</u>
RefSeq Size:	2910 bp
RefSeq ORF:	1794 bp
Locus ID:	10288
Cytogenetics:	19q13.42
Protein Families:	Transmembrane
Gene Summary:	<p>This gene is a member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular immunoglobulin domains, a transmembrane domain, and two to four cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs). The receptor is expressed on immune cells where it binds to MHC class I molecules on antigen-presenting cells and transduces a negative signal that inhibits stimulation of an immune response. It is thought to control inflammatory responses and cytotoxicity to help focus the immune response and limit autoreactivity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (2) uses an alternate in-frame splice site in the central coding region, compared to variant 1. The encoded isoform (2) is shorter, compared to isoform 1. Both variants 2 and 3 encode the same isoform.</p>