

Product datasheet for SC316050

FAM86B1 (NM_001083537) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM86B1 (NM_001083537) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM86B1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316050 representing NM_001083537. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGCGCCCGAGGAGAACGCGGGGACCGAACTCTTGCTGCAGGGTTTGTAGCGCCGCTTCTGGCGGTG
 CGCAGACTGCGCTCCTTCCCCTGGCAGAGCTTAGAGGCAAAGTTAAGAGACTCATCAGATTCTGAGCTG
 CTGCGGGATATTTTGCAGAAGACTGTGAGGCATCCTGTGTGTGAAGCACCCGCCGTCAAGTAT
 GCCTGGTGCTTTCTCTCAGAAGCTCATCAAAAGTCCTCAGGAGGCTCAGTCACACTCTCCAAGAGCACA
 GCCATCATCTCCACGGTACCACAGGCCTGGTCACATGGGATGCCGCCCTCTACCTTGAGAATGGGCC
 ATCGAGAACCCGCGCAGCCTTCATTAACAGGACTGTCCTAGAGCTTGGCAGTGGTGCCGGCTCACAGGC
 CTTGCCATCTGCAAGATGTGCCGCCCGGGCATACATCTTCAGCGACCTCACAGCCGGTCTCGAG
 CAGCTCCGAGGGAATGTCTTCTCAATGGCCTCTCATTAGAGGCAGACATCACTGGCAACTTAGACAGC
 CCCAGGGTGACAGTGGCCAGCTGGACTGGGACGTAGCAATGGTCCATCAGCTCTCTGCCTTCCAGCCA
 GATGTTGTATTGCAGCAGAGCTGCTGTATTGCCAGAAGCCATCGTGTGCTGGTGGGGTCTGCAG
 AGGCTGGTGCCTGCCGGGAGCACAAGCGGGCTCCTGAGGTCTACGTGGCCTTTACCGTCCGCAACCCA
 GAGACGTGCCAGCTGTTACCACCGAGCTAGGCCGGGATGGGATCAGATGGGAAGCGGAAGCTCATCAT
 GACCAGAACTGTTCCCTATGGAGAGCACTTGGAGATGGCAATGCTGAACCTCACACTGTAG
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001083537
Insert Size:	891 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001083537.1</u>
RefSeq Size:	2458 bp
RefSeq ORF:	891 bp
Locus ID:	85002
UniProt ID:	<u>Q8N7N1</u>
Cytogenetics:	8p23.1
MW:	32.9 kDa