

## Product datasheet for **SC316070**

### Striatin3 (STRN3) (NM\_001083893) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Striatin3 (STRN3) (NM_001083893) Human Untagged Clone
Tag:	Tag Free
Symbol:	Striatin3
Synonyms:	PPP2R6B; S/G2NA; SG2NA
Vector:	<u>pCMV6 series</u>



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001083893, the custom clone sequence may differ by one or more nucleotides

```

ATGGACGAGCTTGCCGGAGGCGGTGGTGGCGGCCCGGGATGGCGGCCCTCCCCGGCAG
CAGCAGGGACCTGGGGGAACTGGGCCTTCGCCGGGGGAACGGAGCGGGCGGGCGGC
GGGGTCTCCGGCTCCGAGGGAGCGGGTCCCGCGGCAGGCCCGAGCTGTCCCGGCCG
CAGCAGTACACTATCCCGGGATACTGCACTACATCCAGCAGAGTGGGCTCGGTTCCGAG
ATGGAGCGGGCGCACTGGGAGGTGGAACGGGCCGAAGTGCAGGCCCGGATTGCATTTCTA
CAAGGCGAAAGAAAAGGTCAAGAGAACCCTGAAGAAGGACTTAGTAAGAAGAATAAAGATG
TTAGAGTATGCATTAACAAGAAAGGGCAAAATATCACAAATTAATAATGGCACGGAA
CTGAACCAAGGTGACTTGAAAATGCCAACCTTTGAGTCAGAAGAAACAAAGACACAGAG
GCTCCACAGCACCTCAGAATAGCCAGTTAACGTGGAAGCAAGGCAGACAGCTTTAAGA
CAGTATCTTCAGGAAGTAGGTTATACAGATAACAATATTAGATGTACGGTCTCAGCGGGTA
AGGTCATTACTTGGACTATCTAATTCAGAACCAATGGATCAGTAGAAACAAAGAATTTA
GAACAGATCCTGAATGGAGGTGAATCTCCTAAGCAAAAGGGACAAGAAATAAAAAGTCC
TCTGGTGTGTTCTTGGACGTTCAATTTCTTAGAAAATGCCGATGACAGTGTGAAGAT
GAGGAAAATGACATGATCGAAGGCATCCAGAAGGAAAAGACAAACATCGGATGAATAAA
CATAAAATAGGTAATGAAGGTTTAGCTGCTGACCTAACTGACGATCCTGATACTGAGGAA
GCACTGAAAGAATTTGATTTTTAGTACTGCTGAAGATGGTGAAGGAGCTGGAGAAGCA
CGGAGTTCGGGGGATGGCACAGAATGGGATAAAGATGACCTCTCCCAACTGCTGAGGTT
TGGGATGTAGACCAGGGACTAATAAGTAACTGAAGGAACAGTACAAGAAGGAACGAAAG
GGGAAGAAAGGGTGAAGAGGGCCAACAGGACAAAACCTACGACATGATAGCTGATCTG
GGAGATGATGAGTGCACCATCCCTCAGGAATCATTAACTAGTCTAGGTCAGCCTCT
ACTAGAATGACTGATCATGAAGGTGCAAGAGCAGAGGAAGTGAACCAATAACGTTTCCA
TCTGGAGAGGCAAGTCATTTATTATGGGTTCTGATGATGTTTTGTTAAGTGTACTGGGC
CTTGGAGACCTTGACAGCTTGACGGTAACAAATGATGCAGACTATAGTTATGATTTGCC
GCTAATAAAGATGCCTTCGAAAGACATGGAATCCCAAGTATACACTACGTAGCCATTTT
GATGGAGTACGGGCATTAGCTTTTCATCCTGTAGAACCTGTGCTGGTTACTGCTTCTGAG
GACCATACCCTGAACTTTGGAACCTGCAAAAACAGTTCCTGCCAAAAGAGTGCCTCT
TTAGATGTAGAGCCTATCTACACATTTAGGGCCACATCGGCCCTGTTCTGTCATTAGCT
ATTAGTCTAATGGAGAACAGTGTTTTAGTGGTGGTATTGATGCAACCATCCAGTGGTGG
AATATGCCGAGTCCAGTGTAGATCCATATGATACATATGAGCCAAATGTTCTAGCTGGC
ACTTTAGTTGGTCATACAGATGCAGTTTGGGGTCTTGCTTATAGTGGCATAAAAATCAA
TTACTGTCTTGTTCAGCAGATGGCACTGTTAGGTTATGGAATCCACAAGAAAAATTGCCA
TGTATTTGCACCTTACAATGGAGATAAAAAGCATGGAATACCTACATCAGTTGACTTTATA
GGCTGTGATCCAGCTCATATGGTAACCTCTTCAACACTGGTAGTGACGTAATTTATGAT
TTAGAAACATCACAGTCATTGGTGATACTTTCATCACAGGTAGATTCTGGTTTACAATCT
AATAATCATATCAACAGAGTAGTAAGTCATCCCACTTCCTGTTACAATAACTGCTCAT
GAAGATAGACACATCAAATTTTTGACAATAAACGGGTAATAATGATCCATTCTATGGTA
GCTCACTTGGATGCTGTTACAAGTCTAGCAGTAGATCCTAATGGAATCTATTTGATGTC
GGAAGCCATGACTGTTCCATCAGATTATGGAATTTAGACAGCAAGACATGTGTGCAAGAA
ATAACAGCTCACAGAAAGAAATTGGATGAATCAATTTATGATGTTGCTTCCACTCGTCA
AAAGCATATAGCTAGTGCAGGAGCTGATGCTCTGCCAAAGTATTTGTA

```

**Restriction Sites:** Please inquire

**ACCN:** NM\_001083893

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001083893.1</a> , <a href="#">NP_001077362.1</a>
<b>RefSeq Size:</b>	4233 bp
<b>RefSeq ORF:</b>	2394 bp
<b>Locus ID:</b>	29966
<b>UniProt ID:</b>	<a href="#">Q13033</a>
<b>Cytogenetics:</b>	14q12
<b>Gene Summary:</b>	<p>Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1), also known as the beta isoform.</p>