

Product datasheet for SC316110

GAGE2B (NM_001098411) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GAGE2B (NM_001098411) Human Untagged Clone
Tag:	Tag Free
Symbol:	GAGE2B
Synonyms:	CT4.2; GAGE-2; GAGE-2B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316110 representing NM_001098411. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTTGGCGAGGAAGATCGACCTATCGGCCTAGACCAAGACGCTACGTAGAGCCTCCTGAAATGATT
GGGCCTATGCGGCCGAGCAGTTCAGTGATGAAGTGGAACCAGCAACACCTGAAGAAGGGGAACCAGCA
ACTCAACGTCAGGATCCTGCAGCTGCTCAGGAGGGAGAGGATGAGGGAGCATCTGCAGGTCAAGGGCCG
AAGCCTGAAGCTCATAGCCAGGAACAGGGTCACCCACAGACTGGGTGTGAGTGTGAAGATGGTCCTGAT
GGGCAGGAGATGGACCCGCAAAATCCAGAGGAGGTGAAAACGCCCTGAAGAAGGTGAAAAGCAATCACAG
TGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGCCCGC
```

Restriction Sites:	SgfI-MluI
ACCN:	NM_001098411
Insert Size:	351 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



View online »

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001098411.3</u>
RefSeq Size:	509 bp
RefSeq ORF:	351 bp
Locus ID:	645037
UniProt ID:	<u>Q13066</u>
Cytogenetics:	Xp11.23
MW:	12.8 kDa
Gene Summary:	Antigen, recognized on melanoma by autologous cytolytic T-lymphocytes.[UniProtKB/Swiss-Prot Function]