

Product datasheet for **SC316150**

zinc finger protein 655 (ZNF655) (NM_001085367) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	zinc finger protein 655 (ZNF655) (NM_001085367) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF655
Synonyms:	VIK; VIK-1
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001085367, the custom clone sequence may differ by one or more nucleotides ATGGAGGAAATACCAGCCCAGGAAGCAGCAGGGTCACCAAGGGTCCAGTTTCAGTCTTTG GAGACCCAGTCTGAGTGTCTGTCCCAGAGCCTCAGTTTGTGCAGGACACCGACATGGAA CAGGGACTCACTGGGGCTCCACCTGTTCCCTCAGGTGCCTGCTCTTCCCCGTGAGGGAAGC CCAGGAGACCAGGCAGCTGCGCTCTTGACAGCCAGGTACCAGGAGTTTGTGACATTTCGAG GATGTGGCTGTGCACCTTACTCGAGAGGAATGGGGATACCTGGACCCTGTTTCAGAGGGAC CTCTACAGAGAAGTGATGTTAGAGAATTATGGGAACGTGGTCTCACTGGGCATACTTCTC CGCCTTCCCACCACCCGATTTCATAGTGTGAATTCCTGCCGGCCCTGAGTCATACCCAG GCAAGTGTCTTCTGGAGAAACTTGCCTTACAGCAGGAATCTCCAAGAGATGG CCCAAGTATCGGCTTCCCATCGATATTGCTCGTCCCTGCTCGGAAACTCCTTTCCACGA TTG
Restriction Sites:	Please inquire
ACCN:	NM_001085367
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001085367.1](#), [NP_001078836.1](#)

RefSeq Size: 1435 bp

RefSeq ORF: 546 bp

Locus ID: 79027

UniProt ID: [Q8N720](#)

Cytogenetics: 7q22.1

Protein Families: Transcription Factors

Gene Summary: This gene encodes a zinc finger protein. The zinc finger proteins are involved in DNA binding and protein-protein interactions. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (10) has a distinct 5' UTR and includes alternate exons in the 3' coding region and 3' UTR, compared to variant 7. Both variants 2 and 10 encode the same isoform (b), which has a shorter and distinct C-terminus compared to isoform f.