

Product datasheet for SC316157

PTPMT1 (NM_175732) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PTPMT1 (NM_175732) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTPMT1
Synonyms:	DUSP23; MOSP; PLIP; PNAS-129
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene ORF sequence for NM_175732 edited ATGGCGCCACCGCGCTGCTGGAGGCCGGCCTGGCGGGGTGCTCTTCTACCCGACGCTG CTCTACACCCTGTTCCGCGGGAAGGTGCCGGTTCGGGCGCACCGGACTGGTACCACCGC ATCGACCCACCGTGCTGCTGGGCGCGCTGCCGTTGCGGAGCTTGACGCGCCAGCTGGTA CAGGACGAGAACGTGCGCGGGTGATCACCATGAACGAGGAGTACGAGACGAGGTTCTCTG TGCAACTCTTACAGGAGTGGAAGAGACTAGGAGTCGAGCAGCTGCGGCTCAGCACAGTA GACATGACTGGGATCCCCACCTTGGACAACCTCCAGAAGGGAGTCCAATTTGCTCTCAAG TACCAGTCGCTGGGCCAGTGTGTTTACGTGCATTGTAAGGCTGGGCGCTCCAGGAGTGCC ACTATGGTGGCAGCATACCTGATTCAGGTGCACAAATGGAGTCCAGAGGAGGCTGTAAGA GCCATCGCAAGATCCGGTCATACATCCACATCAGGCCTGGCCAGCTGGATGTTCTTAAA GAGTTCCACAAGCAGATTACTGCACGGCAACAAAGGATGGGACTTTTGCATTTCAAAG ACATGA
Restriction Sites:	Please inquire
ACCN:	NM_175732
Insert Size:	2700 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_175732.1</u> , <u>NP_783859.1</u>
RefSeq Size:	859 bp
RefSeq ORF:	606 bp
Locus ID:	114971
UniProt ID:	<u>Q8WUK0</u>
Cytogenetics:	11p11.2
Protein Families:	Druggable Genome
Gene Summary:	<p>Lipid phosphatase which dephosphorylates phosphatidylglycerophosphate (PGP) to phosphatidylglycerol (PG) (By similarity). PGP is an essential intermediate in the biosynthetic pathway of cardiolipin, a mitochondrial-specific phospholipid regulating the membrane integrity and activities of the organelle (By similarity). Has also been shown to display phosphatase activity toward phosphoprotein substrates, specifically mediates dephosphorylation of mitochondrial proteins, thereby playing an essential role in ATP production (By similarity). Has probably a preference for proteins phosphorylated on Ser and/or Thr residues compared to proteins phosphorylated on Tyr residues (By similarity). Probably involved in regulation of insulin secretion in pancreatic beta cells (By similarity). May prevent intrinsic apoptosis, probably by regulating mitochondrial membrane integrity (PubMed:24709986).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) is the longer transcript.</p>