

Product datasheet for **SC316183**

GPR18 (NM_001098200) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GPR18 (NM_001098200) Human Untagged Clone
Tag:	Tag Free
Symbol:	GPR18
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001098200, the custom clone sequence may differ by one or more nucleotides

```
ATGATCACCCCTGAACAATCAAGATCAACCTGTCCCTTTAACAGCTCACATCCAGATGAATACAAAATTG
CAGCCCTGTCTTCTATAGCTGTATCTTCATAATTGGATTATTTGTTAACATCACTGCATTATGGGTTTT
CAGTTGTACCACCAAGAAGAGAACCACGGTAACCATCTATATGATGAATGTGGCATTAGTGGACTTGATA
TTTATAATGACTTTACCCTTTCGAATGTTTTATTATGCAAAGATGAATGGCCATTTGGAGAGTACTTCT
GCCAGATTCTGGAGCTCTCACAGTGTTTACCCAAGCATTGCTTTATGGCTTCTTGCCTTTATTAGTGC
TGACAGATACATGGCCATTGTACAGCCGAAGTACGCCAAAGAACTAAAAACACGTGCAAAGCCGTGCTG
GCGTGTGGGAGTCTGGATAATGACCCGTGACCAGCACCCTCTGCTACTGCTCTATAAAGACCCAG
ATAAAGACTCCACTCCCGCCACCTGCCTCAAGATTTCTGACATCATCTATCTAAAAGCTGTGAACGTGCT
GAACCTCACTCGACTGACATTTTTTTTCTTGATTCCTTTGTTTCATCATGATTGGGTGCTACTTGGTCATT
ATTCATAATCTCCTTCACGGCAGGACGTCTAAGCTGAAACCCAAAGTCAAGGAGAAGTCCATAAGGATCA
TCATCACGCTGCTGGTGCAGGTGCTCGTCTGCTTTATGCCCTTCCACATCTGTTTCGCTTTCCTGATGCT
GGGAACGGGGGAGAACAGTTACAATCCCTGGGGAGCCTTACCACCTTCCATGAACCTCAGCACGTGT
CTGGATGTGATTCTACTACATCGTTTCAAACAATTTCAAGGCTCGAGTCATTAGTGCATGCTATAACC
GTAATTACCTTCGAAGCATGCGCAGAAAAAGTTCCGATCTGGTAGTCTACGGTCACTAAGCAATATAAA
CAGTGAATGTTATGA
```

Restriction Sites:	NotI-NotI
ACCN:	NM_001098200
Insert Size:	1400 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	The ORF of this clone has been fully sequenced and found to be a perfect match to NM_001098200.1.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001098200.1 , NP_001091670.1
RefSeq Size:	1484 bp
RefSeq ORF:	996 bp
Locus ID:	2841
UniProt ID:	Q14330
Cytogenetics:	13q32.3
Protein Families:	Druggable Genome, GPCR, Transmembrane
Gene Summary:	<p>Receptor for endocannabinoid N-arachidonyl glycine (NAGly) (PubMed:16844083, PubMed:24762058, PubMed:27572937). However, conflicting results about the role of NAGly as an agonist are reported (PubMed:27018161). Can also be activated by plant-derived and synthetic cannabinoid agonists (PubMed:24762058). The activity of this receptor is mediated by G proteins which inhibit adenylyl cyclase (PubMed:16844083). May contribute to regulation of the immune system. Is required for normal homeostasis of CD8+ subsets of intraepithelial lymphocytes (IELs) (CD8alphaalpha and CD8alphabeta IELs)in small intestine by supporting preferential migration of CD8alphaalpha T-cells to intraepithelial compartment over lamina propria compartment, and by mediating their reconstitution into small intestine after bone marrow transplant (By similarity). Plays a role in hypotensive responses, mediating reduction in intraocular and blood pressure (By similarity). Mediates NAGly-induced process of reorganization of actin filaments and induction of acrosomal exocytosis (PubMed:27572937). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Variants 1 and 2 encode the same protein.</p>