

Product datasheet for **SC316191**

CDK10 (NM_052988) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDK10 (NM_052988) Human Untagged Clone
Tag:	Tag Free
Symbol:	CDK10
Synonyms:	ALSAS; PISSLRE
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316191 representing NM_052988. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGAGCCAGATCTGGAGTGCAGCAGATCCGTCTGAAGTGTATTCGTAAGGAGGGCTTCTTCAGC
GTGCCCTCCGGAACACAGGCTGGGACGATGCCGGAGTGTGAAGGAGTTTGAGAAGCTGAACCCGATTGGA
GAGGGTACCTACGGCATTGTGTATCGGGCCCGGGACACCCAGACAGATGAGATTGTCGCACTGAAGAAG
GTGCGGATGGACAAGGAGAAGGATGGCATCCCCATCAGCAGCTTGCAGGAGATCACGCTGCTGCTCCGC
CTGCGTCATCCGAACATCGTGGAGCTGAAGGAGGTGGTTGTGGGAACCACCTGGAGAGCATCTTCTG
GTGATGGGTTACTGTGAGCAGGACCTGGCCAGCCTCCTGGAGAATATGCCAACACCCCTTCTCGGAGGCT
CAGGTCAAGTGCATCGTGTGCAGGTGCTCCGGGGCCTCCAGTATCTGCACAGGAACTTCATTATCCAC
AGGGACCTGAAGGTTTCCAACCTTGTCTATGACCGACAAGGTTGTGTGAAGACAGCGGATTCGGCCCTG
GCCCGGGCCTATGGTGTCCCAGTAAAGCCAATGACCCCAAGGTGGTCACTCTCTGGTACCGAGCCCT
GAAGTGTGTTGGGAACCAACAGCAGACCAACAGCATCGACATGTGGGCTGTGGGCTGCATACTGGCC
GAGCTGCTGGCGCACAGGCTCTTCTCCCGGCACTTCCGAGATCCACCAGATCGACTTGATCGTGCAG
CTGCTGGGCACGCCAGTGAACAACCTGAAGCACAAGTTCCCATGGCTGTGCGAGGCGGGCTGCGCCTG
CTCCGGAAGCAGCCCTACAACAACCTGAAGCACAAGTTCCCATGGCTGTGCGAGGCGGGCTGCGCCTG
CTGCATTCCTGTTTCATGTACGACCCTAAGAAAAGGGCGACGGCCGGGACTGCTGGAGAGCTCCTAT
TTCAAGGAGAAGCCCTACCCTGTGAGCCGGAGCTCATGCCGACCTTTCCCAACCACCGCAACAAGCGG
GCCGCCAGCCACCTCCGAGGGCCAGAGCAAGCGCTGTAACCCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_052988



[View online »](#)

Insert Size:	1083 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_052988.4
RefSeq Size:	1815 bp
RefSeq ORF:	1083 bp
Locus ID:	8558
UniProt ID:	Q15131
Cytogenetics:	16q24.3
Protein Families:	Druggable Genome, Protein Kinase
MW:	41 kDa
Gene Summary:	<p>The protein encoded by this gene belongs to the CDK subfamily of the Ser/Thr protein kinase family. The CDK subfamily members are highly similar to the gene products of <i>S. cerevisiae</i> <i>cdc28</i>, and <i>S. pombe</i> <i>cdc2</i>, and are known to be essential for cell cycle progression. This kinase has been shown to play a role in cellular proliferation and its function is limited to cell cycle G2-M phase. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2009]</p> <p>Transcript Variant: This variant (a) represents the longest transcript and encodes the longest isoform (a).</p>