

Product datasheet for SC316198

CBWD6 (NM_001085457) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CBWD6 (NM_001085457) Human Untagged Clone
Tag:	Tag Free
Symbol:	CBWD6
Synonyms:	CBWD7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316198 representing NM_001085457. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
 ATGTTACCGGCTGTTGGATCTGTGGATGAGGAAGAGGATCCTGCGGAGGAGGATTGTCCTGAATTGGTT
 CCCATTGAGACGACGCAAAGCGAGGAGGAGGAAAAGTCTGGCCTCGGCGCCAAGATCCAGTCACAATT
 ATCACC GG GTATTTAGGTGCTGGGAAGACAACACTTCTGAACATATTTTACAGAGCAACATAGTAAA
 AGAGTAGCGGTCAATTTAAATGAATCTGGGAAGGAAGTGCCTGGAGAAATCCTTAGCTGTCAGCCAA
 GGTGGAGAGCTCTATGAAGAGTGGCTGGAACCTAGAAACGGTTGCCTCTGCTGTTCAAGGACAAT
 GGCCTTAGAGCTATTGAGAATTTGATGCAAAAGAAGGGGAAATTTGATGACATACTGTTAGAGACCACT
 GGATTAGCAGACCCTGGTGCAGTGGCTTCTATGTTTTGGGTTGATGCCGAATTAGGAGTGATATTAC
 CTTGATGGTATCATACTATTGTGGATTCAAAATATGGATTAAACATTTAACAGAGAGAAACCTGAT
 GGCTTATCAATGAAGCTACTAGGCAAGTTGCTTTGGCAGATATCATTCTCATTAAATAAACAGACTTG
 GTCCCAGAAGAAGATGTAAAGAAATTAAGAACGACACTTAGATCCATAAATGGACTAGGACAAATCTTA
 GAAACACAAGATCAAGAGTTGATCTCTCTAATGTATTAGATCTTCATGCCTTTGATAGTCTCTCTGGA
 ATAAGTTTGCAGAAAAAATTCAGCATGTGCCAGGAACACAACCTCACCTTGATCAGAGTATTGTTACA
 ATCACATTTGACGTACCAGGAAATGCAAAGGAAGAATCTTAATATGTTTATTGAGAAATCTTCTGTGG
 GAAAAGAAATGTGAGAAACAAGGACAATCACTGCATGGAGGTCATAAGGCTGAAGGATTGGTGTCAATC
 AAAGACAAATCACAAACAAGTGATTGTCCAGGGTGTCATGAGCTCTGTGATCTGGAGGAGACTCCAGTG
 AGCTGGAAGGATGACACTGAGAGAACAAATCGATTGGTCCTCATTGGCAGAAATTTAGATAAGGATATC
 CTTAAACAGCTGTTTATAGCTACTGTGACAGAAACAGAAAGCAGTGACAAACACATTTCAAAGAAGAT
 CAAGTTTGTACATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

Restriction Sites: SgfI-MluI


[View online »](#)

ACCN:	NM_001085457
Insert Size:	1188 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001085457.1</u>
RefSeq Size:	1719 bp
RefSeq ORF:	1188 bp
Locus ID:	644019
UniProt ID:	<u>Q4V339</u>
Cytogenetics:	9p11.2
MW:	44 kDa