

Product datasheet for SC316199

C5ORF44 (TRAPPC13) (NM_001093756) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C5ORF44 (TRAPPC13) (NM_001093756) Human Untagged Clone
Tag:	Tag Free
Symbol:	C5ORF44
Synonyms:	C5orf44; SHLD3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316199 representing NM_001093756. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGAAGTGAATCCCCCTAAACAGGAGCACCTGCTGGCGCTAAAAGTGATGCGGCTGACTAAGCCTACT
 TTATTCACCAATATCCAGTAACATGTGAAGAGAAAGACTTACCTGGAGATCTCTTTAACCAGCTGATG
 AGAGATGATCCTTCAACCGTTAATGGTGCAGAAGTTTAAATGTTGGGAGAAATGCTGACTTTACCACAG
 AATTTTGGGAATATATTTTGGGAGAGACCTTTCCAGTTATATCAGCGTTCATAATGATAGCAATCAA
 GTTGTAAGACATATTAGTAAAGCTGATCTTCAGACAAGTTCTCAGCGTTTAAATCTTTCAGCCTCC
 AATGCTGCAGTGGCTGAACCTAAACCGGATTGTTGTATTGATGATGTCATACATCATGAAGTCAAAGAA
 ATTGGAACACACATCTTGGTATGTGCTGTGAGTTATACAACCTCAGGCTGGAGAAAAATGTATTTTCA
 AAATCTTCAAATTTTCAAGTTCTCAAACCATTTGGATGTGAAAACCAAATTTTACAATGCAGAGACTGAT
 GAAGTATTTCTGGAAGCCAGATTCAGAAATATGACAACCTCACCTATGTTTATGGAGAAGGTTTCACTG
 GAGCCATCTATTATGTACAATGTAAACAGAAATAAATTCAGTCAGCCAAGCTGGAGAATGTGTGTCTACG
 TTTGGGTCAAGAGCATATTTGCAACCAATGGATACAGCCAGTACTTATACTGCCTAAAGCCAAAGAAT
 GAATTTGCAGAAAAAGCAGGCATCATTAAAGGGAGTAACAGTAATTGGAATTTGGATATAGTATGGAAA
 ACAATCTAGGTGAAAGGGGAAGTTACAGACCAGCCAACCTCAAAGAATGGCTCCAGGTTATGGAGAT
 GTTAGGTTGTCTTTGGAGGCAATACCAGATACCGTAAACCTTGAAGAACCTTTTCATATTACCTGTAA
 ATAACAACTGCAGTGAAAGGACTATGGATCTGGTTTTGGAAATGTGAATACCAATCCATCCACTGG
 TGTGGAATTTCAAGAACAGCTGGGAAAGCTGCATCCAAGTTCTTCGCTCTGTCTTGCCCTTACTCTG
 CTTTCTTCAGTACAGGACTGCAAAGCATCTCTGGCTTAAGACTAACAGACATTTCTAAAGAGAACA
 TATGAATATGATGACATCGCACAAGTCTGTGTGGTATCTTCTGCCATTAAAGTGGAAAGCTGA
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

Restriction Sites: Sgfl-MluI



[View online »](#)

ACCN:	NM_001093756
Insert Size:	1236 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001093756.1</u>
RefSeq Size:	3099 bp
RefSeq ORF:	1236 bp
Locus ID:	80006
UniProt ID:	<u>A5PLN9</u>
Cytogenetics:	5q12.3
MW:	45.9 kDa