

Product datasheet for **SC316205**

hnRNP F (HNRNPF) (NM_001098206) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	hnRNP F (HNRNPF) (NM_001098206) Human Untagged Clone
Tag:	Tag Free
Symbol:	hnRNP F
Synonyms:	HNRPF; mcs94-1; OK/SW-cl.23
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	<p>>NCBI ORF sequence for NM_001098206, the custom clone sequence may differ by one or more nucleotides</p> <pre> ATGATGCTGGGCCCTGAGGGAGGTGAAGGCTTTGTGGTCAAGCTCCGTGGCCTGCCCTGG TCCTGCTCTGTTGAGGACGTGCAGAACTTCCTCTCTGACTGCACGATTCATGATGGGGCC GCAGGTGTCCATTTTCATCTACACTAGAGAGGGCAGGCAGAGTGGTGAGGCTTTTGTGAA CTTGATCAGAAGATGATGTAATAATGGCCCTGAAAAAGACAGGGAAGCATGGGACAC CGGTACATTGAGGTGTTCAAGTCCACAGAACGAGATGGATTGGGTGTTGAAGCACAGT GGTCCCAACAGTGCCGACAGCGCCAACGATGGCTTCGTGCGGCTTCGAGGACTCCCATT GGATGCACAAAGGAAGAAATTGTTCAAGTCTTCTCAGGGTTGGAATTGTGCCAAACGGG ATCACATTGCCTGTGGACCCGAAGCAAGATTACAGGGGAAGCGTTCGTGCAGTTTGCC TCGCAGGAGTTAGCTGAGAAGGCTCTAGGGAACACAAGGAGAGGATAGGGCACAGGTAC ATTGAGGTGTTTAAGAGCAGCCAGGAGGAAGTTAGGTCATACTCAGATCCCCCTCTGAAG TTCATGTCCGTGCAGCGGCCAGGGCCCTATGACCGGCCCGGGACTGCCAGGAGGTACATT GGCATCGTGAAGCAGGCAGGCCTGGAAAGGATGAGGCCTGGTGCTACAGCACAGGCTAC GGGGGCTACGAGGAGTACAGTGGCCTCAGTGATGGCTACGGCTTACCACCGACCTGTT GGGAGAGACCTCAGCTACTGTCTCTCCGGAATGTATGACCACAGATACGGCGACAGTGAG TTCACAGTGCAGAGCACCACAGGCCACTGTGTCCACATGAGGGGCTGCCGTACAAAGCG ACCGAGAACGACATTTACAACCTTCTCTCTCCTCAACCCTGTGAGAGTCCATATTGAG ATTGGCCAGATGGAAGAGTGACGGGTGAAGCAGATGTTGAGTTTGCTACTCATGAAGAA GCTGTGGCAGCTATGTCCAAAGACAGGGCCAATATGCAGCACAGATATATAGAACTCTTC TTGAATTCAACAACAGGGGCCAGCAATGGGGCGTATAGCAGCCAGGTGATGCAAGGCATG GGGGTGTCTGCTGCCAGGCCACTTACAGTGGCCTGGAGAGCCAGTCAGTGAGTGGCTGT TACGGGGCCGGCTACAGTGGGCAGAACAGCATGGGTGGCTATGAC </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001098206
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001098206.1 , NP_001091676.1
RefSeq Size:	2667 bp
RefSeq ORF:	1248 bp
Locus ID:	3185
UniProt ID:	P52597
Cytogenetics:	10q11.21
Gene Summary:	<p>This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins that complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and regulate alternative splicing, polyadenylation, and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has three repeats of quasi-RRM domains that bind to RNAs which have guanosine-rich sequences. This protein is very similar to the family member hnRPH. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR, compared to variant 1. Variants 1-6 encode the same protein.</p>