

Product datasheet for SC316210

C5ORF44 (TRAPPC13) (NM_001093755) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C5ORF44 (TRAPPC13) (NM_001093755) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRAPPC13
Synonyms:	C5orf44; SHLD3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316210 representing NM_001093755. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTACTATAGGGCCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAAGTGAATCCCCTAACAGGAGCACCTGCTGGCGCTAAAAGTGATGCGGCTGACTAAGCCTACT
TTATTCACCAATATCCCAGTAACATGTGAAGAGAAAGACTTACCTGGAGATCTCTTTAACCGCTGATG
AGAGATGATCCTTCAACCGTTAATGGTGCAGAAGTTTTAATGTTGGGAGAAATGCTGACTTTACCACAG
AATTTTGGGAATATATTTTGGGAGAGACCTTTCCAGTTATATCAGCGTTCATAATGATAGCAATCAA
GTTGTAAGACATATTAGTAAAAGCTGATCTTCAGACAAGTTCTCAGCGTTTAAATCTTTCAGCCTCC
AATGCTGCAGTGGCTGAACCTAAACCGGATTGTTGATTGATGATGTCATACATCATGAAGTCAAAGAA
ATTGGAACACACATCTTGGTATGTGCTGTGAGTTATACAACTCAGGCTGGAGAAAAATGTATTTTACA
AAATCTTCAAATTTTCAAGTTCTCAAACCATTTGGATGTGAAAACCAAATTTTACAATGCAGAGAGTGAC
CTCAGTTCTGTGACTGATGAAGTATTTCTGGAAGCCAGATTCAGAAATATGACAACCTCACCTATGTTT
ATGGAGAAGGTTTCACTGGAGCCATCTATTATGTACAATGTAAACAGAATTAATTCAGTCAGCCAAGCT
GGAGAATGTGTCTACGTTTGGGTCAAGAGCATATTTGCAACCAATGGATACACGCCAGTACTTATAC
TGCCTAAAGCCAAAGAATGAATTTGCAGAAAAAGCAGGCATCATTAAAGGAGTAAACAGTAATTGGAAAA
TTGGATATAGTATGAAAAACAACTAGGTGAAAGGGGAAAGTTACAGACCAGCCAACCTTCAAAGAATG
GCTCCAGGTTATGGAGATGTTAGGTTGCTTTGGAGGCAATACCAGATACCGTAAACCTTGAAGAACCT
TTTCATATTACCTGTAATAAACAACACTGCAGCAGTGAAAGGACTATGGATCTGGTTTTGGAAATGTGC
AATACCAATTCATCCACTGGTGTGGAATTTTCAAGAACAGCTGGGAAAGCTGCATCCAAGTTCTTCG
CTCTGTCTTGCCTTACTCTGCTTTCTTCAAGTACAGGACTGCAAAGCATCTCTGGCTTAAAGACTAACA
GACACATTTCTAAAGAGAACATATGAATATGATGACATCGCACAAAGTCTGTGTGGTATCTTCTGCCATT
AAAGTGGAAAGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001093755

Insert Size: 1257 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001093755.1](#)

RefSeq Size: 3120 bp

RefSeq ORF: 1257 bp

Locus ID: 80006

UniProt ID: [A5PLN9](#)
Cytogenetics: 5q12.3
MW: 46.6 kDa