

Product datasheet for SC316222

MAB21L4 (NM_001085437) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	MAB21L4
Synonyms:	C2orf54
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316222 representing NM_001085437. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAAATTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGCCCCCTGCTCTCCCCACCTCAGCCATGGCCGTGCAGGTGCCCTGTGGCACCACCTACCTGCAG
GCCATCCGGTCCCGGGAGGCGCCGCGTGCACGAGCTTCCAGCGCGCAGAGAACGTGCTGCTCACGGTG
CTGGAGCGCGTGCATGCCCTGGACCCCGCTTCATCGTGGACTACTCCCGTGGCTGGAGGCTTCCAG
TTCGCCCTGCGCTCCTCTGAGGACCAATGGACATGGAGGTGCCCTGTGGGTGGATGCCGAGGCTCTA
CTGATTGAGGAACAGAGGCCACCCAGCGGAGGATGGCCTTGAATTATGCCACCTGGGTGTGCCCAGG
GAAGGGGCTGGCCTGGAGCGGTGGACACCGAGGATACCTTCACTGCCTCCTCGGAGGGTGACGCCAAG
TGCCGTGGACACATTGTGCCAGCAAGGTCTGTGTCTCAAGGACCTGCTGGTAGCAGCCATCGTA
CACTGCAAGCACCACAGTCTCATCGCACCAGGTTGCTGAACGCGGCCAGCCTGAGGGAGGAGCAGCTC
CACCTGTCTTGTGGTGTCCAGCGGCTGGAGAACAATCAGCTTCCACGTGGTGGCCGTGGTGAAGG
AAGCTTGGGGCGCCTGCCCTGGAGGGGGTGACGAGATGCCCGGATTCCCTGAGGGCAGCTTGAGAAGG
ATCCTCAGCAAGGGGTGGACCTGGTGCCGGCCAGCGCCAGCTCTGGAGGACCTCCACTGACTACCTG
CTCAGGAGGCTGCTGGGGAGCTGGGCTCCCTGCAGGGTCACTCGCCTGGACAGCCTCTCCATCCTCGAC
CGGGTCAACCACGAGAGCTGGCGTGACAGTGCCAGACTGACGGCTGACCTTTGGCCACCTGAAGATG
GTGCTGCTGTGGGCTCTGTGCTCTTCCCGCGCCCGAGGACTGGGCAAGACTGCAGGGCGCGCTGTAC
CGCCTGCTGGTGGTGTGCTCTGTTGCCTGGCCACGCGGAAGCTGCCCCACTTCTCCACCCGACGCGC
AACCTGCTGCAGGGCAGCGGCTGGACCTTGGTGCCATCTACCAGCGCTGGAGGGCTTCGCCAGCCAG
CCCAGGGCGGCCCTGCGCATCCAGCCACCCACCTGGCCGCGAGCCCCCGCCGCGCATCGGCTCCGGG
CTCAAGGCGCTCCTGCAGCTGCCAGCAGTGACCCCACTTACTGGGCACTGCCTACTTCGACGTCTG
CTGGACAAGTTCAGGTCTTCAACATCCAGGATAAGGACCGGATCTCTGCCATGCAGAGCATCTTCCAG
AAGACCAGGACTCTGGGAGGCGAGGAGAGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

```



Restriction Sites:	Sgfl-MluI
ACCN:	NM_001085437
Insert Size:	1344 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001085437.2
RefSeq Size:	2542 bp
RefSeq ORF:	1344 bp
Locus ID:	79919
UniProt ID:	Q08AI8
Cytogenetics:	2q37.3
MW:	49.5 kDa