

Product datasheet for SC316224

PRAMEF15 (NM_001098376) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRAMEF15 (NM_001098376) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRAMEF15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316224 representing NM_001098376. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGATGAGCATCCGGACTCCACCCAGACTCCTGGAGCTTGCAGGGCGGAGCCTGCTGAGGGACCAA
GCTTTGGCCATGTCCACCCTGGAGGAGCTGCCACAGAACTTTCCCCCACTGTTTCATGGAGGCCCTTC
AGCAGGAGACGCTGTGAGGCCCTGAAGCTGATGGTGCAGGCCTGGCCCTTCGCCGCCTCCCTCTGAGG
CCTCTGATAAAGATGCCTTGTCTGGAGGCCCTCCAAGCTGTGCTCGATGGGCTTGATGCACTGCTTACC
CAAGGGGTTCTGCCAGGAGGTGAAACTCCAAGTGTGGATTTACAGGATGTCTGTGAGAACTTCTGG
ATGGTTTGGTCTGAAGCTATGGCCATGGGTGCTTCTCAATGCCAAGAGGAACAAAAACAGTGCCAG
GACTGTCCAAGGATGAGAGGACGGCAGCCCTTGACTGTGTTCTGTAAGAACTTTGGCTCAAGAACAGGACT
CTGGATGAATACCTCACCTACCTCCTTCTATGGGTCAAGCAGAGGAAAGATTTACTACACCTGTGCTGT
AAGAAGCTGAAAATTTTGGGAATGCCCTTCGCAATATCAGAAGCATCCTGAAAATGGTGAACCTAGAC
TGTATCCAGGAGGTGGAAGTGAATTGCAAGTGGTACTGCCCATCCTGACACAGTTTACCCATACCTG
GGCCACATGAGGAATCTCAGAAGCTCGTTCTCTCCACATGGATGTCTCTGCTACGTTTCCCCAGAG
CAGAAGAAGGAGATTGTTACCCAGTTCACCACTCAGTTCTCAAGCTGCGCTGCCTCCAAAAGCTTTAT
ATGAAGTCTGTTTCTTCTCGAAGGCCACCTGGACCAGCTGCTCAGCTGTCTGAAGACCTCGTTAAAA
GTCCTCACAATACTAAGTGTGTCTTTTGAATCAGACTGAAGCATCTATCCAGTGCCCGAGTATC
AGTCAACTAAAGACCCTGGACCTGAGTGGCATCAGACTGACCAATTATAGTCTTGTGCTCTCCAAATT
CTCCTAGAAAAAGTTGCAGCCACCTTGAGTACCTGGATTTAGATGACTGTGGCATCATAGACTCCCAA
GTCAACGCCATCCTGCCTGCCCTGAGCCGCTGCTTTGAGCTCAACACCTTCAGCTTCTGTGGAATCCC
ATCTGCATGGCCACCCTGGAGAAGCTGTGAGCCACACAATCATACTCAAAAAGTATGTGTGGAGCTG
TATCCTGCCCCCGAGAGAGTTATGGTGTCTGATGGTACTCTGCTGGAGCAGATTTGCTCAAATTAGG
GCTGAGCTGATGAACAGAGTGAGGGACTTAAGGCACCCAAGAGGATCTTGTCTGTACTGACTACTGC
CCTGACTGTGGCAACAGGTCAATTTATGACCTGGAGGCAGATCAATACTGCTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCG



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001098376
Insert Size:	1437 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001098376.2</u>
RefSeq Size:	1862 bp
RefSeq ORF:	1437 bp
Locus ID:	653619
UniProt ID:	<u>Q5VWM5</u>
Cytogenetics:	1p36.21
MW:	55.4 kDa