

## Product datasheet for **SC316225**

### **CHES1 (FOXN3) (NM\_001085471) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	CHES1 (FOXN3) (NM_001085471) Human Untagged Clone
Tag:	Tag Free
Symbol:	FOXN3
Synonyms:	C14orf116; CHES1; PRO1635
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC316225 representing NM\_001085471.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGTCCAGTCATGCCTCCCAGTAAGAAGCCAGAAAGCTCAGGAATTAGTGTCTCCAGTGGACTGAGT
CAGTGTACGGGGCAGCGGTTTCTCCAAGGCCCTTCAGGAAGACGATGACCTCGACTTTTCTGTGCT
GACATCCGATTAGAAGAGGGGGCCATGGAAGATGAAGAGCTGACCAACCTGAACTGGCTGCACGAGAGC
AAGAACTTGCTGAAGAGCTTTGGGGAGTCGGTCTCAGGAGTGTAGCCCGTCCAGGACCTGGACGAT
GACACCCCCCATCCCCTGCCACTCTGACATGCCCTACGATGCCAGGCAGAACCCCAACTGCAAACCC
CCCTACTCCTCAGCTGCCTCATATTTATGGCCATCGAGGACTCTCCAACCAAGCGCCTGCCAGTGAAG
GATATCTACAACGGATCTTGGAACTTTCCGTATTTTGCAAATGCACCTACTGGGTGGAAAACTCA
GTGAGACACAATTTATCATTGAATAAGTGTTTAAGAAAGTGGACAAAGAGAGGAGTCAGAGTATTGGG
AAAGGTGCTTGTGGTGCATAGACCCAGAGTATAGACAAAATCTAATTCAGGCTTTGAAAAGACACCT
TATCACCCACACCCACAGTGTTCATACACCTCCACCTGTCTCAGGCATATCAAAGCACATCAGGT
CCACCCATCTGGCCGGGAGTACCTTCTCAAGAGAAATGGAGCCCTTCTCCAAGATCCTGACATTGAT
GCTGCCAGTGCCATGATGCTTTTGAATACTCCCCCTGAGATACAAGCAGGTTTTCTCCAGGAGTGATC
CAAATGGAGCGCGGGTCTGAGCCGAGGGCTGTTTCTGGCGTGCAGCCGCTGCCAATCACTCCCATT
GGGGTGACAGCGGCCATGAGGAATGGCATCACCAGCTGCCGGATGCGGACTGAGAGTGAAGCCATCTGT
GGCTCCCAGTGGTCCAGCGGAGACCCCAAGGAGGATCACAACTACAGCAGTGCCAAGTCTCCAACGCC
CGGAGCACCTCGCCACACGACTCCATCTCCTCCTCCTCCTCAGCCGACGACCACTATGAGTTT
GCCACCAAGGGGAGCCAGGAGGGCAGCGAGGGCAGCGAGGGGAGCTTCCGGAGCCAGGAGGCCCCAGC
GACACGGGAAGAGGACGACAGGAAGCACAGCCAGGAAGGAGCCCAAGGATTCTCTGGGGGACAGCGGGTAC
GCATCCCAGCACAGAAGCGCCAGCACTTCGCCAAGGCCAGGAAGTCCCAGCGACACACTGCCCCCTC
AAAAAGAGACGCACCGAAAAGCCCCCGAGAGCGATGATGAGGAGATGAAAGAAGCGGCAGGGTCCCTC
CTGCACTTAGCAGGGATCCGGTCTGTTTGAATAACATCACCAATCGGACGGCAAAGGGGCAGAAAGAG
CAAAGGAAACCACAAAAATTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001085471

<b>Insert Size:</b>	1473 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001085471.1</a>
<b>RefSeq Size:</b>	7856 bp
<b>RefSeq ORF:</b>	1473 bp
<b>Locus ID:</b>	1112
<b>UniProt ID:</b>	<a href="#">O00409</a>
<b>Cytogenetics:</b>	14q31.3-q32.11
<b>Protein Families:</b>	Druggable Genome, Transcription Factors
<b>MW:</b>	53.8 kDa
<b>Gene Summary:</b>	<p>This gene is a member of the forkhead/winged helix transcription factor family. Checkpoints are eukaryotic DNA damage-inducible cell cycle arrests at G1 and G2. Checkpoint suppressor 1 suppresses multiple yeast checkpoint mutations including mec1, rad9, rad53 and dun1 by activating a MEC1-independent checkpoint pathway. Alternative splicing is observed at the locus, resulting in distinct isoforms. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).</p>