

Product datasheet for **SC316235**

HMGCS1 (NM_001098272) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HMGCS1 (NM_001098272) Human Untagged Clone
Tag:	Tag Free
Symbol:	HMGCS1
Synonyms:	HMGCS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC316235 representing NM_001098272.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCCTGGATCACTTCCTTTGAATGCAGAAGCTTGCTGGCCAAAAGATGTGGGAATTGTTGCCCTTGAG
ATCTATTTTCCTTCTCAATATGTTGATCAAGCAGAGTTGGAAAAATATGATGGTGTAGATGCTGGAAG
TATACCATTGGCTTGGGCCAGGCCAAGATGGGCTTCTGCACAGATAGAGAAGATATTAACCTCTTTGC
ATGACTGTGGTTCAGAATCTTATGGAGAGAAAACCTTTCTATGATTGCATTGGGCGGCTGGAAGTT
GGAACAGAGACAATCATCGACAAATCAAAGTCTGTGAAGACTAATTTGATGCAGCTGTTGAAGAGTCT
GGGAATACAGATATAGAAGGAATCGACACAATAATGCATGCTATGGAGGCACAGCTGCTGTCTTCAAT
GCTGTTAACTGGATTGAGTCCAGCTCTTGGATGGACGGTATGCCCTGGTAGTTGCAGGAGATATTGCT
GTATATGCCACAGGAAATGCTAGACCTACAGGTGGAGTTGGAGCAGTAGCTCTGTAATTGGCCAAAT
GCTCCTTAATTTTTGAACGAGGGCTTCGTGGACACATATGCAACATGCCTATGATTTTTACAAGCCT
GATATGCTATCTGAATATCTATAGTAGATGGAAAACCTCCATACAGTGTACCTCAGTGCATTAGAC
CGCTGCTATTCTGTCTACTGCAAAAAGATCCATGCCCAGTGGCAGAAAAGAGGGAAATGATAAAGATTTT
ACCTTGAATGATTTTGGCTTCATGATCTTCACTCACCATATTGTAACCTGGTTCAGAAATCTCTAGCT
CGGATGTTGCTGAATGACTTCCTTAATGACCAGAATAGAGATAAAAAATAGTATCTATAGTGGCCTGGAA
GCCTTTGGGGATGTTAAATTAGAAGACCTACTTTGATAGAGATGTGGAGAAGGCATTTATGAAGGCT
AGCTCTGAACTCTTCAGTCAGAAAACAAAGGCATCTTTACTTGTATCAAATCAAATGGAATATGTAC
ACATCTTCAGTATATGGTCCCTTGCATCTGTTCTAGCACAGTACTCACCTCAGCAATTAGCAGGGAAG
AGAATTGGAGTGTCTTCTATGGTCTGGTTGGCTGCCACTCTGTACTCTTTAAAGTCACACAAGAT
GCTACACCGGGTCTGCTCTTGATAAAATAACAGCAAGTTTATGTATCTTAAATCAAGGCTTGATTCA
AGAACTGGTGTGGCACCAGATGCTTCGCTGAAAACATGAAGCTCAGAGAGGACACCCATCATTTGGTC
AACTATATTTCCCAAGGTTCAATAGATTCACTCTTTGAAGGAACGTGGTACTTAGTTAGGTGGATGAA
AAGCACAGAAGAATTACGCTCGGCGTCCCACTCCAAATGATGACACTTTGGATGAAGGAGTAGGACTT
GTGCATTCAAACATAGCAACTGAGCATATTCGAAGCCCTGCCAAGAAAAGTACCAAGACTCCCTGCCACA
GCAGCAGAACCTGAAGCAGCTGTCATTAGTAATGGGGAACATTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:

NM_001098272

Insert Size:	1563 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001098272.2
RefSeq Size:	5450 bp
RefSeq ORF:	1563 bp
Locus ID:	3157
UniProt ID:	Q01581
Cytogenetics:	5p12
Protein Families:	Druggable Genome
Protein Pathways:	Butanoate metabolism, Metabolic pathways, Synthesis and degradation of ketone bodies, Terpenoid backbone biosynthesis, Valine, leucine and isoleucine degradation
MW:	57.3 kDa
Gene Summary:	This enzyme condenses acetyl-CoA with acetoacetyl-CoA to form HMG-CoA, which is the substrate for HMG-CoA reductase.[UniProtKB/Swiss-Prot Function]