

## Product datasheet for **SC316239**

### zinc finger protein 655 (ZNF655) (NM\_001085368) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	zinc finger protein 655 (ZNF655) (NM_001085368) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF655
Synonyms:	VIK; VIK-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC316239 representing NM\_001085368.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGGAAATACCAGCCAGGAAGCAGCAGGGTCACCAAGGTCAGTTTCAGTCTTTGGAGACCCAG
CTGAGTGTCTGTCCCCAGAGCCTCAGTTTGTGCAGGACACCCGACATGGAACAGGGACTCACTGGGGGA
TTTCCAATTTCCAAGCCTGATGGGATCTCCAGCTGGAACAGGATCTACAGGTCTTTGATCTGAAACT
AAGACTAGAGAAGTCTTAAGAGATGACTGCTCAGATGGAGAGACCAGAGAAGAGAACAAGCTGTTGATT
CCTAAGCAGAAAAATTCGGAAGAAGTGCATTCATACAAAGTGAGAGTAGGAAGACTCAAACAGATATT
ACCCAAGTTCCTGAGACTAGAGAAGTGTATAAGTCTGAGGACAGATTAGAAAGACTTCAGGAAATCTA
AGGAAATTTCTGTACCTGGAGAGAGAGTTTAGGCAAATAACAATCAGCAAGGAAACCTTCACCAGTGAG
AAGAACAATGAATGTCATGAACCCGAAAAAAGCTTCAGTCTGGACTCTACTATTGATGCAGATCAGAGA
GTTCTTAGAATACAGAATACCGATGACAATGATAAGTATGACATGAGCTTCAACCAGAATTCAGCCTCT
GGTAAACATGAACACTTAAATCTAACAGAGGATTTTCAGAGTAGTAATGTAAGGAAAGCTTAATGGAT
CTCTCCACCTTAATAAATGGGAGAGCATCCCTAACACTGAGAAATCCTATAAATGTGATGTATGTGGG
AAAATTTCCATCAGAGCTCAGCCCTTACTAGACATCAGAGAATCCATACTAGAGAGAAGCCCTACAAA
TGTAAGAATGTGAAAAGTCTTTCAGTCAGAGCTCAAGTCTTAGTCGACATAAAAAGAATACACACTAGA
GAAAAACCTTACAAATGTGAAGCATCTGATAAATCCTGTGAAGCGTCTGATAAATCCTGTAGTCCAAGC
TCAGGCATAATTCAGCATAAGAAAATTCACACCAGAGCCAAATCTTACAAATGTAGCAGTTGTGAAAAGA
GTCTTCAGTCGTAGTGTCCACCTTACTCAACATCAGAAAATTCACAAAGAGATGCCCTGTAAGTGTACT
GTATGTGGCAGTGACTTCTGCCATACTTCATACCTACTTGAACATCAGAGGGTCCATCATGAAGAGAAA
GCCTATGAGTATGATGAATATGGGTTGGCCTATATTAACAACAAGGAATTCATTTAGAGAAAAGCCC
TATACGTGTAGTGAATGTGGAAGAAAGACTTCAGATTGAATTCACATCTTATTCAGCATCAAAGAATTCAC
ACAGGAGAGAAAAGCACATGAATGTAATGAATGTGGAAGGCTTTTCAGTCAAACCTCATGCCTTATTCAG
CATCACAAAATGCATAGGAAAAGAGAAATCGTATGAATGTAATGAGTATGAGGGCAGTTTCAGTCATAGC
TCAGATCTTATCCTGCAACAAGAAGTCTCACCAGACAGAAAAGCCTTTGATTGTGATGTATGGGAAAAG
AACTCCAGTCAGAGAGCACATCTAGTTCAACATCAGAGCATTTCATACCAAAGAGAAGCTATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

Sgfl-Mlul

**Plasmid Map:**



**ACCN:**

NM\_001085368

<b>Insert Size:</b>	1581 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001085368.1</a></u>
<b>RefSeq Size:</b>	4690 bp
<b>RefSeq ORF:</b>	1581 bp
<b>Locus ID:</b>	79027
<b>UniProt ID:</b>	<u><a href="#">Q8N720</a></u>
<b>Cytogenetics:</b>	7q22.1
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	61.4 kDa
<b>Gene Summary:</b>	<p>This gene encodes a zinc finger protein. The zinc finger proteins are involved in DNA binding and protein-protein interactions. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (8) differs in the 5' UTR compared to variant 7. Variants 7 and 8 encode the longest isoform (f).</p>