

Product datasheet for **SC316249**

PNMA8B (NM_020709) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PNMA8B (NM_020709) Human Untagged Clone
Tag:	Tag Free
Symbol:	PNMA8B
Synonyms:	PNMAL2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC316249 representing NM_020709.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGGCCATGAGCCTTTTGCAGGACTGGTGCCGGAGCCTGGACGTGGACGCGCACAGGGCCCTGCTGGTC
ACCGGCATCCCGAGGGCCTGGAGCAGGACAGCTCGAAGCCGCTCTGCAGCCGACCCCTCTGCCCTG
GGCAGTTTCAGGTTGCGACACATGAAGGCTTTGATGAACGAGAAGGCCAGGCCGCTTGGTGGAGTTT
GTGGAGGACGTCAATCACGCTGCCATTCCAGGGAGATCCAGGCAAGGATGGGGTCTGGAGGGTCTG
TGAAGGACCGTGCAGGACACGAGGGTCTGAGGCAGATGAGACGCTGCTGCTGGATGACGGGCC
ACGAGGCCGCGGAGGCTGGGACCCCGGGGAGGCACCCACCCCTCCCGCTTCGGAGACGAGGCCAG
GATTCTGGGGAGGTAACAGGCGAGGCTGGCTCGCTTCTGGGGCAGCCAGGAACCAAGGAGGGCCGT
CGGGTTCGAGAAACAGAACGACGCAACAGGTTGACCCAGAAGGGCAAGAAGAGAAGCCGAGGAGGA
CGGCCGTCTGCTCCGCGAGGAGTGAGGCCGAGGACTTCCGACGAGAGCTGGGCATCGTGATCGAG
GAGATCGACAGGGCGACCTGAGCGGAGAAGAGGACCAGAGCGCGCTGTACGCCAGCTGCAGGCCGCT
GCCAGGGAGCTGGTTAGGCAGTGGGCGCCCTGCAACTCCAGGGGGAAGAAGACGGTCCCGCGAGTTC
TTGGCTCTGGTCACCGTCACCGACAAATCGAAGAAAGAAGAGGCAGAGAAGGAGCCAGCTGGGGCCGAA
TCCATCCGCTTGAACACCAAGAAGACAAAAATGGTGTCCCGACTTAGTGGCCCTGCTGGCTGTGAGA
GACACCCCGGACGAGGAGCCGCTGGACAGCGACACTTCGAGAGCGACTCGCAGGAAAGTGGGGACCA
GAAACAGAGGAGTTGGATAATCTGAGTTCTGGCCATTGTGGCTATACCGACCCGTGGACCCCTGG
GCCCGGGAGGAGATGTTGAAAATCGCTTCTGTTATCGAGTCGCTGGGCTGGAGCGACGAGAAAGACAAG
CGAGACCCCTCCGACAGGTCTTGTCCGTCATGTCCAAGGACACTAACGGGACCCGCGTGAAGGTGAA
GAGGCGGGCCGCGAGGTGGACGCCGTGGTCTGCGCAAGGCCGGGATGACGGGGACCTCCGGGAGTGC
ATTTCCACCTTGGCGCAGCCGGATCTCCCTCCCGAGGCGAAGAAGGCTGGGCGTGGCCTCTTCGGGGGC
TGGAGCGAGCACCGTGAGGACGAAGGGGTCTTCTGGAGCTGGTGGCGCTCTGGCTGCCAGGACATG
GCGGAGGTGATGAAGGAGGAAAAAGAAACGCCTGGGAAGGCGGGAAGTACAAATACCCAAAGGCAAA
CTGGGGGAGGTATTGGCGCTCTGGCCGCCCGGGAACATGGGGTCCAACGAGGGGTGGAGGAGGCT
TCGACGAAACAGTCCGAGGAGGAGTGGAGGACACCGAGAGCGAGGCGTGGAGCCGAGGACAGGGCA
TCCAGGAAGCCCGGGCAAGAGGGCGCGACGGCCCCAGGGGCTGACTCCGGCCGGCGCGCTCCC
ACCGCTTCCGGGGCCCGAAAACCCGCGGGCGGCCGAGGCCGAGGCCGGGGCGTCACTCCGAGAAG
AAAGCCGGGAGCCGGGGCTCGGCCAGGACGACCGCGCAGGAAGCAGGAAGAAGAAGGGGAGCGCGGC
GCCGGGGCCCATGCCAGGCGAGGCGAGGCCAAGGGCCAGGCGCCACTGGATCCAAGGCCGCGCGGG
AAGAAGGCCCGTGGGGCCGGAGGCTGCCCCCTAAATGCCGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_020709

Insert Size: 1908 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_020709.2](#)

RefSeq Size: 4739 bp

RefSeq ORF: 1908 bp

Locus ID: 57469

UniProt ID: [Q9ULN7](#)

Cytogenetics: 19q13.32

MW: 68.6 kDa