

Product datasheet for **SC316266**

MYSM1 (NM_001085487) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MYSM1 (NM_001085487) Human Untagged Clone
Tag:	Tag Free
Symbol:	MYSM1
Synonyms:	2A-DUB; 2ADUB; BMFS4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC316266 representing NM_001085487.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGGCTGAAGAGGCGGATGTGGATATCGAAGGGGACGTGGTAGCGGGCGGGGGCCACAGCCAGGA
AGTGGTGAAAATACAGCATCAGTTTTACAAAAAGATCACTATCTTGATTCATCTTGAGAACAGAGAAT
GGCCTTATTCTTGGACCTTGGATAACACCATCAGTGAAGAGAACAGAGCTGTTATTGAGAAAATGTTG
TTGAAGAAGAATATTATTTATCTAAAAAATCACAAACCGAAAAAGTCTGGCTTGATCAAAAAGGAAGAT
GATAAAAAATACATGAAGAGTCTGCAGAAAAACAGCAAAAATCATGGTACACTCTCCTACAAAACCAGCC
AGTTACTCAGTAAAGTGGACGATAGAAGAAAAAGAGCTGTTTGAACAAGGGCTGGCTAAATTTGGCCGA
AGATGGACCAAAATTTCAAAGCTAATTGGAAGCCGCACTGTTTTACAAGTGAAGAGTTATGCAAGACAG
TATTTTAAAAATAAGGTCAAATGCGGTCTGGATAAAGAAAACCAAAATCAGAAGACCGGCCATAATCTT
CAAGTTAAAAATGAAGATAAAGGGACAAAGGCATGGACACCATCATGTTTAAAGGGGACGTCTGATCCC
AACTTGAATGCTGTAATAATGAAAAGTTATCTGATGATGAAGAAGTAGACATCACAGATGAGGTGGAC
GAGTTGCTTCTCAAACACCCCGAGAAGAATTCTAGCAGTGATCTCTTGTAGACTTTCCTAATAGTAAA
ATGCATGAAACCAATCAAGGAGAATTCATTACTTCTGACAGCCAGGAAGCTCTCTTTTCTAAGTCTTCC
AGGGGCTGTCTTCAAATGAAAAGCAAGATGAAACACTTTCAAGCTCAGAAAATTACACTGTGGACTGAG
AAACAGAGCAATGGTGACAAAAATCAATTGAATTAATGACCAGAAAATTAATGAATTGATTAATAAAC
TGCAACAAGCATGATGGAAGGGGAATAATAGTTGATGCCAGGCAGTTGCCTTCTCCAGAGCCTTGTGAA
ATTCAGAAAAATTTGAATGATAATGAAATGCTTTTTTATTCTTGCCTTGGTGGTGGAGGAAAGCCATGAG
GAAGAAGAGCTTAAGCCACCAGAACAGGAAATAGAAAATAGATAGAAAATATCATTCAAGAAGAAAAAAA
CAAGCAATTCCTGAGTTTTTTGAGGGGCCCAAGCTAAAACACCCAGAACGCTATTTGAAAATAGAAAAT
TATATTTTGGATCAATGGGAGATATGCAAACCAAAATACTTAAATAAGACCTCAGTACGTCTGGCCTG
AAGAAGCTGTGGAGATGTTAATTGTATTGGACGGATTATACATACCTCGAATTGATAGGACAAATCAAT
TTTGGATGTGAACAGGCTGTGTATAATAGGCCACAAAACAGTTGACAAAAGTACGAATCAGAGACAGAAAA
GATGCAGTAGAAGCATACCAACTGCCAGCGTCTGCAGTCTATGCGTACAAGGAGAGCTAGGGTCCGA
GACCCATGGGAAACTGGTGTGATGCAAAGGACTTAGAAGGACAAAACGTTTGGCATCTCTCTGCTGAG
GAGTTGGCAAAAAGAAGAGAAGAGGAAAAAGGCAGACCTGTTAAATCTTAAAAGTGCCAAGACCAACA
AAAAGCTCGTTTATCCCTTCCAAGTACCTTGTAAATTTTTTAAAGTGAAGAAAAGCAGGAGCCATTT
CAGGTGAAAGTGGCTTCAAGCACTTTTAAATAATGGATTTGCATGCTCATGTTTCTATGGCAGAAGTG
ATTGGTCTGTTAGGAGGAAGATACTCAGAAGTTGATAAAGTAGTTGAAGTCTGTGCAGCAGAACCATGT
AACAGTCTGAGTACAGGACTACAGTGTGAGATGGATCCTGTATCACAACACAGGCCTCAGAAAACCTTG
GCTGTTAGAGGCTTCAAGTGTATTGGATGGTATCATTCTCATCCTGCTTTTGTATCCTAATCCTTCTTA
CGAGATATTGACACACAAGCTAAATACCAGAGTACTTCTCCAGAGGAGGTGCAAAGTTCATTGGGATG
ATTGTTAGTCCCTATAATCGAAATAATCCCTTACCATATTCTCAGATTACCTGCCTGTTATAAGTGAG
GAAATTAGCCAGATGGCTTTATCGCTTACCTTACAAATTTGAAGTACAGCAGATGTTAGAAGAACCT
CAGTGGGATTAGTATTTGAAAAGACAAGATGGATAATAGAAAAATACAGGCTCTCCCATAGCAGCGTC
CCATGGATAAAAATCTTTCGCCGGGATTCTGACCTGACTTGTGTTGCAGAAAATTTGGAGTGTATGAGG
AAGACTCTGAGCAAAGTGACCAATTGCTTATGGCTGAAGAATCTTGACTGAAATAGAAAATTTGTTT
CTTTCCAATTATAAAAAGCAACCAAGAGAATGGAGTAACCGAAGAGAAGTGTACAAAGGAATTGTTAATG
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001085487

Insert Size: 2487 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001085487.2](#)

RefSeq Size: 7768 bp

RefSeq ORF: 2487 bp

Locus ID: 114803

UniProt ID: [Q5VVJ2](#)

Cytogenetics: 1p32.1

Protein Families: Protease

MW: 95 kDa

Gene Summary: Metalloprotease that specifically deubiquitinates monoubiquitinated histone H2A, a specific tag for epigenetic transcriptional repression, thereby acting as a coactivator. Preferentially deubiquitinates monoubiquitinated H2A in hyperacetylated nucleosomes. Deubiquitination of histone H2A leads to facilitate the phosphorylation and dissociation of histone H1 from the nucleosome. Acts as a coactivator by participating in the initiation and elongation steps of androgen receptor (AR)-induced gene activation. Required for correct regulation of hematopoiesis and lymphocyte differentiation (PubMed:28115216, PubMed:26220525). [UniProtKB/Swiss-Prot Function]