

Product datasheet for **SC316268**

CTNND1 (NM_001085467) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CTNND1 (NM_001085467) Human Untagged Clone
Tag:	Tag Free
Symbol:	CTNND1
Synonyms:	BCDS2; CAS; CTNND; p120; p120(CAS); p120(CTN); P120CAS; P120CTN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC316268 representing NM_001085467.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCAGGAGCCGGGCGAGATTGTGGAGACCTACACGGAGGAGGATCCTGAGGGAGCCATGTCTGTAGTC
TCTGTGGAGACCTCAGATGATGGGACCACTCGGCGCACAGAGACCACGGTCAAGAAAGTAGTGAAGACT
GTGACAACACGGACAGTACAGCCAGTCGCTATGGGACCAGACGGTTGCCTGTGGATGCTTCATCAGTT
TCTAACAACTATATCCAGACTTTGGGTCGTGATTTCCGCAAGAATGGCAATGGGGAGCTGGTCCCTAT
GTGGGGCAAGCTGGCACTGCTACCCTTCTAGGAACCTTCCACTACCCTCCTGATGGTTATAGTCGCCAC
TATGAAGATGGTTATCCAGGTGGCAGTGATAACTATGGCAGTCTGTCCCGGGTGACCCGATTGAGGAG
CGGTATAGGCCACAGATGGAAGGCTACCGGGCACCTAGTAGACAGGATGTGTATGGGCCCAACCCAG
GTTCCGGTAGGTGGGAGCAGCGTGGATCTGCATCGCTTTCATCCAGAGCCTTATGGGCTAGAGGATGAC
CAGCGTAGTATGGGCTATGATGACCTGGATTATGGTATGATGTCTGATTATGGCACTGCCGTGGACT
GGGACACCCTCTGACCCTCGTCGGCGCCTCAGGAGCTATGAAGACATGATTGGTGAGGAGTGCCATCG
GATCAATACTACTGGGCTCCTTTGGCCACAGATGAGCGAGGAAGTTAGCAAGCTTGGATAGCCTGCCG
AAAGGAGGGCCTCCACCTCCTAATTGGAGACAGCCAGAGCTGCCAGAGGTGATCGCCATGCTTGGATTC
CGCTTGGATGCTGTCAAGTCCAATGCAGCTGCATACCTGCAACACTTATGCTACCGAATGACAAGGTG
AAGACTGACGTGCGGAAGCTCAAGGGCATCCCAGTACTGGTGGGATTGTTAGACCATCCCAAAAAGGAA
GTGCACCTTGGAGCCTGTGGAGCTCTCAAGAATATCTCTTTGGACGTGACCAGGATAACAAGATTGCC
ATAAAAACTGTGATGGTGTGCCCTGGCTTGTGCGATTGCTTCAAAGGCTCGTGATATGGACCTTACT
GAAGTTATTACCGAACCCTGTGGAATCTTTCATCCCATGACTCAATCAAATGGAGATTGTGGACCAT
GCACTGCATGCCCTTGACAGATGAAGTGATCATTCTCATTCTGGTTGGGAGCGGGAACCTAATGAAGAC
TGTAAAGCCACGCCATATTGAGTGGGAATCGGTGCTCACCAACACAGCTGGCTGCCTTAGGAATGTAAGC
TCAGAGAGGAGTGAAGCTCGCCGAAACTCGGGAATGTGATGGTTTAGTTGATGCCCTCATTTCATT
GTTCAGGCTGAGATTGGGCAGAAGGATTAGACAGCAAGCTTGTAGAGAAGTGTGTTGCTTCTTCGG
AACTTATCATATCAAGTTCACCGGGAGATCCCACAGGCAGAGCGTTACCAAGAGGCAGCTCCCAATGTT
GCCAACAACTACTGGCCACATGCTGCCAGTTGCTTTGGGGCAAGAAGGGCAAAGGGAAAAAACCTATA
GAGGATCCAGCAAACGATACAGTGGATTTCCCTAAAAGAACGAGTCCAGCTCGAGGCTATGAGCTCTTA
TTTCAGCCAGAGGTGGTTCGGATATACATCTCACTTCTTAAGGAGAGCAAGACTCCTGCCATCCTAGAA
GCCTCAGCTGGAGCTATCCAGAACTTGTGTGCTGGGCGCTGGACGTATGGTCGATACATCCGCTCTGCT
CTGCGTCAAGAGAAGGCTCTTCTGCCATAGCTGACCTCCTGACTAATGAACATGAACGGGTGGTAAA
GCTGCATCTGGAGCACTGAGAACTGGCTGTGGATGCTCGCAACAAAGAATTAATTGGTAAACATGCT
ATTCTAACTTGGTAAAGAATCTGCCAGGAGGACAGCAGAACCTCCTTGGAAATTTCTCTGAGGACACT
GTCATCTCTATTTTGAACACTATCAACGAGGTTATCGCTGAGAAGTGGAGGCTGCCAAAAAGCTTCGA
GAGACACAGGGTATTGAGAAGCTGGTGTGATCAACAAATCAGGGAACCGCTCAGAAAAAGAAGTTCGA
GCAGCAGCACTTGTATTACAGACAATCTGGGGATATAAGGAATCGCGAAGCCACTGGAAAAAGAAGGA
TGAAGAAATCAGACTTTCAGGTGAATCTAAACAATGCTTCCCGAAGCCAGAGCAGTCATTTCATATGAT
GATAGTACTCTCCCTCTCATTGACCGGAACCAAAAATCAGATAAGAAACCTGATCGGGAAGAAATTCAG
ATGAGCAATATGGGATCAAAACAAAAATCACTAGATAACAATTTCCACACCAATGAGAGAGGAGAC
CACAATAGAACACTGGATCGATCGGGGATCTAGGCGACATGGAGCCATTGAAGGGAACAACACCCTTG
ATGCAGAAGATTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

Plasmid Map:


ACCN: NM_001085467

Insert Size: 2499 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001085467.1](#)

RefSeq Size: 5850 bp

RefSeq ORF: 2499 bp

Locus ID: 1500

UniProt ID: [O60716](#)

Cytogenetics: 11q12.1

Protein Families:	Druggable Genome
Protein Pathways:	Adherens junction, Leukocyte transendothelial migration
MW:	92.7 kDa
Gene Summary:	<p>This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length natures of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene. [provided by RefSeq, Dec 2010]</p> <p>Transcript Variant: This variant (11) lacks multiple alternate segments in the 5' UTR and in the coding region, compared to variant 1. This variant shares the alternate name of hp120ctn-3A with variants 12 and 13. Variants 11, 12, and 13 encode the same protein (isoform 3A), which is shorter than isoform 1ABC.</p>