

## Product datasheet for **SC316272**

### CTNND1 (NM\_001085464) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CTNND1 (NM_001085464) Human Untagged Clone
Tag:	Tag Free
Symbol:	CTNND1
Synonyms:	BCDS2; CAS; CTNND; p120; p120(CAS); p120(CTN); P120CAS; P120CTN
Vector:	<u>pCMV6 series</u>



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001085464, the custom clone sequence may differ by one or more nucleotides

```

ATGCAGGAGCCGGGGCAGATTGTGGAGACCTACACGGAGGAGGATCCTGAGGGAGCCATG
TCTGTAGTCTCTGTGGAGACCTCAGATGATGGGACCACTCGGCGCACAGAGACCACGGTC
AAGAAAGTAGTGAAGACTGTGACAACACGGACAGTACAGCCAGTCGCTATGGGACCAGAC
GGGTTGCCTGTGGATGCTTCATCAGTTTCTAACAACTATATCCAGACTTTGGGTCGTGAT
TTCCGCAAGAATGGCAATGGGGACCTGGTCCCTATGTGGGGCAAGCTGGCACTGTACC
CTTCCTAGGAACCTCCACTACCCTCCTGATGGTTATAGTCGCCACTATGAAGATGGTTAT
CCAGGTGGCAGTGATAACTATGGCAGTCTGTCCCGGGTGACCCGATTGAGGAGCGGTAT
AGGCCCAGCATGGAAGGCTACCGGGCACCTAGTAGACAGGATGTGTATGGGCCCAACCC
CAGGTTCCGGTAGGTGGGAGCAGCGTGGATCTGCATCGCTTTCATCCAGAGCCTTATGGG
CTAGAGGATGACCAGCGTAGTATGGGCTATGATGACCTGGATTATGGTATGATGTCTGAT
TATGGCACTGCCCGTCCGACTGGGACACCCTCTGACCCTCGTCGGCGCCTCAGGAGCTAT
GAAGACATGATTGGTGGAGAGTGCCATCGGATCAATACTACTGGGCTCCTTTGGCCAG
CATGAGCGAGGAAGTTTAGCAAGCTTGGATAGCCTGCGCAAAGGAGGGCCTCCACCTCCT
AATTGGAGACAGCCAGAGCTGCCAGAGGTGATCGCCATGCTTGGATTCCGCTTGGATGCT
GTCAAGTCCAATGCAGCTGCATACCTGCAACACTTATGCTACCGCAATGACAAGGTGAAG
ACTGACGTGCGGAAGCTCAAGGGCATCCCAGTACTGGTGGGATTGTTAGACCATCCCAAA
AAGGAAGTGCACCTTGGAGCCTGTGGAGCTCTCAAGAATATCTCTTTGGACGTGACCAG
GATAACAAGATTGCCATAAAAACTGTGATGGTGTGCCTGCCCTTGTGCGATTGCTTCGA
AAGGCTCGTGATATGGACCTTACTGAAGTATTACCGAACCCCTGTGGAATCTTTCATCC
CATGACTCAATCAAAAATGGAGATTGTGGACCATGCACTGCATGCCTTGACAGATGAAGTG
ATCATTCTCATTCTGGTTGGGAGCGGGAACCTAATGAAGACTGTAAGCCACGCCAATT
GAGTGGGAATCGGTGCTACCAACACAGCTGGCTGCCTTAGGAATGTAAGCTCAGAGAGG
AGTGAAGCTCGCCGAAACTTCGGGAATGTGATGGTTAGTTGATGCCCTCATTTTCATT
GTTCAGGCTGAGATTGGGCAGAAGGATTCAGACAGCAAGCTTGTAGAGAAGTGTGTTGC
CTTCTTCGGAACTTATCATATCAAGTTCACCGGGAGATCCACAGGCAGAGCGTTACCAA
GAGGCAGCTCCCAATGTTGCCAACAACTGGGCCACATGCTGCCAGTTGCTTTGGGGCC
AAGAAGGGCAAAGGGAAAAAACCTATAGAGGATCCAGCAAACGATACAGTGGATTCCCT
AAAAGAACGAGTCCAGCTCGAGGCTATGAGCTCTATTTTCAGCCAGAGGTGGTTCGGATA
TACATCTCACTTCTTAAGGAGAGCAAGACTCCTGCCATCTAGAAGCCTCAGCTGGAGCT
ATCCAGAACTTGTGTGCTGGGCGCTGGAGCTATGGTCGATACATCCGCTCTGCTCTGCGT
CAAGAGAAGGCTCTTCTGCCATAGCTGACCTCCTGACTAATGAACATGAACGGGTGGT
AAAGCTGCATCTGGAGCACTGAGAAACCTGGCTGTGGATGCTCGCAACAAAGAATTAATT
GGTAAACATGCTATTCCTAACTTGGTAAAGAACTGCCAGGAGGACAGCAGAAGCTCCTCT
TGGAATTTCTCTGAGGACACTGCATCTCTATTTTGAACACTATCAACGAGGTTATCGCT
GAGAAGCTGGAGGCTGCCAAAAAGCTTCGAGAGACACAGGGTATTGAGAAGCTGGTGTG
ATCAACAAATCAGGGAACCGCTCAGAAAAAGAAGTTCGAGCAGCAGCACTTGATTACAG
ACAATCTGGGGATATAAGGAACCTCGGGAAGCCACTGAAAAAGAAGGATGGAAGAAATCA
GACTTTCAGGTGAATCTAAACAATGCTTCCCGAAGCCAGAGCAGTCATTATATGATGAT
AGTACTCTCCCTCTCATTGACCGGAACCAAAAATCAGATAAAGAACTGATCGGGAAGAA
ATTGAGATGAGCAATATGGGATCAAACACAAAATCACTAGATAACAACACTTCCACACCA
AATGAGAGAGGAGACCACAATAGAACACTGGATCGATCGGGGGATCTAGGCGACATGGAG
CCATTGAAGGGAACAACCCCTTGGATGCAGGACGAGGGGCAGGAATCTCTGGAGGAAGAG
TTGGATGTGTTGGTTTTGGATGATGAGGGGGGCCAAGTGTCTTACCCTCCATGCAGAAG
ATT
    
```

**Restriction Sites:** Please inquire  
**ACCN:** NM\_001085464

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001085464.1</a></u> , <u><a href="#">NP_001078933.1</a></u>
<b>RefSeq Size:</b>	5864 bp
<b>RefSeq ORF:</b>	2586 bp
<b>Locus ID:</b>	1500
<b>UniProt ID:</b>	<u><a href="#">O60716</a></u>
<b>Cytogenetics:</b>	11q12.1
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Adherens junction, Leukocyte transendothelial migration
<b>Gene Summary:</b>	<p>This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length natures of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene. [provided by RefSeq, Dec 2010]</p> <p>Transcript Variant: This variant (8) has multiple differences and uses a downstream start codon, compared to variant 1. The resulting protein (isoform 3AB) has a shorter N-terminus, compared to isoform 1ABC. This variant is also called hp120ctn-3AB.</p>