

Product datasheet for **SC316284**

CTNND1 (NM_001085458) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CTNND1 (NM_001085458) Human Untagged Clone
Tag:	Tag Free
Symbol:	CTNND1
Synonyms:	BCDS2; CAS; CTNND; p120; p120(CAS); p120(CTN); P120CAS; P120CTN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC316284 representing NM_001085458.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGACGACTCAGAGGTGGAGTCGACCGCCAGCATCTTGGCCTCTGTGAAGGAACAAGAGGCCAGTTT
GAGAAGCTGACCCGGGCGCTGGAGGAGAACGGCGCCACGCTCGGCGCAGCTGGAACCGTCCGGGTC
TACCACAAGATGCCAACCCACTCATGGCCAACGGCACACTCACCCGCCGCATCAGAACGGCCGGTTT
GTGGGCGATGCTGACCTTAAAGACAGAAATTTTCAAGTTTAAACTCAACGGACCCAGGATCACAGT
CACCTTCTATATAGCACCATCCCCAGGATGCAGGAGCCGGGGCAGATTGTGGAGACCTACACGGAGGAG
GATCCTGAGGGAGCCATGTCTGTAGTCTCTGTGGAGACCTCAGATGATGGGACCACTCGGCGCACAGAG
ACCACGGTCAAGAAAGTAGTGAAGACTGTGACAACACGGACAGTACAGCCAGTCGCTATGGGACCAGAC
GGGTTGCCTGTGGATGCTTCATCAGTTTCTAACAATAATCCAGACTTTGGTCTGTGATTTCCGCAAG
AATGGCAATGGGGACCTGGTCCCTATGTGGGCAAGCTGGCACTGCTACCTTCTAGGAACTCCAC
TACCCTCTGATGTTATAGTCGCCACTATGAAGATGGTTATCCAGGTGGCAGTGATAACTATGGCAGT
CTGTCCCGGGTGACCCGATTGAGGAGCGGTATAGGCCACAGCATGGAAGGCTACCGGGCACCTAGTAGA
CAGGATGTGTATGGGCCCAACCCAGGTTCCGGTAGGTGGGAGCAGCGTGGATCTGCATCGCTTTCAT
CCAGAGCCTTATGGGCTAGAGGATGACCAGCGTAGTATGGGCTATGATGACCTGGATTATGGTATGATG
TCTGATTATGGCACTGCCCGTCGGACTGGGACACCCTCTGACCCTCGTCGGCGCCTCAGGAGCTATGAA
GACATGATTGGTGAAGGAGTCCCATCGGATCAATACTACTGGGCTCCTTTGGCCAGCATGAGCGAGGA
AGTTTAGCAAGCTTGGATAGCCTGCGCAAAGGAGGGCCTCCACCTCCTAATTGGAGACAGCCAGAGCTG
CCAGAGGTGATCGCCATGCTTGGATTCCGCTTGGATGCTGTCAAGTCCAATGCAGCTGCATACCTGCAA
CACTTATGCTACCGCAATGACAAGGTGAAGACTGACGTGCGGAAGCTCAAGGGCATCCCAGTACCTGGTG
GGATTGTTAGACCATCCCAAAAAGGAAGTGCACCTTGGAGCCTGTGGAGCTCTCAAGAATATCTCTTTT
GGACGTGACCAGGATAACAAGATTGCCATAAAAACTGTGATGGTGTGCCTGCCCTTGTGCGATTGCTT
CGAAAGGCTCGTGATATGGACCTTACTGAAGTTATTACCGGAACCTGTGGAATCTTTCATCCCATGAC
TCAATCAAAATGGAGATTGTGGACCATGCACTGCATGCCTTGACAGATGAAGTGATCATTCTCATTCT
GGTTGGGAGCGGGAACCTAATGAAGACTGTAAGCCACGCCATATTGAGTGGGAATCGGTGCTCACCAAC
ACAGCTGGCTGCCTTAGGAATGAAGCTCAGAGAGGAGTGAAGCTCGCCGAAACTTCGGGAATGTGAT
GGTTTAGTTGATGCCCTCATTTCATTGTTGAGGCTGAGATTGGGAGAAAGGATTCAGACAGCAAGCTT
GTAGAGAACTGTGTTTGCCTTCTCGGAACTTATCATATCAAGTTCACCGGGAGATCCCACAGGCAGAG
CGTTACCAAGAGGCAGCTCCAATGTTGCCAACATACTGGGCCACATGCTGCCAGTTGCTTTGGGGCC
AAGAAGGGCAAAGATGAGTGGTTCTCCAGAGGGAAAAACCTATAGAGGATCCAGCAAACGATACAGTG
GATTTCCCTAAAAGAACGAGTCCAGCTCGAGGCTATGAGCTCTTATTTACGCCAGAGGTGGTTCCGGATA
TACATCTCACTTCTTAAGGAGAGCAAGACTCCTGCCATCCTAGAAGCCTCAGCTGGAGCTATCCAGAAC
TTGTGTCTGGGCGCTGGACGTATGGTCGATACATCCGCTCTGCTCTGCTCAAGAGAAGGCTCTTCT
GCCATAGCTGACCTCCTGACTAATGAACATGAACGGGTGGTGAAGCTGCATCTGGAGCACTGAGAAAC
CTGGCTGTGGATGCTCGCAACAAAGAATTAATTGGTAAACATGCTATTCTAAGTGGTAAAGAATCTG
CCAGGAGGACAGCAGAACTCCTTGGAAATTTCTGAGGACACTGTCATCTCTATTTTGAACACTATC
AACGAGGTTATCGCTGAGAACTTGGAGGCTGCCAAAAAGCTTCGAGAGACACAGGATATTGAGAAGCTG
GTGTTGATCAACAAATCAGGGAACCGCTCAGAAAAAGAAGTTCGAGCAGCAGCACTTGTATTACAGACA
ATCTGGGGATATAAGGAACTGCGGAAGCCACTGGAAAAAGAAGGATGGAAGAAATCAGACTTTCAGGTG
AATCTAAACAATGCTTCCGAAGCCAGAGCAGTCAATCATATGATGATAGTACTCTCCCTCATTGAC
CGGAACCAAAAAATCAGATAAGAAACCTGATCGGGAAGAAATTCAGATGAGCAATATGGGATCAAACACA
AAATCACTAGATAACAATAATCCACACCAATGAGAGAGGAGACCACAATAGAACACTGGATCGATCG
GGGATCTAGGCGACATGGAGCCATTGAAGGGAACAACACCCTTGTATGCAGGACGAGGGGCAAGGAACT
CTGGAGGAAGAGTTGGATGTGTTGGTTTGGATGATGAGGGGGCCAAGTGTCTTACCCCTCCATGCAG
AAGATTTAG
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

Plasmid Map:



ACCN: NM_001085458

Insert Size: 2907 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	NM_001085458.1
RefSeq Size:	6363 bp
RefSeq ORF:	2907 bp
Locus ID:	1500
UniProt ID:	O60716
Cytogenetics:	11q12.1
Protein Families:	Druggable Genome
Protein Pathways:	Adherens junction, Leukocyte transendothelial migration
MW:	108.2 kDa
Gene Summary:	<p>This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length natures of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene. [provided by RefSeq, Dec 2010]</p> <p>Transcript Variant: This variant (1) is the longest variant and it encodes the longest protein (isoform 1ABC). This variant has also been called hp120ctn-1ABC.</p>