


```

AGCAAGATTAGAGAGTTTGAGGTGGAAACAGAACGGCTGAATAGCCGGATTGAGCACCTC
AAATCCCAAAATGACCTCCTGACCATAACCTTGGAGGAATGTAAGCAATGCCGAGAGG
ATGAGCATGCTGGTGGGAAAATACGAATCCAATGCCACAGCGCTGAGGCTGGCCTTGCA
TACAGCGAGCAGTGCAATCGAAGCCTACGAACCTCCTCTGGCGCTGGCAGAGAGTGAGCAG
AGCCTCATCTGGGGCAGTTCGAGCGGGCGGTGGGGTCTCCCCTGGAGACCACTGCG
GGGGATGAAAACATCACTCAGATGCTCAAGCGAGCTCATGACTGCCGGAAGACAGCTGAG
AACGCTGCCAAGGCCCTGCTCATGAAGCTGGACGGCAGCTGTGGGGGAGCCTTTGCCGTG
GCCGGCTGCAGCGTGCAGCCCTGGGAGAGCCTTCTCCAACAGCCACACCAGCACAAAC
AGCTCCACAGCCAGTAGTTGCGACACCGAGTTCCTAAAGAAGACGAGCAGAGGCTGAAG
GATTATATCCAGCAGCTCAAGAATGACAGGGCTGCGGTCAAGCTGACCATGCTGGAGCTG
GAAAGCATCCACATCGATCCTCTCAGCTATGACGTCAAGCCTCGGGGAGACAGCCAGAGG
CTGGATCTGAAAACGCAGTGCTTATGCAGGAGCTCATGGCCATGAAGGAGGAGATGGCC
GAGTTGAAGGCCAGCTCTACCTACTGGAGAAAGAGAAGAAGGCCCTGGAGCTGAAGCTG
AGCACGGGGAGGCCAGGAGCAGGCCTACCTGGTGCACATTGACACCTGAAGTCCGAG
GTGGAGGAGCAGAAGGAGCAGCGGATGCGATCCCTCAGCTCCACCAGCAGCGGCAGCAA
GACAAACCTGGCAAGGAGTGTGCTGATGCTGCCTCCCCAGCTCTGTCCCTAGCCGAACTC
AGGACAACGTGCAGCGAGAATGAGCTGGCTGCGGAGTTACCAACGCCATTTCGTGAGAA
AAGAAGTTGAAGGCCAGAGTTCAAGAGCTGGTGAAGTGCCTTGGAGAGACTACCAAGAGC
AGTGAATCCGACATCAGCAATCTGCAGAGTTCGTGAATGATCTAAAGCGGGCCAACAGC
AACCTGGTGGCTGCCTATGAGAAAGCAAAGAAAAGCATCAAACAACTGAAGAAGTTA
GAGTCGCAGATGATGGCCATGGTGGAGAGACATGAGACCCAAGTGAGGATGCTCAAGCAA
AGAATAGCTCTGCTAGAGGAGGAGAAGTCCAGGCCACACACCAATGAAACTTCGCTT

```

- Restriction Sites:** Please inquire
- ACCN:** NM_001085377
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
- RefSeq:** [NM_001085377.1](#), [NP_001078846.1](#)
- RefSeq Size:** 8679 bp
- RefSeq ORF:** 3060 bp
- Locus ID:** 4163

UniProt ID: [P23508](#)

Cytogenetics: 5q22.2

Protein Families: Druggable Genome

Gene Summary: This gene is a candidate colorectal tumor suppressor gene that is thought to negatively regulate cell cycle progression. The orthologous gene in the mouse expresses a phosphoprotein associated with the plasma membrane and membrane organelles, and overexpression of the mouse protein inhibits entry into S phase. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.