

Product datasheet for **SC316359**

JKAMP (NM_001098625) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	JKAMP (NM_001098625) Human Untagged Clone
Tag:	Tag Free
Symbol:	JKAMP
Synonyms:	C14orf100; C24orf100; CDA06; HSPC213; HSPC327; JAMP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316359 representing NM_001098625. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCATGCCCTGGACTTTATTGTGGGAAGACCCTATTATTTAAAAATGGCTCAACTGAAATATATGGA
GAATGTGGGGTATGCCCAAGAGGACAGAGAACGAATGCACAGAAATATTGTCAGCCTTGACAGAAATCT
CCTGAACCTTATGATTGGCTCTATCTTGGATTTATGGCAATGCTTCCTCTGGTTTTACATTGGTTCTTC
ATTGAATGGTACTCGGGAAAAAGAGTTCCAGCGCACTTTCCAACACATCACTGCATTATTTGAATGC
AGCATGGCAGCTATTATCACCTACTTGTGAGTGATCCAGTTGGTGTCTTTATATTCGTTTCATGTCGA
GTATTGATGCTTCTGACTGGTACACGATGCTTTACAACCAAGTCCAGATTACGTTACCACAGTACAC
TGTAATCATGAAGCCGCTACCCACTATATACCATTGTATTTATCTATTACGCATTCTGCTTGGTATTA
ATGATGCTGCTCCGACCTCTTCTGGTGAAGAAGATTGCATGTGGGTTAGGGAAATCTGATCGATTTAAA
AGTATTTATGCTGCACCTTACTTCTTCCCAATTTAAACCGTGCTTCAGGCAGTTGGTGGAGGCCCTTTTA
TATTACGCCTTCCCATACATTATATTAGTGTTATCTTTGGTACTCTGGCTGTGTACATGTCTGCTTCT
GAAATAGAGAAGTGTATGATCTTCTGGTCAGAAAAGAAAAGACTATTGTTCTTTCAGCCACTGGTTA
CTTCATGCCTATGGAATAATCTCCATTTCCAGAGTGGATAAAGTTGAGCAAGATTGGCCCTTTTGGCT
TTGGTACCTACACCAGCCCTTTTACTTGTTCAGTCAAAAATTTACCGAACCTCAAGGATACTCTCA
GAAGGAGCCAATGGACACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001098625

Insert Size: 918 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001098625.2](#)

RefSeq Size: 2411 bp

RefSeq ORF: 918 bp

Locus ID: 51528

UniProt ID: [Q9P055](#)

Cytogenetics: 14q23.1

Protein Families: Transmembrane

MW: 34.6 kDa

Gene Summary: May be a regulator of the duration of MAPK8 activity in response to various stress stimuli. Facilitates degradation of misfolded endoplasmic reticulum (ER) luminal proteins through the recruitment of components of the proteasome and endoplasmic reticulum-associated degradation (ERAD) system (By similarity).[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) uses an alternate in-frame splice site in the 5' coding region, compared to variant 1. The encoded isoform (2) is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.