

Product datasheet for SC316368

GTF2H2C (NM_001098728) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GTF2H2C (NM_001098728) Human Untagged Clone
Tag:	Tag Free
Symbol:	GTF2H2C
Synonyms:	GTF2H2C_2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316368 representing NM_001098728. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGGATGAAGAACCTGAAAGAACTAAGCGATGGGAAGGAGGCTATGAAAGAACATGGGAGATTCTTAAA
 GAAGATGAATCTGGATCACTTAAAGCTACAATAGAAGACATTCTATTCAAGGCAAAGAGAAAAAGAGTA
 TTTGAGCACCATGGACAAGTTCGACTTGGAAATGATGCGCCACCTTTATGTGGTAGTAGATGGATCAAGA
 ACAATGGAAGACCAAGATTTAAAGCCTAATAGACTGACGTGTACTTTAAAGTTGTTGGAATACTTTGTA
 GAGGAATATTTTATCAAAATCCTATTAGTCAGATTGGAATAATTGTAACCTAAGAGTAAAAGAGCTGAA
 AAATTGACTGAACTTTCAGGAAACCCAAGAAAACATATAACGTCTTTGAAGGAAGCTGTGGATATGACC
 TGCCATGGAGAGCCATCTCTTTATAATTCCTAAGCATGGCTATGCAGACTCTAAACACATGCCTGGA
 CATACAAGTCGAGAAGTACTAATCATCTTTAGCAGCCTTACAACCTGCGATCCATCTAATATTTATGAT
 CTAATCAAGACCCCTAAAGGCAGCTAAAATTAGAGTATCTGTTATTGGATTGTCTGCAGAAGTTCCGCTT
 TGCACTGTACTTGCTCGTGAACTGGTGGCAGTACCATGTTATTTTAGATGAAAGCCATTACAAAGAG
 TTGCTCACACATCATCTTAGTCCTCCTGCTAGCTCAAGTTCTGAATGCTCACTTATTCGTATGGGA
 TTTCTCAGCACACCATTGCTTCTTTATCTGACCAGGATGCAAAACCCTTTTTCAGCATGGCGCATTTG
 GATGGCAATACTGAGCCAGGGCTTACATTAGGAGGCTATTTCTGCCACAGTGTGGGCAAAGTACTGT
 GAGCTACCTGTTGAATGTAATCTGTGGTCTTACTTTGGTGTCTGCTCCCACTTGGCAGGCTTAC
 CATCATTTGTTTCTTTGGATGCTTTTCAAGAAATTCCTAGAGAATAATGGAGAAAAGATTTTGT
 TATGGATGTGAGGGGAATTGAAAGACCAACATGTTTATGTTTGTGCTGTGTGCCAAAATGTTTCTGT
 GTGGACTGTGATTTTTGTTTATGATTCTCTACACTGTTGCCCTGGCTGATTTCATAAGATTCCAGCT
 CCTTCAGGTGTT**TGA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

Restriction Sites: SgfI-MluI


[View online »](#)

ACCN:	NM_001098728
Insert Size:	1188 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001098728.1</u>
RefSeq Size:	2106 bp
RefSeq ORF:	1188 bp
Locus ID:	728340
UniProt ID:	<u>Q6P1K8</u>
Cytogenetics:	5q13.2
Protein Families:	Druggable Genome
MW:	44.5 kDa