

Product datasheet for **SC316405**

CR1L (NM_175710) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CR1L (NM_175710) Human Untagged Clone
Tag:	Tag Free
Symbol:	CR1L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC316405 representing NM_175710.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGCCTCCCGTCCGTCCTCGAGCGTCCCTTCCCTCCCGCGCTTTCCTGGGTTGCTTCTGGCGGCC
CTGGTGTGCTGCTGCTCCTCCTTCCGATCAATGCAATGTCCCGAATGGCTTCCATTTGCCAGGCCT
ACCAACCTAACTGACTTTGAGTTTCCATTGGGACATATCTGAACTATGAATGCCGCCCTGGTTAT
TCCGGAAGACCGTTTTCTATCATCTGCCTAAAAAACTCAGTCTGGACAAGTGCTAAGGACAAGTGCAAA
CGTAAATCATGTCGTAATCCTCCAGATCCTGTGAATGGCATGGCACATGTGATCAAAGACATCCAGTTC
AGATCCCAAATTAATATTTCTTGTCTAAAGGATACCGACTCATTGGTTCCTCGTCTGCCACATGCATC
ATCTCAGGCAACACTGTCATTTGGGATAATAAAACACCTGTTTGTGACAGAATTATTTGTGGGCTACCC
CCCACCATCGCCAATGGAGATTTACTAGCATCAGCAGAGAGTATTTTCACTATGGATCAGTGGTGACC
TACCCTGCAATCTTGAAGCAGAGGGAAAAAGGTGTTTGAGCTTGTGGTGAGCCCTCCATATACTGC
ACCAGCAAAGATGATCAAGTGGGCATCTGGAGTGGCCAGCCCTCAGTGCATTATACCTAACAAATGC
ACGCCTCCAAATGTGGAAAATGGAATATTTGGTATCTGACAACAGAAGCTATTTTCTTAAATGAAGTT
GTGGAGTTTAGGTGTCAGCCTGGCTTTGGCATGAAAGGGCCCTCCCATGTGAAGTGCCAGGCCCTGAAC
AAATGGGAGCCAGAGTTACCAAGCTGCTCCAGGTATGTGAGCCACCTCCAGATGTCTGATGCTGAG
CGTACCCAAAGGGACAAGGACAACTTTTCACCCGGGCAGGAAGTGTCTACAGCTGTGAGCCCGGCTAC
GACCTCAGAGGATCTACGTATTTGCACTGCACACCCAGGGAGACTGGAGCCCTGCAGCCCCCAGATGT
GAAGTGAATCCTGTGATGACTTCTGGCCAACTTCTAATGGCCATGTGCTATTTCCACTTAATCTC
CAGCTTGGAGCAAAGTGGATTTTGGTTTGTGATGAAGGATTTCAATTAAGGCAGCTCTGCTAGTTAC
TGTGTTTTGGCTGGAATGAAAGCCTTTGGAATAGCAGTGTCCAGTGTGTGAACGTAATCATGTGAA
ACTCTCCAGTTCAGTGAATGGCATGGTGCATGTGATCACAGACATCCATGTTGGATCCAGAATCAAC
TATTCTTGTACTACAGGGCACCGACTCATTGGTCACTCATCTGCTGAATGTATCCTCTCGGGCAACTACT
GCCATTGGAGCATGAAGCCACCAATTTGTCAACAAATCTTTTGTCCAAATCCTCCAGCTATCCTTAAT
GGGAGACACACAGGAACTCCCCTTGGAGATATTCCTATGGAAAAGAAGTATCTTACACATGTGACCCC
CACCCAGACAGGGATGACCTTCAACCTCATTGGGGAGACCCATCCGCCGACAAGTGAACCTCAT
GGGAATGGGGTTTGGAGCAGCCCTGCCCTCGCTGTGAACTTCTGTTGGTGCTGGTTCACATGATGCT
CTTATAGTTGGTAAGTTTATGAAGTGTTGCTGAGGAATTCTGTCATCTTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:	NM_175710
Insert Size:	1710 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_175710.1</u>
RefSeq Size:	1837 bp
RefSeq ORF:	1710 bp
Locus ID:	1379
UniProt ID:	<u>Q2VPA4</u>
Cytogenetics:	1q32.2
Protein Families:	Druggable Genome
MW:	62.7 kDa