

Product datasheet for **SC316414**

LRRC43 (NM_001098519) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC43 (NM_001098519) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRC43
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC316414 representing NM_001098519.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAGGCGTCGTACGAGTCCGAGTCCGAGTCCGAGTCTGAGGCCGGCCTGGGACTCAGCGGCCCGGG
ACCGGGACCGTGAGCGCGCCGTCGCGAGCACTTGCAGGAGCTGTGTCTGCGCGAGTTCCCGTGCAGT
GCCGGCAGCTGGAATAAGTCGCGCTTTCTTCTCAAACCTTGCGCAACTTGAGGGAGCTTGCCCCAGA
GAGGAGGATGTGGTGAGCCCCGAGAGGAGACGGTGGAGGCCCTGCTGGGCCTGGTCCGACAGCCGCAC
TCCCCCTGGGCTCTGCTGAACAACCTCGAATGCAGAAGACAGTTTCTGAGAGAATTGGCCATCCGGAAC
CCGCTGACGATCACAGACACCTTCTTCTACTCCTACTTCCGGTCCCTGCGGGTAATAGACAAGAAGGTC
ACCCTGGTGGATAAAGACCTCCTGAAATTTCTAAAGCTGGAGGAGTTGGTACTGAGCGCAATCGAATC
AAGGAGGTGGATGCCACCAATCTGCCCCCACACTCAAGGTGCTGGAGCTCTACGGCAATGAGATCAGC
AGCATGGAGTGTCTGTGTGCCACCCACCCGCCGGCCTGCAGCACTTGGGGTTAGGCCACAACAACCTT
CTAGGCCCTTGAAAGTCTCTACGTACCGCTAATCACTGGCCAACCTCGTCTCCCTGGACCTGGGC
TTCAACGACCTGACAGACCTGCAGAGCATGGTACCAGCCTGAGGACCCTCCGGCACCTGCGACTCCTG
GTGCTGCAGGGAAACCCACTGGCCTTGGTGCCCTACTACCGCGCCTCACCATCGACAGCCTGGCCAG
CTCTGCGTGTGGACGACATCACCGTGTCTCCCAATGAGAAGCATCTTCCGGGGGCTCAGCCTCAAT
GGCGATCTCTTGGCACAGGAGGCGCAGTTTGTGGTGACCATCGGAAACATCAGAGGAGTCTGGACACC
TCTGTCTTAGACCCGGAACCCAGGCCCGAAGGCCCTTTCATCACTTACAACATTTACGTGACCTATGAT
TTTGTGAAAGATGAAGAAGGCGAAATGAATGAGTCCGCGGGCGTCTGGCCGAGATCGTCAAGCCCTCT
CCCAGCTTAGAATTATTAGTTGAGGAATCTCCTGAAGAGTCTGGAAGAGCTCATCGAAGACATTGTT
GAAGAGGTACTGAAGAGGTGGAAGGCTCTGGAGTCTGAGTGGAGGAGTCAAGAGAGTCCGAGCTG
TCTGTCTCGGGCCCTCGACCATCTTGCAGATGCCGAGGGCCTTGCAGAAGAGCTGGCCAAGTTG
AGGCTGCGTATAGATCCCCGGCTCTGCCGTCCCCAGGGACTGTCTTTCAGCACTGCCACAAGCCC
TGGGCTGAGGTACATCCCCTGCAGTTACGAGATGCAGCACTCTCTCAGGGACCTGGTCCCCTGAAGGCC
TTCCTGTGGCGGGGACCACCGTGACCATCGTGGAGGAGAAGATTCTCTCCTGGCCTGTGGTGTACCT
GCTGTTGACAGTCCCCTGTCTGCCAAGAAAGGAAAGGGGAGAAAGACAAGAAAGGAAGGAGAAAGAC
AGGACGGGAAAGGAGAGAAAGAGCCGGCAAGGAGTGAAGGTGCTGAAGAAGAAGAAAGAGCCGCC
AAGGAGTCCGGCAGGACCCCCCATCTCCAGGTGCTGGGCCGGGGCCTGGTGTCTGGAGCCCTG
CTCGCCGGGAGCCCTGGTGTCCACCGTGTGCAACTTCGGCGTGGTCCGCACATTGACATCTGACAGG
CTGACGTTGGCCAGGATTCAAAGAAGATTAAGAAAGTTGCCAAAAAAGAAAAGCCGAAAGCCGTGATT
CCGATCTACGAAGGCGATTACCACCTGAGCCCTGACCGTAGAGGTGCAGATCCAGCTGAACCAAGTGC
CGCTCGGCGGAGGAGGCTCTGCGCATGTTGCCCGTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001098519

Insert Size: 1971 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001098519.1</u>
RefSeq Size:	2045 bp
RefSeq ORF:	1971 bp
Locus ID:	254050
UniProt ID:	<u>Q8N309</u>
Cytogenetics:	12q24.31
MW:	73 kDa