

Product datasheet for **SC316432**

IRAG1 (NM_001098579) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IRAG1 (NM_001098579) Human Untagged Clone
Tag:	Tag Free
Symbol:	IRAG1
Synonyms:	IRAG; JAW1L; MRV11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC316432 representing NM_001098579.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGAATGGATTGACTTGTCTTTTGGCATCTCCCCAGCCTGTGGAGCCCAGGCTCTTGGAGCATC
TTTGGGGCTGACGACGGAGGTTCCGGGCACACGTGGCCACTCCAGCAGGAGGCTGCCATGCCCCAC
ATTCCCGAGGACGAGGAGCCCCCGGAGAGCCACAGGCAGCCAGAGCCCTGCCGCGCAAGTCTCTCT
GCCGACGAGATCTTGCAGTCCAACCTCCACGATTGTCCTGACTGGGGATGCCACTTACCAGAAGGA
GAAACCGACAAAAACCTGGCCAACAGAGTTCACAGTCCCCACAAGAGGCTTTCTACCGACACTTGAAG
GTGTCCACTGCCTCCCTGACATCTGTGGACCCCGGGGCACATCATTGACCTGGTGAATGACCAGCTG
CCAGACATCAGCATCTCAGAGGAGGACAAGAAGAAAAACCTGGCGCTGCTGGAAGAAGCCAAGTTGGTG
AGTGAGCGATTCTGACCCGCCGTGGGAGGAAGTCCAGGAGCAGCCCGGAGACTCCCCATCAGCTGTT
TCCCCGAACCTCAGCCCCAGCGTTCTCCTACATCTCTCGGAGCAACTCACTTACAGTCCCCACCCCG
CCAGGTTTGGATGTGTGCAGTGGCCCGCATCCCCTCTGCTGGAGCACCACCACAGCAGAAGGGGGAT
GAGGCCGACGCTCTTACCTCACCTGGCGAGCCTAACGTCCCCAAAGGGCTAGCTGACAGGAAGCAG
AATGACCAGAGGAAAGTGTCTCAGGGCAGGCTGGCTCCTCGTCTCTCCAGTTGAGAAGTCCAAGAG
ATTGCAATAGAACAAAAGGAAAACCTTCGATCCCCTCCAGTACCCCGAGACCACACCCAAAGGCTAGCT
CCTGTTACAAACAGCAGTGGGAAAATGGCCCTGAACAGCCCTCAGCCTGGCCCGTGGAGAGCGAGCTG
GGGAAGCAGCTCTTGA AACGGGCTGGGAGGGCAGCCCTCTGCCGAGAAGTCCAACCCAGGATGCGGCA
GGAGTGGTCCCCAGCCTCCAGGGGAGAGGCCAGCTGGAGAGCCGATGGGGCCCGAGGCTGGCTCC
AAAGTGAGCTTCCACCCACTGTGTCCCGCCCCCGCTGCTGCGAGGGCTCTCTGGGACAGTGGCCCT
GAAGAACCTGGCCCCGGCTGCAGAAAAGTGTGGCAAGCTGCCACTGGCAGAGGAAGAAAAGTCTTT
GCAGGCAAGGCCCGGGCAAGCTGGCCAAGGCCCTGGTCTCAAAGACTTTCAGATACAAGTGCAGCCC
GTGCGGATGCAGAACTGACCAAGCTCCGAGAGGAGCACATCCTGATGAGAAATCAGAACTTAGTGGGG
CTCAAGCTTCCAGACCTTAGTGAAGCAGCTGAGCAGGAAAAAGGGCTTCTTCTGAACTCTCCCCAGCT
ATTGAGGAAGAAGAGTCAAAGAGTGGCTTAGATGTGCCTAATATTTCTGATGTGCTGCTGCGCAAA
CTGCGGGTCCACAGGAGTCTCCCTGGAAGTGCCCTCCACTCACTGAAAAGGAAGTTGAGAACGTGTTT
GTGCAACTGTCTTGGCCTT TAGAAATGACAGCTACACTCTGGAATCTAGAATTAACCAGGCTGAAAGG
GAACGCAACCTGACAGAGGAGAACACTGAGAAAAGAACTGAAAACCTCAAAGCTTCCATTACGCTCTCA
GCTTCACTCTGGCACCCTGTGAGCACC GGAAACCTACCAGAAGTTGCTGGAGGACATCGCTGTCTTG
CACCGCCTGGCTGCCCGCTCTCCAGCCGAGCTGAGGTGGTAGGGCCGTCCGCCAGGAAAAGCGCATG
TCGAAAGCAACGGAAGTATGATGACAGTATGTGGAGAATCTAAAGAGGACGTATGAGAAGGACCATGCG
GAGCTCATGGAGTTTAAAAAGCTTGCAAAATCAGAATTCAGCCGAGCTGTGGCCCTCTGAAGATGGG
GTCCCTCGCACGGCACGGTCCAATGTCCCTCACGCTGGGAAAGAATATGCCTCGCCGGAGGGTCAGCGTT
GCTGTGGTTCTAAGTTAATGCCCTGAATCTGCCTGGCCAAACTCCCAGCTCATCATCCATTCCCTCC
TTACCAGCCTGTGCGAATCACCAATGGGAAAGGCAGCTACCTGTCACTTACAGCACTGCCTGCACTT
TTGGAAAATGAAAGACAAATGGGGACCCAGATTGTGAAGCCTCTGCTCCTGCGCTGACCCTGAGCTGC
CTGGAGGAGCTTAGTCAGGAGACCAAGGCCAGGATGGAGGAAGAAGCCTACAGCAAGGGATTCCAAGAA
GGTCTAAAGAAGACCAAAGAACTTCAAGACCTGAAGAGGAGGAGGAAGAACAAGAGTGAAGTCTCT
GAGGAACCTGAAGAGGTAGAAGAACTGAGGAAGAGGAAAAGGGCCCAAGAAGCAGCAAACTTGAAGAA
TTGGTCCATTTCTACAAGTCAATGATCCCAAACCTGTGTGAGCACTGGCAAGTACTGGAATGATGGCT
GCAGTGTGCTGGTCTTACTGTTGTGCTGGGGCTCTACAATTCATAACTCTTGTGCAGAGCAGGCT
GATGGGCCCTTGAAGATCCACTTGTCTGCGCAGCCAGAGGGACTCCTGGTGGAGCTCAGGACTCCAG
CATGAGCAGCCTACAGAGCAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001098579

Insert Size: 2715 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001098579.2](#)

RefSeq Size: 6045 bp

RefSeq ORF: 2715 bp

Locus ID: 10335

UniProt ID: [Q9Y6F6](#)

Cytogenetics: 11p15.4

Protein Families: Transmembrane

Protein Pathways: Vascular smooth muscle contraction

MW: 98 kDa

Gene Summary: This gene is similar to a putative mouse tumor suppressor gene (Mrvi1) that is frequently disrupted by mouse AIDS-related virus (MRV). The encoded protein, which is found in the membrane of the endoplasmic reticulum, is similar to Jaw1, a lymphoid-restricted protein whose expression is down-regulated during lymphoid differentiation. This protein is a substrate of cGMP-dependent kinase-1 (PKG1) that can function as a regulator of IP3-induced calcium release. Studies in mouse suggest that MRV integration at Mrvi1 induces myeloid leukemia by altering the expression of a gene important for myeloid cell growth and/or differentiation, and thus this gene may function as a myeloid leukemia tumor suppressor gene. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene, and alternative translation start sites, including a non-AUG (CUG) start site, are used. [provided by RefSeq, May 2011]

Transcript Variant: This variant (1) encodes isoform a, also known as MRVI1A. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.