

Product datasheet for **SC316433**

Epac1 (RAPGEF3) (NM_001098531) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Epac1 (RAPGEF3) (NM_001098531) Human Untagged Clone
Tag:	Tag Free
Symbol:	Epac1
Synonyms:	bcm910; CAMP-GEFI; EPAC; EPAC1; HSU79275
Mammalian Cell Selection:	Neomycin
Vector:	<u>PCMV6-Neo</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001098531, the custom clone sequence may differ by one or more nucleotides

```

ATGAAGTGGGCTGGCCAGGTGAGAGCTGCTGGCAGGTGGGCTGGCTGTGGAGGATAGCCCAGCTCTGG
GAGCACCGCGGGTGGGAGCCCTCCCTGACGTGGTGCCGGAGGGGACACTACTCAACATGGTGTGAGAAG
GATGCACCGGCCCGAAGCTGCTCCTACCAGCTGCTGCTGGAGCACCAGCGTCCGAGCTGCATCCAGGGG
CTGCGCTGGACACCACTACCAACAGCGAGGAGTCCCTGGATTTAGCGAGAGCCTGGAGCAGGCCCTCCA
CAGAGCGGGTGTCTCAGGGCTGGGAGGCAGCTGCATCGGCATCTGCTGGCCACCTGCCAAACCTCATCCG
AGACCGGAAGTACCACCTTAGGCTCTATCGGCAGTGTCTGCTGGCCGGGAGCTGGTGGATGGGATCTTG
GCCCTGGGACTTGGGTCCATTCGCGAGCCAAGTTGTGGGAATCTGCCAGGTGTCTGGATGAAGGTG
CCCTCTGCCATGTGAAACACGACTGGGCCTTCAGGACCGAGATGCCAATTCTACCGGTTCCCGGGCC
CGAGCCCAGCCCGTGGAACTCATGAGATGGAGGAGGATTGGCCGAAGCTGTGGCCCTGCTCTCCAG
CGGGGGCTGACGCCCTGCTCACTGTGGCACTTCGAAAGCCCCAGGTGACGCGCACGGATGAAGAGCTGG
ACCTCATCTTTGAGGAGTGTGCACATCAAGGCTGTGGCCACCTCTCCAACCTCGTGAAGCGAGAATT
AGCGGCTGTTCTGCTCTTTGAACCACACAGCAAGGCAGGACCGTGTGTTGACGCCAGGGGGACAAGGGC
ACTTCGTGGTACATTATCTGGAAGGGATCTGTCAACGTGGTGACCCATGGCAAGGGGCTGGTGACCAACC
TGCATGAGGGAGATGATTTTGGACAGCTGGCTCTGGTGAATGATGACCCCGGGCAGCCACCATCATCTT
GCGAGAAGACAAGTGCATTTCTGCGTGTGGACAAGCAGGACTTCAACCGTATCATCAAGGATGTGGAG
GCAAAGACCATGCGGCTGGAAGAATGGCAAAGTGGTGTGGTGTGGAGAGAGCCTCTCAGGGCGCCG
GCCCTTCCCGACCCCAACCCAGGCAGGAACCGGATACAGTGTGCTGGCACCCAGAGAAGATCCT
AGAGCTTCTGTTGAGGGCCATGGGACCAGATTCCAGTGTCTATGACCAACAGAGACATCTCAGCGAC
TTCTCTGACCCACAGGGTCTTCATGCCAGCGCCCAACTCTGCGCTGCCCTTCTGACCCATCCATG
TGGAGCCTGCGGGTGGCAGCGAGCAGGAGCGCAGCACTACGTCTGCAACAAGAGGCAGCAGATCTTGGC
GCTGGTACGCGAGTGGGTGGCCCTGTATGGCTCCATGCTCCACACTGACCCGTGGCCACCAGCTTCTC
CAGAAACTCTCAGACCTGGTGGGAGGACACCCGACTCAGCAACCTGCTGAGGGAGCAGTGGCCAGAGA
GGCGGGGATGCCACAGGTTGGAGAATGGCTGTGGGAATGCATCTCCTCAGATGAAGGCCCGGAATTGCC
TGTTTGGCTCCCCAACAGGACGAGCCCTTCTGGCAGCAGCTGTGCCATCCAAGTTGGGGATAAAGTC
CCCTATGACATCTGCCGGCCAGACCACTCAGTGTGACCCCTGCAGCTGCCTGTGACAGCCTCCGTGAGAG
AGGTGATGGCAGCGTTGGCCAGGAGGATGGCTGGACCAAGGGGAGGTGCTGGTGAAGGTCAATTCTGC
AGGTGATGCCATTGGCCGCAGCCAGATGCCCGTGGTGTGGCCACATCTCTGGGGCTCAATGAGCGTCTC
TTTGTGTCAACCCACAGGAAGTGCATGAGCTGATCCCACACCCTGACCCAGCTGGGGCCACTGTGGGCT
CTGCTGAGGGGCTGGACCTGGTGTGAGTGCCAAGGACCTGGCAGGCCAGCTGACGGACCAGACTGGAGCCT
CTTCAACAGTATCCACCAGGTGGAGCTGATCCACTATGTGCTGGGCCCCCAGCATCTGCGGGATGTCAAC
ACCGCCAACCTGGAGCGCTTCATGCGCCGCTTCAATGAGCTGCAGTACTGGGTGGCCACCAGCTGTGTC
TCTGCCCCGTGCCCGGCCCGGGCCAGCTGCTCAGGAAGTTCATTAAGCTGGCGGCCACCTCAAGGA
GCAGAAGAATCTCAATTCCTTCTTTGCCGTGATGTTGGCCTCAGCAACTCGGCCATCAGCCGCCTAGCC
CACACTGGGAGCGGCTGCCTCACAAGTCCGGAAGCTGTAAGTCCGCCCCGAGAGGCTGCTGGATCCCT
CATGGAACCACCGGTATACCGACTGGCCCTCGCCAAGCTCTCCCTCCTGTCATCCCCTCATGCCCCCT
TCTTCTCAAAGACATGACCTTCAATTCATGAGGAAACCACACACTAGTGGAGAATCTCATCAACTTTGAG
AAGATGAGAATGATGGCCAGAGCCGCGGATGCTGCACCACTGCCGAAGCCACAACCTGTGCCTCTCT
CACCACCTCAGAAGCCGAGTTTCCACCTCCACGAGGACAGCCAGGTGGCGAGGATTTCCACATGCTCGGA
GCAGTCCCTGAGCACCCGAGTCCAGCCAGCACCTGGGCTTATGTCCAGCAGCTGAAGGTCATTGACAAC
CAGCGGGAACCTCTCCGCTCTCCGAGAGCTGGAGCCATGA
    
```

Restriction Sites: Please inquire
ACCN: NM_001098531

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<p>NM_001098531.1, NP_001092001.1</p>
RefSeq Size:	<p>6239 bp</p>
RefSeq ORF:	<p>2772 bp</p>
Locus ID:	<p>10411</p>
Cytogenetics:	<p>12q13.11</p>
Protein Pathways:	<p>Leukocyte transendothelial migration, Long-term potentiation</p>
Gene Summary:	<p>Guanine nucleotide exchange factor (GEF) for RAP1A and RAP2A small GTPases that is activated by binding cAMP. Through simultaneous binding of PDE3B to RAPGEF3 and PIK3R6 is assembled in a signaling complex in which it activates the PI3K gamma complex and which is involved in angiogenesis. Plays a role in the modulation of the cAMP-induced dynamic control of endothelial barrier function through a pathway that is independent on Rho-mediated signaling. Required for the actin rearrangement at cell-cell junctions, such as stress fibers and junctional actin.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (a).</p>