

## Product datasheet for SC316517

### FOXD4L2 (NM\_001099279) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	FOXD4L2 (NM_001099279) Human Untagged Clone
Tag:	Tag Free
Symbol:	FOXD4L2
Synonyms:	forkhead box D4-like 2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316517 representing NM_001099279. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACCTGCCAAGAGCTGAGCGCCTCGCTCCACACCGCAGCGCAGCCTCCGGGACTCCGATGGGGAA
GACGGTAAAATCGATGTCCTGGGAGAGGAGGAAGATGAAGACGAGGTGGAAGACGAGGAGGAGGAGCG
AGACAGCAGTTCCCTAGAGCAGTCGCTCCAGCCGGGGCTGCAGGTGGCCCGGTGGGGCGGGTTGCGCTT
CCCCGAGAGCACATCGAGGGCGGGCGGGCCCGAGTGACCCTTCAGAGTTTGGCACCAAGTTCAGGGCA
CCGCCAAGGTCTCGGGCGCCTCTGAAGATGCCCGGCAGCCGGCAAAGCCCCCTACTCGTACATCGCG
CTCATCACCATGGCCATCTGCAAAACCGCACAAAGCGCCTCACGCTCAGCGGCATCTGCGCCTTATT
AGTGGCCGCTTCCCATACTACCGCCGCAAGTTCGCCCGCTGGCAGAACAGCATCCGCCACAACCTCTCG
CTGAACGACTGCTTCTGTTAAGATCCCCCGCAGCCGGGCCACCCAGGCAAGGGCAACTACTGGAGCCTG
GACCCCGCCTCCCAGGACATGTTTCGACAATGGCAGCTTTCTCCGGCGTAGGAAGCGTTTCAAGCGCCAC
CAACTGACCCCGGGAGCCACCTGCCCCACCCCTTCCCTTACCTGCTGCACACGCGCCCTGCACAAC
CCCCACCCAGGCCCTCTGCTTGGGGCCCTGCCCGCCGAGCCAGTCCCGGGGGCTACCCCAACACC
GCCCCCGGAGACGCCCTTACGCTCTGCTGCACCCGCATCCTCTTCGCTACCTACTGCTCTCGGCCCGC
GTCTATGCCGGGGCACCAAGAAAGCAGAAGGCGGGACCTGGCGACCCCGGCACCCTTCCCGTGTGC
AGCCCTCACTTGGTCTCAGCCTTGGGAGGAGGCAAGGGTCTGGCGTCGCCACCCGGGAGGCGGATGCT
TCTTTTCAGCATTGAGAGTATTATGCAAGGGGTGAGGGGAGCGGGTACAGGGGGCTGCGCAGAGTTTGT
CCCCGACCGCGTGGAGCTACTGCCACCTGCTCCAGCGACCATCAAGCCTGTTGCATCCCCAAACCGCTG
CCCTTTGCTGCAAGTGTCCGCCCGCTGCTGCTCGGACAATTTGCAGCAATAGCAGCAGCATCAGG
AGGACTGCGCCAACGGCTGCGCTCCACCAAGGGCGCGGTGCTGGGCGGGCACCTGTCGGCTCTGCGG
CGCTGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001099279

**Insert Size:** 1251 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001099279.1](#)

**RefSeq Size:** 3109 bp

**RefSeq ORF:** 1251 bp

**Locus ID:** 100036519

Cytogenetics: 9p12

MW: 45.9 kDa