

## Product datasheet for **SC316538**

### ZBTB34 (NM\_001099270) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZBTB34 (NM_001099270) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZBTB34
Synonyms:	ZNF918
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC316538 representing NM\_001099270.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGACAGCAGCAGTTTTATTTCAGTTTGTATGTGCCCGAGTACAGCAGCACCGTTCTGAGCCAGCTAAAC
GAATCCGCTGCAGGGAAACTATGTGACATCATTGTACACATTCAGGGTCAGCCATTCGAGCCCAC
AAAGCAGTCCTTGTGCCAGCTCCCATATTTCCGGGACCATTACAGCGTTAAGTACCATGAGTGGCTTG
TCAATATCAGTGATTAATAATCCCAATGTGTTTGAGCAGTTGCTTTCTTTTTGTTACACTGGAAGATG
TCCTTGACAGTGAAGGATGTTGTGAGTTTTCTGACTGCAGCCAGCTTTCTTCAGATGCAGTGTGTCATT
GACAAGTGCACGCAGATCCTAGAGAGCATCCATTCCAAAATCAGCGTTGGAGATGTTGACTCTGTTACC
GTCGGTGTGAAGAGAATCCCGAGAGTCGAAACGGAGTGAAGACAGCAGCTTCTTTGCCAACCCAGTG
GAGATCTCTCCTCCATATTGCTCTCAGGGACGGCAGCCACCAGCAAGCAGTGACCTCCGGATGGAGACG
ACCCCAGCAAAGCTTTGCGCAGCCGCTTACAGGAGGAGGGGCACTCAGACCGGGGAGCAGTGGGAGC
GTTTCTGAATATGAGATTCAGATAGAGGGAGACCATGAGCAAGGAGACCTATTGGTGAGGGAGAGCCAG
ATCACCGAGGTGAAAGTGAAGATGGAGAAGTCCGACCGGCCAGCTGTTCCGACAGCTCCTCCTGGGT
GACGATGGGTACCACACCGAGATGGTTGATGGGAAACAAGTTGTGGCAGTGAATGTGGGCTCCTATGGT
TCTGTGCTCCAGCACGCATACTCCTATTCCAAGCAGCCTCACAGCCAACCAATGTATCAGAAGCTTTT
GGAAGTTTGAGTAATCCAGCCCATCCAGGTCCATGCTGAGCTGTTTCCGAGGAGGGCGTCCCGCCAG
AAGCGGGCTTTGTCTGTCCACCTGCACAGTGACCTGCAGGGCCTGGTGCAGGGCTCTGACAGTGAAGCC
ATGATGAACAACCCGGGTATGAGAGCAGTCCCGGGAGAGGAGTGCAGAGGGCATTGGTACCCGTAC
AATGAGAGGTTGATCTGTATTTACTGTGAAAGTCCTTCAACCAGAAAGGAAGCCTTGATAGGCACATG
CGACTCCATATGGGAATCACCCCTTTGTGTGCAAGTTCTGTGGGAAGAAGTACACACGGAAGACCAA
CTGGAGTACCACATCCGGGGCCATACAGATGATAAACCATTCGCTGTGAGATCTCGGGGAAGTGCTTT
CCATTCCAAGGTACCCTCAACCAGCACTTGCAGAAAAACCAACCAGGCGTTGCTGAAGTCAGGAGTCGC
ATTGAGTCCCGGAGAGAACAGATGTGTACGTGGAACAGAAACTAGAAAATGACGCATCGGCCTCAGAG
ATGGGCCATAGATCCCGGATGGAATTCACACAGTGTCTGATGCTCCCGATTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001099270

<b>Insert Size:</b>	1503 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001099270.1</a></u>
<b>RefSeq Size:</b>	6577 bp
<b>RefSeq ORF:</b>	1503 bp
<b>Locus ID:</b>	403341
<b>UniProt ID:</b>	<u><a href="#">Q8NCN2</a></u>
<b>Cytogenetics:</b>	9q33.3
<b>MW:</b>	55.5 kDa
<b>Gene Summary:</b>	May be a transcriptional repressor.[UniProtKB/Swiss-Prot Function]