

Product datasheet for **SC316544**

KIAA1430 (CFAP97) (NM_020827) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIAA1430 (CFAP97) (NM_020827) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIAA1430
Synonyms:	hmw; KIAA1430
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC316544 representing NM_020827.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATCAGTTTGGAGATATATTAGAAGGTGAAGTGGACCATTCTTTCTTTGACAGTGACTTTGAAGAA
GGAAAGAAATGTGAACTAACTCAGTTTTTGACAAGCAAAATGATGACCCAAAGGAAAGATAGATAAA
GATACAAAAAATGTAAATTCGAACACTGGAATGCAAAACAGAAAATTATCTTACTGAGAAGGGAAT
GAAAGAAACGTGAAATTTCCCCAGAACACCCCGTAGAGAATGATGTTACACAACTGTAAGTTCTTTC
TCATTGCCAGCCTCTTCAAGATCAAAAAATTGTGTGATGTACAAACAGGACTTAAATACACGTGTCC
ATTCCAAATAGAATTCCTCAAAATGTAAAAGAAGGTGAAGATGATTACTACACAGATGGAGAGGAAAGC
AGTGATGATGGGAAGAAATACCATGTGAAGTCCAAGTCCGCTAAACCATCTACTAACGTTAAAAAAGC
ATAAGGAAAAAGTATTGCAAAGTTAGCTCCTCTTCTCCTCCTCTTTATCTTCTCATCTTCAGGTTCA
GGTACAGATTGTTAGATGCAGGTCTGATAGCCATCTATCTGATTCTGCTCCGTCATCTAAGTCATCT
AAGAAACATGTATCTGGTATAACCTCCTGTACCAAAACACAAGTATAAATCAGGAATAAAATCGACA
GAAACACAGCCTTCAAGTACTACACCAAAATGTGGCCACTACCCTGAGGAGTCTGAAGATACTGTGACT
GACGTAAGTCCCTTATCAACTCCAGACATTAGCCCTCTTCAGTCTTTTGAAGTGGGCATAGCAAATGAT
CAAAAAGTGAAAATTAAGGCAAGAAAATGTGAGCCAAGAAATATGAAGATGTTGAGGATTTGAAA
AATAATTCAAAATATTTGAAAGCAGCAAAAAAGGGAAGAAAAACATGAGCCTGATGTCTCCTCAAG
TCGTCTTCAGTGTTAGACTCCAGTTTAGACCACAGACATAAACAGAAAGTCTTACATGACACAATGGAT
CTGAATCATCTCTTGAAAGCTTTTCTGCAATTAGATAAAAAAGGACCACAAAAACATCACTTTGATCAG
CCTTCAGTAGCACCCGGGAAAACTACTCTTTCACAAGAGAAGAGGTGAGACAGATCGATCGGAAAAAT
CAGAGGCTTTTGAAAGAACTGTCAAGACAGGCGGAAAGCCGGGAAGCAAAAGTACAATTCCTAGATCG
GCTGATCATCCCCAAAGTTATATCACAGTGCTCTCAACAGACAGAAGGAACAACAAAGGATTGAGAGA
GAAAACCTTGCTTTATTGAAAAGCTTGAGGCCGTGAAACCAACAGTTGGTATGAAACGTTTCAGAACAA
CTGATGGACTATCATCGCAATATGGGCTATCTCAACTCATCACCATTGTCAAGACGGGCCAGATCCACT
CTTGCCCAATATAGCCCATTAAGAGCTTCCAGGACATCCAGTGCTACGAGTGGTCTCAGTTGTAGGAGT
GAGCGATCAGCGTTGACCCCTCCAGTGCCACCCTCGAAGAAGACCTAAACCCCTAATGTCCGTACA
GCTTGTTATTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_020827

Insert Size: 1599 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_020827.2](#)

RefSeq Size: 4882 bp

RefSeq ORF: 1599 bp

Locus ID: 57587

UniProt ID: [Q9P2B7](#)

Cytogenetics: 4q35.1

MW: 59.5 kDa