

Product datasheet for **SC316556**

FDXACB1 (NM_138378) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FDXACB1 (NM_138378) Human Untagged Clone
Tag:	Tag Free
Symbol:	FDXACB1
Synonyms:	hCG_2033039
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC316556 representing NM_138378.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGGCCCTCGGCGCCTCTGTGGTTGGGAGGGGAATTTCTCCTTCGCCCGCGCTCTGAGCGAAACC
CTGGATCAGAGCACTCAACTTACCGCCACCTGCCCTCCAGCGCCGCGAGTTGGCTCGGGATCCACTG
GCCTGGGAGAATCTGCAGTGCCTGCGCGAGCGAGGTATCGATGTACGTTTCGGTGTGGAAGTGCACCCAG
CTGGCAGATGTCTTTGAACTGCACGAGAGAGAATTTGATCAAATTTATTTTCATCTTCCCGCATTGTGGA
CGCAAAGCTGGCGTAGCTAAGAACAGGGAAGTGTGCAAAATTTTCCAAAGCTGTGCAGACGTTCTT
GCAGAGGAAGGAGAAGTCCACGTGGCATTATGTAGAGGACAAGGTGGAAGTCTGCGGATAAGCCCCAG
AGAGAATGGCACAACAGTTGGCAAGTGGTTGCCATGCCAGCCCTGGGGGGCTCATTTTAAGCGACGTG
TATCCATTCAGCTGTAAAGCTGTGGCAGGTACAAGTGCAGTGGATATAGGAGTCAAGATAAGTCCTTT
CATGTAGAAGGTGCTTTGAACCATATCTTACCAGGAGCTTACCTTTGAAGTTCTCAACCCAGAATC
TTCAGGATCAAAGTGGTAACAGTGGTTTCTTTCCAGAACCAAGCAAGTGTAGGAAAGTTGAAC
AGGGGTTTCTGGAAGCACCTTCATGTATCCTATCAAAACCATAAATGAGAAACTATTGCTGAATTA
GGCAAAGTTTCCCGCTAAAAAGGCTGAAGTGTCTACCTTTGCTGCCACAGGAAGGTACCAAGTGT
CTTCTTTCTGGAATTGTGACTTTCTGTGAGTGTCTTTTGGATTAGTCTCCATGAAGATAACTCAAT
TCTGAGTCCCTGACTGGTGGGACATCACAAGATGTGGAAGACTTCTAGTGTCAATTTTCAAGACTTAGC
CTTCTCAAGAATCCTGGAAGAGATGGTAAAGAAGAAGCTGTGAAGGAACCTGTGGCCAGGCCAAGATC
TGCCTTAGACCTTCTCTCCTAGTGCATGTTTCAAGATGTCATCGAAGTACCAGACTTCTCTCAGGTTCT
CTGCACATCCTCAGTGGACCTGTCTTTCAGAAGTGCCACATTTTGCCTTTCACAATGCCAGCATTTTCAT
GAAACTTTATTTATCCTTGGGGTTAATCAAAATCTGAAGGATGGCTGTCTTCAATCACTGCTGGATCAT
CTGAAGGGCATTCTAGATAGCCTGCTGACCCAGACATTGCCGGAGAGCTCTAAGCTGAGCAGTTTATGTC
AAATTTGTCTTCAAGTCAATGGAAGGATTATGATTCTGTGAAGACTCATAATTTTAGCCAGAT
TGTAAGGATCTAATTATTGGGTCTGTTATCAGTCTGCCACTAGTGTATACATAAAGACCAAGTGT
TTTGTGTTGTGTCTATGAAGTGGACTTATTAGCCATGCTTGTCTGGTGTATCTCTGACTGGAGAATG
TTGTGGACGTTTGATAACCGTTTCTGAAAAATTTTGTCCCTGGCAAAATAGAACCCTTTAAAGTCAT
TCTCTGTATCCTCCATGTTATGTGATGATGTTAGTTTTTGGATAGATCAGAAGAAAGGATTGATGAA
CTAGAGTTTCACTGTGGCCGAGCAGTGTCTCAGGACACTATTATATCCATACAGTTTCTTAGCCGT
TTCCAGCATCAAAGACTCAACAGGTGAGTCTCTGCTATAGATTGACCTACCAGACCTGTGACAAGGCC
CTACCCAGCAGCAAGTAGCATCAATGCAGTCCCAGTTTAGGAAGGAGATTCAACAACACCTATATGTT
ATACCTCGGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_138378

Insert Size: 1875 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_138378.2</u>
RefSeq Size:	3081 bp
RefSeq ORF:	1875 bp
Locus ID:	91893
UniProt ID:	<u>Q9BRP7</u>
Cytogenetics:	11q23.1
MW:	70.4 kDa
Gene Summary:	This gene encodes a protein which contains a ferredoxin-fold anticodon-binding domain which is contained in a subset of phenylalanyl tRNA synthetases. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011] Transcript Variant: This variant (1) represents the protein-coding transcript.