

## Product datasheet for **SC316564**

### SEC14L5 (NM\_014692) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SEC14L5 (NM_014692) Human Untagged Clone
Tag:	Tag Free
Symbol:	SEC14L5
Synonyms:	PRELID4B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC316564 representing NM\_014692.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGTGCAAAGATACCAAGTCTCCTGTCCGAGTCTACAAGTACCCGTTTGTAGCTGGTTCATGGCGCCCTAC
GAGAAGCGTTTCCCACGTGCCACAGATCCCAGTCTTCCCTGGCAGCGAGGTCTTGCAGGAGTCCCGC
AGCCCGGACGGGGCTGTGCACGTGGTGGAGCGGAGCTGCCGGCTGCGCGTGGACGCCCGCGGCTGCTG
CGGAAGATCGCAGGTGTTGAGCAGTGGTCTTCTGTCAGACAAACATCTTGAAGTGAAGGAGAGGACG
CTCCTCATCGAAGCGCACAATGAGACCTTCGCAACCGCGTGGTGGTGAACGAGCACTGCAGCTACACG
GTCCACCCTGAGAATGAAGACTGGACTTGTTCGAGCAGTCTGCCTCACTGGACATTCGGTCTTTCTTT
GGCTTTGAAAATGCCTGGAGAAGATCGCCATGAAGCAGTACACCGCCAACGTCAAGAGGGGGAAGGAG
GTGATTGAGCATTACCTGAATGAGCTCATCTCCAGGGTACCTCGCACATTCGCGCTGGACGCCTGCC
CCAGTCCGTGAGGAGGATGCCCGCAACCAGGCTGGACCGAGGGACCCAGCTCCCTGGAGGCCACGGG
CCCCGTAGCACCTGGGGCCCGCTTGGAGGCGGTGAGTATGGACGGGGACAAGCTGGATGCGGACTAC
ATTGAGAGGTGCCTGGGCCACCTCACGCCATGCAGGAGAGCTGCCTGATCCAGCTTCGGCACTGGTTA
CAGGAGACCCACAAGGCAAGATTCCCAAAGATGAGCACATCCTTCGGTTCCTGCGGGCTCATGACTTC
CACCTGGACAAGGCCCGGAAATGTGCGCCAGTCTTGGAGTGGCGCAAGCAGCACCAGGTGGATCTC
CTCCTTACAGACCTGGCAACCCCTGCCCTGCTGGAGGAGTTCTATGCAGGGGGCTGGCATTACCAGGAC
ATAGATGGCCGCCCTCTACATCTCCGCTGGGCCAGATGGACACCAAAGGCTTGATGAAGGCCGTG
GGGGAGGAGGCGCTGTGCGGCATGTTCTCCGTCAACGAGGAAGGACAGAAGCGGTGTGAGGGGAGC
ACAAGGCAGTGGGCCGTCCCATCAGCTCCTGGACCTGCCTGTAGACCTGGAGGGACTCAACATGCGG
CACCTGTGGCGGCCGGGGTGAAGGCCCTGCTGCGGATGATTGAGGTGGTTGAGGACAATTACCCAGAG
ACCCTGGGTGCGGTGCTCATCGTGCAGACCCCGGAGTCTTCCCGTGTCTGGACACTGATCAGCCCC
TTCATCAATGAGAACACCAGGCGGAAGTTCCTCATCTACAGTGGCAGCAACTACCAGGGACCCGGAGGC
CTTGTTGACTATCTGGATAGAGAAGTGATCCCTGACTTCCTTGGGGGAGAGAGTGTGTGTAATGTCCCC
GAAGGAGGGCTGGTCCCCAAGTCCCTCTACATGACAGAAGAGGAGCAGGAGCACACGGACCAGCTGTGG
CAGTGGAGTGAGACCTACCATTACGCCAGCGTGTCCGCGGAGCCCCACGAGGTGGCCGTGGAGATC
CTGGAAGGAGAGTCGGTCATCACCTGGGACTTTGACATCCTGCGAGGGGACGTGGTGTTCAGCCTGTAC
CACACCAAGCAGGCGCCAGGCTGGGCGCCGGGAACCGGGGACCAGGGCCAGCGGGCAGCTGATCGAC
AAAGGCTGGTCTGGGCAGGATTACAGCCGTGTGGAGGCTCCCCTTGTCTGCCGGGAGGGGGAGAGC
ATCCAGGGCTCCCATGTGACCCGGTGGCCCGGCTCTACCTGCTCCAGTGGCAAATGCACAGCCCCCCC
AGCAGCGTGGCCTGCAGCCTCCCGGTGTGGACGATGTCCTGACGGCTCTGCACAGCCCCGGCCCAAG
TGCAAACCTTCTACTACTGTGAGGTGCTGCCTCTGAGGACTTCAGGGGCTCCATGTCCAGCCTGGAA
TCCTGCACCAGCGGCTTCTCCAGCTCAGCGCCGCCACCTCGTCTCCTCCTCCGGCCAGTCTCATAGC
AGCTCCCTGGTCTCCAGATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_014692

**Insert Size:** 2091 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_014692.1](#)

**RefSeq Size:** 6456 bp

**RefSeq ORF:** 2091 bp

**Locus ID:** 9717

**UniProt ID:** [O43304](#)

**Cytogenetics:** 16p13.3

MW: 78.9 kDa