

Product datasheet for **SC316572**

C8orf80 (NUGGC) (NM_001010906) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C8orf80 (NUGGC) (NM_001010906) Human Untagged Clone
Tag:	Tag Free
Symbol:	C8orf80
Synonyms:	C8orf80; HMFN0672; SLIP-GC
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC316572 representing NM_001010906.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGGCAGAAACGAAGGATGTTTTGGCCAGGAACCGCATCCAGTTGAAGATGATTTATATAAAGAACGA
ACGAGAAAAAGAAGGAAATCAGATCGAGACCAGCGTTCCGAGCATTTCCCTCCATGGAGCAGAGTGCT
CTTAAGGAATATGAAAAATTGGAATCACGGACCAGAAAGGTTTTGAGCAACACTTATCAGAAACTTATT
CAGTCTGTCTTCCTGGATGACAGCATCCCTAATGGAGTCAAGTATCTCATAAATAGGCTTCTTGCCCTTG
ATTGAAAAGCCGACAGTGGACCAATCTACATTGCATTATTTGGAAGCACTGGGGCTGGGAAGAGCTCC
CTGATCAATGCCATCATCCAGCAAGCAATGTTTCTACAGTGTCTGGAGAAAGCATATGACTTCTGCTGC
ATTGTACAAGTGAGCTCTGGCTGCTGTGTCAGTATGAGGCCAAAATCCACCTTCTGTCTGACCAGGAG
TGGAGGGAGGAGCTGAAGAACTGACCAAACTCCTGCATAGGACGGAGGAGCTGAGCAGAGAAGAGGCA
GATGCGTGAACAGGGATGAGGCAGTGGAGGAAGCCACCTGGAAGCTACAATGATTTATGGAATGGG
GCAGAGAGTAAGAACTATGAGGAGTTACTGAGGGCGAAGCCAAAAGGAAGATCCCCACCTCCAGAGTC
ATCACCTCAAGGCGGAAGAGGCAGAAGAGCTGTCCATCAAGCTGGACCCCTACATCCGCACACAGAGG
AGAGATTGGGATGGAGAGGCCGCTGAGATGCGCATCTGGCCCTTGATCAAACATGTGGAAGTGACACTT
CCCAATCCGACCTGATCCCAGAAGGGGTGCTGCTGGTGACATCCAGGCACAGGCGACTTCAACAGC
AAGAGGGACGAGATGTGGAAGAACATTGACAAGTGCTCAGTGATCTGGTGATCAGCGACATAGAG
CGAGTTTCTGGGGGCAAGCCACGAAGACCTTCTGAATGAGAGCATCAAAGCCTGCCAGCGGGGCTTC
TGTAGGGACGTGGCCCTGGTGGTCACCAAGATGGACAACTCCACTTGCCAGAATACCTAAGGGAAGA
AAGGCAGGAAATCAAGCTATTCAGAGTCAGCGAGAGGCTGTTTTAGAAAGAAATGAAATGATAAACTA
CAAAGGACTAGAATTCTCAAGGAAAACTGAAGAGAAAACTGCCAGCTGACTTCAAGGTTCTGGAAGCC
TCAGATCTGGTGTATACAGTCAGCGCCAGGAGTACTGGCAGCAGGCTCTCCTCACCGAGGAGAAACG
GAAATCCCTAAGTTAAGAGAATACATCAGGAAGAGCCTCTTGACAAGAAGAAGAGGACAGTGACCAAG
TATGTGACTGAAGCCTTTGGCCTGTTGCTCCTCACAGATAGTTTCAACTCCACGCAAAACCTGCCGAAT
GAACACTTGACATGAGTGTCTGCGGAGATTTGCGGAAGAGAAGGTTGAGCTGCTGGAGAAGGCCATC
GCACAGTGCTTCGCTGCATGGAGCAGCCTCTGCAAGAAGGGGTGAGGACCGCCAGGACTTCTTACCGC
TGCATCCTCAGAGCATGCTTGGTGAAGGTAAGGAAACCAAGGTTTTATCAGACCTGAAAGCTGTT
TGCTGAAAAATGGCATCTATGCCTCCAGGACTCTGGCGAGAATTGATCTAAATGAAGCCCTCACTCAG
CCCGTCTATGACCAGATCGACCCTGTTTTTGAAGCATTTTTAGGACGGGGAAGCCCACTGGTTCAGCT
CTGATGCCTCACATAGATGCTTTTAAGCAGTCCCTGCAGGAGAAAATGACAGAAATTGGGATAAGAAGT
GGCTGGAATATGATAGCTGCAAAAAAATTTCTGATCCAGGAGATAAGTGCCATCCTCGGGGGCTG
GAGGACCACATCCTCAGAAGGAAGAGGAGGATCTACGAGTCCCTCACTGCCTCTGTCCAGAGTGACCTG
AAGCTCTGTTACGAAGAGGCAGCTCAGATCACGGGCAAAAAGCGTGTGAGCGGATGAAAGATGCCATC
AGAAGAGGAGTGGACCGGAGGTGGCTGAGGGCATGTTTGAAAGGGCCAGGAAAGGATGCAGCACCAG
TTTCAGCAGCTGAAGACTGGAATCGTGGAGAAGGTGAAGGGCAGCATCACCCTATGCTGGCCCTTGCT
TCGTCCCAGGGGGATGGCCTCTACAAGGAGCTTGACAGTGTGGGAGTGAATACAAGGAGATGGAGAAG
CTGCACAGAAGCCTGAGGGAGGTGCGGGAAGTGCACGGCTGAGGAAGGGCATGCAAGAATTCTCTCTA
AGGGCATCCCCAGCAAGGCTGGCCCCCGGGACATCACTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001010906
Insert Size: 2391 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001010906.1</u>
RefSeq Size:	3912 bp
RefSeq ORF:	2391 bp
Locus ID:	389643
UniProt ID:	<u>Q68CJ6</u>
Cytogenetics:	8p21.1
MW:	91.1 kDa
Gene Summary:	Nuclear GTPase found in germinal center B-cells, where it may inhibit function of the activation-induced cytidine deaminase AICDA (PubMed:19734146). Reduces somatic hypermutation in B-cells which may enhance genome stability (By similarity). [UniProtKB/Swiss-Prot Function]