

Product datasheet for SC316633

ST20 (NM_001100879) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ST20 (NM_001100879) Human Untagged Clone
Tag:	Tag Free
Symbol:	ST20
Synonyms:	HCCS-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316633 representing NM_001100879. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGCGCGATCTCGGCTCACTGCAACCTCTGTCTCCAGGTTTCAGGAAAATGGCTTTGTAAAGAAGCTT
 GAGCCTAAATCTGGCTGGATGACTTTTCTAGAAGTTACAGGAAAGATCTGTGAAATGCTCTTCTGTCCCT
 GAAGCAATACTGTTGACCAGAAAGGACACTCCATATTGTGAAACCGGCCTAATTTTCTGACTCTTACG
 AAAACGATTGCCAACACATACTTCTACTTTAA
 ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001100879
Insert Size:	240 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001100879.1</u>
RefSeq Size:	1025 bp
RefSeq ORF:	240 bp
Locus ID:	400410
Cytogenetics:	15q25.1
MW:	9 kDa
Gene Summary:	<p>May act as a tumor suppressor. Promotes apoptosis of cancer cells.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript. Variants 1, 2 and 3 encode the same protein.</p>