

Product datasheet for **SC316668**

CYBC1 (NM_001100408) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CYBC1 (NM_001100408) Human Untagged Clone
Tag:	Tag Free
Symbol:	CYBC1
Synonyms:	C17orf62; CGD5; Eros
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316668 representing NM_001100408. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTACCTGCAGGTGGAGACCCGCACCAGCTCCCGCTCCATCTGAAGAGGGCTCCAGGCATCCGGTCC
TGGTCCCTGCTGGTTGATAGCCTGGGCTGGAAGCTCTTCTACGTCACAGGCTGCCTGTTTGTGGCTGTG
CAGAACTGGAGGACTGGGAGGAAGCCATCTCGACAAGAGCACAGGGAAGGTTGTTTTGAAGACGTTT
AGCCTCTACAAGAAGCTGCTGACTCTTTTCAGAGCTGGCCACGACCAGGTGGTGGTCTGCTCCATGAT
GTCCGTGATGTGAGCGTGGAGGAGGAGAAGGTCCGGTACTTCGGGAAAGGCTACATGGTGGTCTCCGG
CTTGCGACGGGCTTCTCCACCCCTCACGCAGAGTGCAGTCATGGCCACCGCAGTGATGTGGAAGCC
ATCGCCAAGCTCATCACCAGCTTCTGGAGCTGCACTGCCTTGAGAGCCACAGAGCTGTCTCAGAGC
AGCGACAGTGAGGCCGGTGACCCGCAAGCCAGAGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites:	Sgfl-MluI
Plasmid Map:	<input type="checkbox"/>
ACCN:	NM_001100408
Insert Size:	522 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001100408.2
RefSeq Size:	2125 bp
RefSeq ORF:	522 bp
Locus ID:	79415
UniProt ID:	Q9BQA9
Cytogenetics:	17q25.3
Protein Families:	Transmembrane
MW:	19.4 kDa
Gene Summary:	<p>Necessary for a stable expression of the CYBA and CYBB subunits of the cytochrome b-245 hetrodimer. Controls the phagocyte respiratory burst and is essential for innate immunity. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) lacks an alternate in-frame exon, compared to variant 1. The resulting protein (isoform b) has the same N- and C-termini but is shorter compared to isoform a.</p>