

Product datasheet for **SC316718**

RALYL (NM_001100391) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RALYL (NM_001100391) Human Untagged Clone
Tag:	Tag Free
Symbol:	RALYL
Synonyms:	HNRPCL3
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF:	>OriGene sequence for NM_001100391 edited <pre>GTAGAAGCGGCAGAGACCGCATATTCGGAGGGGCTGGCTCCAGGCCACGCGAGCCGGAG CTGGAGACATTTCCAGAGCCTTGAGGATGGCACCAGGCCCTCGCACGCTCAACTCCTGCTC CTCCTCTTCGCAATGAGTAAACGACGCAGATTAAGCAAGGAGAGCCAATCATGACTGGC AAAACACAGACCAGCAACGTACCAATAAGAATGACCCCAAGTCCATCAACTCCCGTGTT TTCATCGGCAATCTAAATACGGCAATTGTCAAGAAAGTTGACATTGAAGCCATTTTTTCA AAGTATGGAAAAATAGTTGGATGTTCCGTTACAAAGGTTATGCATTTGTACAGTACATG AGTGAGCGACATGCAAGAGCTGCAGTGGCTGGAGAAAATGCCAGAGTCATCGCCGGCCAA CCACTTGATATCAACATGGCAGGAGAGCCAAACCATACAGACCAAAACCTGGAAACAAG AGGCCCTTTCTGCACTTTACAGACTTGAATCAAAGGAACCTTTCTGTCTGTTGGCGGT TATGTCTTTGACTATGATTACTACAGAGATGATTTCTACAATCGTTATTTGATTACCAC GGGCGTGTGCCTCCACCTCCCGTGCAGTAATTCGCTGAAGCGTCCAGAGTGGCAGTC ACAACGACTCGCAGGGGAAAGGAGTCTTTTCATGAAAGGTGGATCGAGATCTACTGCC AGTGGGTCAACAGGTTCTAAATTGAAATCAGATGAGTTACAGACCATCAAGAAAGAATTA ACCCAGATCAAACTAAATTTGACTCCTTGCTAGGGCGCTGGAGAAGATTGAGAAACAG CAGAAGCGGAGGCAGAAAGCTCAGAAGAAGCAATTGGAAGAGAGTCTAGTGCTGATCCAA GAGGAATGTGTGCAGAGATTGCAGATCACTCTACAGAGGAGCCTGCTGAAGGAGGGCCA GATGCCGATGGAGAAGAGATGACAGATGGGATAGAGGAGGACTTCGATGAAGATGGGGGT CATGAGCTGTTTCTACAGATAAAGTGATCTGAAATAACGCATGATGCCACAAAGCAGAAA AGAGAACTGTGACAACCCCCAGAAATGTGAAAGGAGGTTTCTTACTGGACAGCAGCATC TTTGGTTCAATTTATATAAAACCCAAATAAATAAATGGACAGTATTGCTCAGTTTTAG AAATTCATTTCTCTATGTTTTAAGCTGTACAATTGTCAGGTTTTATGGTTTAAATTG TAAATGTGTTTTCCCTTTGCTAATTATG: TTTTTTTTTTCAGTCTTAAATGTGAAAGG CATTTATGAATGTAAGGGAACACTATATACAATGTATATTTGTAAAGCTATTTTTTA TGATTAGCATGTTTCACTGTTGATCATATATAAAGTCAGGTGATTTGCAATTCTGTATT TAAAGCTTATTTCCAACAATGTCATGTAAGAAAAGATGCATCTTATGCTAGTTTTTATAA TTTATTTATAATTTATAGTTTAAAGTACTTCAGATCATAATGATAAAATACTTGAAAAAG TTATATTTCTGCCCTGTATAAGCACCTTTTTATTAATAAAGAATGCAGATATTTAGAT GTGATATAATAGTTAAAGAACTGTTGGTTTGATCTGTGATTAAGTTGAGCATGCTCCGCT CTACTGAACTAAATGATCCAATTATTACTTCAGTCTGGGTATGAGATTCCATGGACAAGT AAGGACTAGATTGCCAAGGAAAAGACTGTCTTGCCTTGATCCAAAAGTTTAAATTAGT GCATACATCATGTCATTTACCTCCTGTTCTAGGAACTCTCATTCCCAAGCATTGCCA GTGTTTTCCAGATAATCTTAGCTGTTGTCTGTGCTGTGGAAATGGAAGAAACCATCTTC ACAGACTGTAGGAGAATTCAACATATAATTTCTTAATAAATACTGTTTCTTTAAACAA AAAAAAAAAAAAA</pre>
Restriction Sites:	Please inquire
ACCN:	NM_001100391
Insert Size:	2100 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001100391.1](#), [NP_001093861.1](#)

RefSeq Size: 1995 bp

RefSeq ORF: 915 bp

Locus ID: 138046

UniProt ID: [Q86SE5](#)

Cytogenetics: 8q21.2