

Product datasheet for **SC316774**

PHACTR2 (NM_001100166) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PHACTR2 (NM_001100166) Human Untagged Clone
Tag:	Tag Free
Symbol:	PHACTR2
Synonyms:	C6orf56
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001100166, the custom clone sequence may differ by one or more nucleotides

```

ATGGACAACGCCGTTGACGGACTGGACAAAGCTTCTATAGCAAACCTCAGATGGCCCCACA
GCAGGTTCCAAACACCTCCCTTCAAAGAAAGGGGAAACTATCCACCATTGGTAAAATC
TTAAGCCTTGAAATGGAGGAAAAAGAAGACCAGCGACAAATTTAGAGAAACCTCGGCA
GTATTAGAAAGGAAGATATCCACACGACAAAGTAGAGAGGAGCTGATAAGAAGGGGAGTG
CTTAAGGAATTGCCTGATCAAGATGGAGATGTAACAGTTAACTTTGAAAATTCAAACGGG
CACATGATACCCATCGGAGAGGAATCTACCCGAGAGGAAAATGTAGTAAAGTCTGAAGAA
GGTAATGGCTCTGTATCTGAAAAACACCACCTCTGGAGGAACAGGCAGAAGATAAGAAA
GCTGGCTCCTCTATTCAAAAAACAACCTGGCTCTAAAGCATCAGTTGCCCATCCACT
TCATCCACCTCATCTCGTCCCAAAGCTTCAAAGGAGACAGTTTCTAGCAAAGCAGGGACA
GTGGGGACCACCAAGGGCAAGAGAAAAACTGACAAGCAGCCAATAACTTCTCACCTGTCC
TCAGACACAACAACCTTCTGGCACATCCGACCTGAAAGGAGAGCCTGCAGAGACCAGAGTG
GAGAGTTTCAAACCGAACAGACTGTCCCTGGAGCTGAGGAGCAGAACACAGGCAAAATTC
AAGTCCATGGTCCCTCCACCCCTGTGGCTCCAGCACCTTCTCCTCTGGCCCCCCTCTC
CCTCTTGAGGATCAGTGCATTAAGTGCCTCAGACACTCCAGTTGTCCTCGTCAGCGTTGGA
GCTGACCTGCCGCTCTCTGCCTTAGACCAAGTCAGCTTCTTTGGGCTGAAGAGCCGAGC
AACAGAACCACCTCTACTCAGGCACTGGCTTAAGTGTTAACAGAGAAAATGCAAAATGT
TTCACAACCAAAGAGGAGCTGGGGAAGACAGTGCCTCAGCTACTGACTCCTGGGCTGATG
GGCGAATCTTCAGAATCCTTTAGTGCCTCAGAAGATGAAGGCCACAGGGAATACCAAGCC
AATGACTCTGACTCGGACGGGCTATCTTGTACACCGATGATGAGGACGAAGACGAAGAT
GAGGATGGCAGTGGAGAAAGTGCTTTGGCAAGTAAAATACGCCGAGGGATACTCTTGCT
ATCAAACCTGGCAACAGACCATTAAGAAAGAACTAGAGGACAAAACATCTTGACGCGT
ACATCTGAAGAAGAGAGGAGGAAATCCGACAACAATTTGGAACCAAGTTAGTCAGGAGG
CTGAGCCAGAGGCCCAACTGAAGAATTAGAGCAAAGAAACATCCTAAAACAAAAGAAT
GAAGAAGAGGAACAAGAAGCAAAAATGGAACCTTAACCGCAGACTCAGCAGAAAGCTCAGC
CTGAGACCCACAGTGGCAGAGCTACAAGCTCGAAGGATCCTGCGATTTAACGAGTATGTA
GAAGTCACGGATTCTCCTGACTATGACCGCCGAGCAGACAAGCCCTGGGCCAGGTTAAACA
CCTGCAGACAAGGCAGCAATAAGGAAAGAACTAAATGAATTTAAAAGCACAGAAATGGAA
GTTTCATGAAGAGAGTCGACAGTTTACAAGGTTTCATCGTCCA

```



[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_001100166
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001100166.1</u> , <u>NP_001093636.1</u>
RefSeq Size:	9300 bp
RefSeq ORF:	1665 bp
Locus ID:	9749
UniProt ID:	<u>O75167</u>
Cytogenetics:	6q24.2