

Product datasheet for **SC316789**

GRB10 (NM_005311) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GRB10 (NM_005311) Human Untagged Clone
Tag:	Tag Free
Symbol:	GRB10
Synonyms:	Grb-10; GRB-IR; IRBP; MEG1; RSS
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene ORF within SC316789 sequence for NM_005311 edited (data generated by NextGen Sequencing)

```

ATGGCTTTAGCCGGCTGCCAGATTCTTTTTGCACCATCCGTACTIONACCAGGACAAGGTGGAGCAGACAC
CTCGCAGTCAACAAGACCCGGCAGGACCAGGACTCCCCGCACAGTCTGACCGACTTGCGAATCACCAGGA
GGATGATGTGGACTGGAAGCCCTGGTGAACGATATGAATGCATCCCTGGAGAGCCTGTACTCGGCCTGC
AGCATGCAGTCAGACACGGTGCCCTCCTGCAGAATGGCCAGCATGCCCGCAGCCAGCCTCGGGCTTCAG
GCCCTCCTCGGTCATCCAGCCACAGGTGTCCCAAGGCAGAGGGTGCAGCGCTCCAGCCTGTGCACAT
CCTCGCTGTGAGGCGCTTCAGGAGGAAGACCAGCAGTTTAGAACCCTCATCTCTGCCGCCATCCCCAAT
CCTTTTCTGAACTCTGTGGCCCTGGGAGCCCCCTGTGCTCACGCCGGTTCTTTACCTCCGAGCCAGG
CCGCCGAAAGCAGGATGTTAAAGTCTTAGTGAAGATGGGACAAGCAAAGTGGTGGAGATTCTAGCAGA
CATGACAGCCAGAGACCTGTGCCAATTGCTGGTTACAAAAGTCACTGTGTGGATGACAACAGCTGGACA
CTAGTGGAGCACCACCCGACCTAGGATTAGAGAGGTGCTTGAAGACCATGAGCTGGTGGTCCAGGTGG
AGAGTACCATGGCCAGTGAGAGTAAATTTCTATTAGGAAGAATTACGAAAATACGAGTCTTTAAAAA
TCCCATGAATTTCTCCAGAACAGATGGTTACTTGGTGGCAGCAGTCAATGGCAGTCAAACCCAGCTT
TTGCGAATTTCTGAACTCCAGTAGTTGTCTGAAATCAAGGGTTTTGCATGTGAAAGAGCTGGAA
AGAAATCATGAAAAAGCTGTATGTGTGTTGCGGAGATCTGGCCTTTATTGCTCCACCAAGGAACTTC
AAAGGAACCCAGACACCTGCAGCTGCTGGCCGACTGGAGGACAGCAACATCTTCTCCCTGATCGCTGGC
AGGAAGCAGTACAACGCCCTACAGACCACGGGCTCTGCATAAAGCCAAACAAAGTCAGGAATGAAACTA
AAGAGCTGAGGTTGCTCTGTGCAGAGGACGAGCAAACCAGGACGTGCTGGATGACAGCGTTCAGACTCCT
CAAGATGGAATGCTCCTTTACCAGAATTACCGAATCCCTCAGCAGAGGAAGGCCCTGCTGTCCCGTTC
TCGACGCCAGTGGCAGTGTCTCCGAGAACTCCCTCGTGCAATGGATTTTCTGGGCAACAGGACGCG
TGATAGAGAATCCGGCGGAGGCCAGAGCGCAGCCCTGGAGGAGGGCCACGCCTGGAGGAAGCGAAGCAC
ACGGATGAACATCCTAGGTAGCCAAAGTCCCTCCACCCTTCTACCCTAAGTACAGTGATTACAGGACA
CAGCACTGGTTTACGGGAGGATCTCCAGGAGGAATCCACAGGATCATTAAACAGCAAGGGCTCGTGG
ATGGGCTTTTTCTCCTCCGTGACAGCCAGAGTAATCCAAAGGCATTTGTACTCACACTGTGCATACCA
GAAAATAAAAATTTCCAGATCTTACCTTGGCAGGACGACGGGCAGACGTTCTTCAGCCTAGATGACGGG
AACACCAAATCTCTGACCTGATCCAGCTGGTTGACTTTTACCAGCTGAACAAAGGAGTCTGCCTTGCA
AACTCAAGCACCCTGCATCCGAGTGGCCTTATGA
    
```

Clone variation with respect to NM_005311.4
315:a=>g

- Restriction Sites:** Please inquire
- ACCN:** NM_005311
- Insert Size:** 2800 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_005311.3 , NP_005302.3
RefSeq Size:	4713 bp
RefSeq ORF:	1785 bp
Locus ID:	2887
UniProt ID:	Q13322
Cytogenetics:	7p12.1
Domains:	RA, SH2, PH
Protein Families:	Druggable Genome
Gene Summary:	<p>The product of this gene belongs to a small family of adapter proteins that are known to interact with a number of receptor tyrosine kinases and signaling molecules. This gene encodes a growth factor receptor-binding protein that interacts with insulin receptors and insulin-like growth-factor receptors. Overexpression of some isoforms of the encoded protein inhibits tyrosine kinase activity and results in growth suppression. This gene is imprinted in a highly isoform- and tissue-specific manner, with expression observed from the paternal allele in the brain, and from the maternal allele in the placental trophoblasts. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Oct 2010]</p> <p>Transcript Variant: This variant (1), also known as hGrb10gamma and hGrb10zeta, encodes the longest isoform (a).</p>